



Quickie M6™
May 2025
 EIF1

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ Prov: _____ PC: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ Prov: _____ PC: _____
 Ship To Phone: () _____ Fax: () _____

✓ Option # Definition Price ✓ Option # Definition Price

1. FRAME

Base Model

<input type="checkbox"/>	115WM21	Quickie M6	\$4,811.00
--------------------------	---------	------------	------------

2. SIZE / SEATING

Frame Width

<input type="checkbox"/>	115W13	22"	NC
<input type="checkbox"/>	115W14	24"	NC
<input type="checkbox"/>	115W15	26"	NC
<input type="checkbox"/>	115W16	28"	NC
<input type="checkbox"/>	115W17	30"	NC

Frame Length

<input type="checkbox"/>	115FL9	Regular (18"-19" depth)	NC
<input type="checkbox"/>	115FL18	Long (18"-22" depth)	NC

Seat Depth

<input type="checkbox"/>	115SD11	18"	NC
<input type="checkbox"/>	115SD12	19"	NC
<input type="checkbox"/>	115SD13	20" (Long frame only)	\$308.00
<input type="checkbox"/>	115SD14	21" (Long frame only)	\$308.00
<input type="checkbox"/>	115SD33	22" (Long frame only)	\$308.00

Seat Height

<input type="checkbox"/>	115HT65	17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster	\$560.00
<input type="checkbox"/>	115HT11	18.5" Only available w/ 5" casters	NC
<input type="checkbox"/>	115HT14	20"	NC

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	115JC20	JAY Bundled Cushion Cushion Name or Part # _____ Width: _____ Depth: _____	Price List
--------------------------	---------	--	------------

3. CASTERS

Castors

<input type="checkbox"/>	115CW6	5" Soft Roll	\$126.00
<input type="checkbox"/>	115CW36	8" x 2" Pneumatic	NC
<input type="checkbox"/>	115CW37	8" x 2" Pneumatic w/Airless Insert	\$393.00

Seat Height Matrix

Caster	Front/Rear
5"	17" (Different Axle & Fork)
5"	18.5"
8"	20"

4. REAR WHEELS

Rear Wheel Type

<input type="checkbox"/>	115RW2	24" Spoke Wheel	NC
<input type="checkbox"/>	115RW29	24" Spinergy Wheel	\$1,878.00

Rear Wheel Tire

<input type="checkbox"/>	115RT1	Pneumatic Tire	NC
<input type="checkbox"/>	115RT2	Pneumatic Airless Insert Tire	\$247.00
<input type="checkbox"/>	115RT25	Primo Orion Poly	\$271.00

Handrim

<input type="checkbox"/>	115HR1	Aluminum Anodized	NC
<input type="checkbox"/>	115HR4	Plastic Coated	\$146.00
<input type="checkbox"/>	115HR7	Omit Handrim	NC

Adjustable Axle Plate

<input type="checkbox"/>	115AL2	Standard (for 18.5", 20" Seat Ht.)	NC
<input type="checkbox"/>	115AL28	Low (only for 17" Seat Ht.)	NC

Axle

<input type="checkbox"/>	115AX1	Quick Release	NC
--------------------------	--------	---------------	----

Rear Tire & Wheel Matrix

Tire	Rear Wheel	
	Spoke wheel	Spinergy wheel
Pneumatic	24	24
Pneumatic w/ Airless Insert	24	24
Primo Orion Poly	24	24

5. WHEEL LOCKS

Wheel Lock				Extension Handle			
<input type="checkbox"/>	115WL1	High Mount Push Wheel Lock	NC	<input type="checkbox"/>	115EH1	6" Pair	\$73.00
<input type="checkbox"/>	115WL12	Do Not Mount Wheel Lock	NC				
<input type="checkbox"/>	115WL10	Omit Wheel Lock	NC				

6. HANGERS AND FOOTREST

Front Hangers				Extension Tubes			
<input type="checkbox"/>	115H1	60° Swing Away	NC	<input type="checkbox"/>	115E1	Short (13"-15.5")	NC
<input type="checkbox"/>	115H146	70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$193.00	<input type="checkbox"/>	115E2	Medium (15.5"-18")	NC
<input type="checkbox"/>	115H132	80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$188.00	<input type="checkbox"/>	115E5	Long (18"-20.5") NA w/ 5" caster	NC
<input type="checkbox"/>	115H41	ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>	\$1,043.00				

Lower Leg Length Range Matrix

Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"
	Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"
	Long	16.00" - 16.50"	16.00" - 17.75"	NA
70°	No Ext Tube	NA	9.75" - 13.50"	NA
80°	No Ext Tube	NA	10.75" - 14.75"	NA
ELR	Short	NA	NA	13.00" - 17.00"
	Medium	NA	NA	14.25" - 17.00"

All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle

Footplates

Heel loop or leg strap not included, See Footplate Options Section

<input type="checkbox"/>	115F1	Composite	NC
<input type="checkbox"/>	115F25	Adult Angle Adj	\$168.00
<input type="checkbox"/>	115F5	Aluminum	NC

Footplate Options

<input type="checkbox"/>	115HL1	Heel Loop	NC
--------------------------	--------	-----------	----

7. BACKREST

Backrest				Back Height Preset			
<input type="checkbox"/>	115BT34	8° w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i>	NC	<input type="checkbox"/>	115B91	16-3/4"	NC
				<input type="checkbox"/>	115B96	17-7/8"	NC
				<input type="checkbox"/>	115B39	19"	NC
				<input type="checkbox"/>	115B98	20-1/8"	NC

8. ARMREST

Armrest			
<input type="checkbox"/>	115AA33	Dual Post Height Adj - Desk 10"	\$381.00
<input type="checkbox"/>	115AA34	Dual Post Height Adj - Full 14"	\$381.00

9. KEY OPTIONS & ACCESSORIES

Rear Anti-Tip				Positioning Belt			
<input type="checkbox"/>	115AT1	Rear Anti-Tip	\$121.00	<input type="checkbox"/>	115A115	Belt 2" Aircraft Buckle Long	\$139.00
<input type="checkbox"/>	115AT3	Omit Rear Anti-Tip	NC	<input type="checkbox"/>	115A116	Belt 2" Aircraft Buckle Padded Long	\$204.00
				<input type="checkbox"/>	115A142	Belt 1.5" Auto Buckle	\$97.00

10. COLOR

Standard Frame Colors				Personalized Frame Color			
<input type="checkbox"/>	115C211	Celestial Sapphire	NC	<input type="checkbox"/>	115C112	Desert Camouflage	\$368.00
<input type="checkbox"/>	115C212	Desert Sage	NC	<input type="checkbox"/>	115C43	Mossy Oak Camouflage	\$368.00
<input type="checkbox"/>	115C213	Desert Tan	NC	<input type="checkbox"/>	115C42	Stars & Stripes	\$368.00
<input type="checkbox"/>	115C214	Gloss Black	NC	<input type="checkbox"/>	115C108	Zebra (<i>Choose base color</i>) NA w/ <i>Matte Black and Gloss Black.</i>	\$368.00
<input type="checkbox"/>	115C215	Matte Black	NC				
<input type="checkbox"/>	115C216	Metallic Shadow	NC				
<input type="checkbox"/>	115C217	Midnight Blue	NC				
<input type="checkbox"/>	115C218	Mystic Purple	NC				
<input type="checkbox"/>	115C219	Petal Pink	NC				
<input type="checkbox"/>	115C220	Pink Power	NC				
<input type="checkbox"/>	115C221	Red Pepper	NC				
<input type="checkbox"/>	115C222	Royal Blue	NC				
<input type="checkbox"/>	115C223	Salute Green	NC				
<input type="checkbox"/>	115C224	Scarlet Noir	NC				
<input type="checkbox"/>	115C225	Silver Ore	NC				
<input type="checkbox"/>	115C226	Solar Orange	NC				
<input type="checkbox"/>	115C227	Starlight White	NC				
<input type="checkbox"/>	115C228	Storm Gray	NC				
<input type="checkbox"/>	115C229	Turbo Green	NC				
<input type="checkbox"/>	115C230	Velocity Yellow	NC				



Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8

Customer Service: 800-263-3390 Email: cscanada@sunmed.com www.sunrisemedical.ca © 5/2025

MK-100063-CAN Rev 7