

ZIPPIE Frame Colors

| | | | |
|-----------------|-------------------|--------------------|---------------|
| Gloss Black | Matte Black | Metallic Shadow | Silver Ore |
| Starlight White | Storm Gray | Petal Pink | Desert Sage |
| Desert Tan | Salute Green | Celestial Sapphire | Midnight Blue |
| Royal Blue | Scarlet Noir | Red Pepper | Mystic Purple |
| Pink Power | Turbo Green | Velocity Yellow | Solar Orange |
| Mossy Oak Camo | Stars and Stripes | Desert Camo | Zebra* |

IMPORTANT NOTE:

The color options shown are for reference only. They are not intended as accurate representations.

Refer to the QUICKIE & ZIPPIE Color Brochure or a QUICKIE & ZIPPIE Color Ring for representations of ZIPPIE colors.

Print processes cannot accurately represent paint finishes. Actual colors may vary.

Colors and patterns may vary based on type of chair or frame material.


*Zebra Kolorfusion™ is available on any frame color. Shown on Pink Power.



ZIPPIE Q300 M Mini Available Shroud Colors

| | | |
|------------|--------------|-------------------|
| White | Matte Black | Candy Apple Red |
| Candy Blue | Candy Purple | Green Apple |
| Pearl Pink | Silver | Stars and Stripes |

ZIPPIE Q300 M Mini Available Accent Colors



A comprehensive guide
for selecting the right
ZIPPIE® pediatric wheelchair.

ZIPPIE® Pediatric
WHEELCHAIR
REFERENCE GUIDE

| | ZIPPIE X'CAPE® | ZIPPIE® GS / GS SE | ZIPPIE ZONE® | ZIPPIE SPHYNX® |
|-------------|----------------|--------------------|--------------|------------------------|
| FRAME TYPE: | FOLDING | FOLDING or RIGID | RIGID | FOLDING STATIC TILT |



| KEY FEATURES: | <ul style="list-style-type: none">• XLOCK® locking crossbrace for rigid performance and one-step push-to-lock and pull-to-fold features• Z-finity® System for adjusting lower extremity length and angles• VersaRail™ seat rail for shared placement of components | <ul style="list-style-type: none">• Built-in growth system• Available in reverse configuration• Proven design and the highest quality construction for trouble-free durability | <ul style="list-style-type: none">• Ultra lightweight• Growth capabilities• Easy adjustments• Cool style | <ul style="list-style-type: none">• Variable tilt can be set at 10°, 20°, 30°• Quick-adjust recline back with 85° to 100° range• Ultra compact folded package for easy transport |
|---|--|--|---|--|
| HCPCS CODE | E1238 | Rigid: E1237 Folding: E1238 | E1235, E1237 | E1232, E1234, E1161 |
| Transit Option | Yes ¹ | Yes | Yes | Yes ¹ |
| Tilt Range | _____ | _____ | _____ | 10°, 20°, 30° |
| Chair Weight* | 19 lbs. w/o footrests | Folding: 25 lbs. w/o footrests Rigid: 29 lbs. w/o footrest | 14 lbs. | 28 lbs. w/o footrests |
| Weight Capacity ^{w/o seating} | 165 lbs. | 165 lbs. | 165 lbs. | 250 lbs. |
| Overall Width* | 15" to 26.5" | 17" to 27" | 15" to 25" | 21" to 26" |
| Seat Width | 8" to 16" (grows to 18") | 10" to 18" | 8" to 16" | 10" to 18" |
| Seat Depth | 8" to 18" (grows to 20") | 10" to 20" | 8" to 18" | 10" to 19" (Includes 5" depth growth standard) |
| Back Angle Adjustment | Fore (-) 8° to Aft (+) 24° | Fore (-) 5° to Aft (+) 30° | Fore (-) 17° to Aft (+) 8° | 85° to 100° |
| Front Seat-To-Floor Height | 13" to 18" | 14" to 20" | 14" to 19" | 16.5" to 21" |
| Rear Seat-To-Floor Height | 11" to 18" | 12.5" to 20" | 13" to 19" | Tilt angle 10° = 16" Tilt angle 20° = 14" Tilt angle 30° = 11.5" |
| Footrest / Frame Angle | 80° Swing-In/Out, 80° Fixed Front, Z-finity System: Adjustable ELR: 0° to 90° | Swing-In/Out: 60°, 70°, 80° HD Lift Off: 70°, 80°, 90° Swing-Away: 90° ELR, ALR: 0° to 90° | 80° | Adjustable, removable: 0° to 120° |
| Back Type | Fixed Height, Height Adjustable, Dynamic | Fixed Height and Height Adjustable | Folding Angle Adjustable Double Lock-Down | Folding Angle Adjustable Double Lock-Down |
| Back Height | Fixed Height: 16", 20", 23" Ht Adjustable: 11" to 24", 12" to 24" Dynamic: 12" to 24" | Fixed Height: 19", 22", 25" Ht Adjustable: Rigid: 16" to 22", 19" to 25" Folding: 17" to 23", 20" to 26" | 6" to 8" through 15" to 18" | Fixed Height: 24" or 26" |
| Center of Gravity Adjust. | 3" | 3" | 2.5" | N/A |
| Rear Wheel Sizes | 12", 16", 20", 22", 24" | 12", 16", 18", 20", 22", 24" | 18", 20", 22", 24" | 12" |
| Caster Sizes | 3", 4", 5", 6" | 3", 4", 5", 6", 8" | 3", 4", 5", 6" | 6" |
| True Fit or One-Time Free Growth Kit | Yes | GS: Yes GS SE: No | No ² | No |
| Warranty ^{Refer to Owner's Manual for details} | Frame: Lifetime Parts: 1 Year | Frame: Lifetime Parts: 1 Year | Frame: Lifetime Parts: 1 year | Frame: Lifetime Parts: 1 year |

| ZIPPIE® TS / TS SE | ZIPPIE IRIS® | ZIPPIE VOYAGE® |
|--------------------|-----------------------------------|--------------------------------|
| FOLDING | TILT-IN-SPACE Folding or Rigid | EARLY INTERVENTION Stroller |



| | | |
|--|--|---|
| <ul style="list-style-type: none">• Tilt-in-space up to 45°• Lightweight welded aluminum frame construction ideal for portability• Designed to accept the widest possible variety of seating solutions | <ul style="list-style-type: none">• Patented rotation-in-space technology ensures a smooth, easy-to-control tilting motion• Lightweight frame• Optional XLOCK® folding crossbrace for convenient transport | <ul style="list-style-type: none">• Tilt-in-space range of 35°, optional 40° of recline with Advanced Seating and 90° of recline with Moderate Seating• Personalized seating with extensive built-in adjustments• Functional accessories for child and caregivers |
| Folding: E1234 Rigid: E1233 | Rigid: E1233 Folding: E1234 | E1232 without seating: E1234 |
| Yes | Yes ¹ | Yes |
| 45° | 40°: -5° to 35°, 0° to 40°, +5° to 45° 55°: -5° to 50°, 0° to 55°, +5° to 60° | 35°: -5° to 30° |
| 29 lbs. w/o footrests | Folding: 22.8 lbs. ** Rigid: 26.6 lbs. ** | 28 lbs. |
| 165 lbs. | 225 lbs. ◇ | 75 lbs. |
| 17" to 26" | 17.5" to 25.5" | 24" |
| 10" to 18" | Rigid: 10" to 18" Folding: 10" to 16" (grows to 18") | Adv.: 6" to 13" Mod.: 9.5" to 12" |
| 13" to 20" | 10" to 18" | Adv.: 5" to 13" Mod.: 8" to 13" |
| Fore (-) 5° to Aft (+) 30° | Fore (-) 5° to Aft (+) 30° | Adv.: 85° to 140° Mod.: 90° to 180° |
| 15.5" to 20" | 13" to 19" w/drop seat: 11" to 17.5" | _____ |
| 14" to 19" | 13" to 19" w/drop seat: 11" to 17.5" | _____ |
| Swing-In/Out: 60°, 70°, 80°, 90° HD Lift Off: 60°, 70°, 80° ELR, ALR: 0° to 90° | Swing-In/Out: 60°, 70°, 80°, 90° HD Lift Off: 70°, 80°, 90° ELR, ALR: 0° to 90 Contracture Z-finity System: Adjustable | 0° to 90° |
| Fixed Height and Height Adjustable | Fixed Height, Height Adjustable, Reclining, Dynamic, MONO Backrest System® | Height Adjustable, Reclining |
| Fixed Height: 19", 22", 25" Ht Adjustable: Rigid: 16" to 22", 19" to 25" Folding: 17" to 23", 20" to 26" | Fixed Height: 19", 22", 25" Ht Adjustable: 16" to 22", 19" to 25" Reclining: 21", 24" MONO: 18.7", 21.7" | Adv.: Short: 7" to 13" Adv.: Tall: 10" to 17" Mod.: 16" to 24" (includes head support) |
| None | 6" | _____ |
| 12", 16", 18", 20", 22", 24" | 12", 16", 20", 22", 24" | 12" |
| 3", 4", 5", 6", 8" | 4", 5", 6", 8" | 7" |
| TS: Yes TS SE: No | IRIS: Yes IRIS SE: No | No |
| Frame: Lifetime Parts: 1 Year | Frame: Lifetime Parts: 1 Year | Frame: Lifetime Parts: 1 Year |

| ZIPPIE Q300® M Mini |
|-----------------------|
| MID-WHEEL DRIVE POWER |



| | |
|----------------------|---|
| KEY FEATURES: | <ul style="list-style-type: none">• Our smallest power wheelchair base for superior maneuverability• Optional 12" power seat elevate, power tilt, and manual recline• Pediatric-specific seating with quick seat width and depth growth adjustments |
| HCPCS CODE | K0848 K0856 |
| Transit Tested | Yes ¹ |
| Drive Wheel | Mid |
| Max Speed | 6.0 mph |
| User Weight Capacity | 165 lbs. |
| Electronics | PG Drives VR2 or R-net |
| Suspension | 6-wheel |
| Turning Radius | 17" |
| Overall Base Width | 20.5" |
| Overall Base Length | 33" w/o footrests |
| Drive Wheels | 12" or 14" |
| Casters | 7" |
| Seat Width | 10" to 16" |
| Seat Depth | 10" to 19" |
| Seat-to-Floor Height | 17" to 19" |
| Tilt | Power |
| Recline | Manual |
| Seat Elevate | 12" |
| Batteries | Group 22 |

* Actual measurements may vary according to the chair configuration and accessories.
** Without footrests and rear wheels

1. Conforms with ANSI/RESNA WC-19 standards. Please refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at www.sunrisemedical.com.

2. Includes 3" of built-in frame depth growth and 2" of built-in seat sling depth growth.

◇ WARNING: Addition of the Z-Finity Footrest System to this device limits the maximum allowable user weight to 165 lbs. (75 kg)