

Customized control for every unique journey

No two wheelchair users are alike, and their controls shouldn't be either. That's why SWITCH-IT offers a wide range of alternative drive controls, switches, and buttons—empowering users to navigate, adjust seating, and manage daily functions their way.

Customer				Ship To			
Account				Name			
Account #			Date	Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City			Province
Buyer				PC	Phone		
Mark For				Email / Fax			

* Pricing provided by quote only. Please send all quote requests to cscanada@sunmed.com.

Joysticks

Switch-It MicroGuide Proportional Light Touch Joystick

1 Switch-It MicroGuide Proportional Joystick

<input type="checkbox"/> MG	MicroGuide Joystick Only, DB9 Connector (displacement-style light-touch joystick) MicroGuide requires roughly 40g of activation force for joystick displacement
-----------------------------	--

2 Switch-It MicroGuide Joystick Mounting Packages (select joystick above)

<input type="checkbox"/> MA-SA-MS	Armrest Mount, MicroSeries Connector
<input type="checkbox"/> BT-SA	BulletTray, Armrest Mount
<input type="checkbox"/> SD-SA	SpaceDisc, Armrest Mount
<input type="checkbox"/> MA-MID-BP-MS	Midline Mount (BodyPoint) MicroSeries Connector
<input type="checkbox"/> BT-MID-BP	BulletTray Hand Tray, Midline Mount (BodyPoint)
<input type="checkbox"/> SD-MID-BP	SpaceDisk Hand Tray, Midline Mount (BodyPoint)
<input type="checkbox"/> MA-MID-GAT-MS	Midline Mount (Gatlin), MicroSeries Connector
<input type="checkbox"/> BT-MID-GAT	BulletTray Hand Tray, Midline Mount (Gatlin)
<input type="checkbox"/> SD-MID-GAT	SpaceDisk Hand Tray, Midline Mount (Gatlin)
<input type="checkbox"/> MA-BIB-TH-MS	Bib Mount (Therafin), MicroSeries Connector
<input type="checkbox"/> MA-SVA-CHIN-MS-WH	S/A Chin Mount, MicroSeries Clamp, Whitmyer Connector

Switch-It VersaGuide and VersaGuide EZ Proportional Joysticks

1 Switch-It VersaGuide, VersaGuide EZ, and Flex Switches

<input type="checkbox"/> VG	VersaGuide Joystick Only, DB9 Connector (versatile proportional compact joystick) VersaGuide requires roughly 240g of activation force for joystick displacement
<input type="checkbox"/> VG-EZ	VersaGuide EZ Joystick Only, DB9 Connector (medium-force proportional joystick) VersaGuide requires roughly 120g of activation force for joystick displacement
<input type="checkbox"/> VG-FLEX	Flex Switches with 6" Gooseneck Only (2) Recommended to be ordered in combination with VG or VG-EZ

2 Switch-It VersaGuide Joystick Mounting Packages (select joystick above)

<input type="checkbox"/> MA-SA-VG	Armrest Mount, VersaGuide Connector
<input type="checkbox"/> MA-MID-BP-VG	Midline Mount (BodyPoint), VersaSeries Connector
<input type="checkbox"/> MA-MID-GAT-VG	Midline Mount (Gatlin), VersaSeries Connector
<input type="checkbox"/> MA-BIB-TH-VG	Bib Mount (Therafin), VersaSeries Connector
<input type="checkbox"/> MA-SWA-CHIN-VG-WH	S/A Chin Mount, VersaSeries Clamp, Whitmyer Connector

Switch-It MicroGuide, VersaGuide, and VersaGuide EZ Mode Switch Packages

<input type="checkbox"/>	ML-SWA-WH	Micro Light Switch, 12" Rigid HDR S/A Mount, Whitmyer Connector	
<input type="checkbox"/>	ML-FXE-SWA-WH	Micro Light Switch, 12" Flexible HDR S/A Mount, Whitmyer Connector	
<input type="checkbox"/>	PIKO-SWA-WH	1" Single Switch, 12" Rigid HDR S/A Mount, Whitmyer Connector	
<input type="checkbox"/>	PIKO-FXE-SWA-WH	1" Single Switch, 12" Flexible HDR S/A Mount, Whitmyer Connector	
<input type="checkbox"/>	ES-SWA-WH	Egg Switch, 12" Rigid S/A Mount, Whitmyer Connector	

See Switches section for additional unmounted switch selections.

Head Controls

Switch-It DualPro Proportional Head Array System

Switch-It DualPro Articulating Facial Pads	Switch-It DualPro Round Pads
Long (approximately 6") Lateral Pads	Round (2" diameter) Spot Pads
DualPro Proportional Head Array* with Long "Facial Pads"	DualPro Proportional Head Array* with Round "Spot" Pads
*Make specific DualPro and mounting hardware configuration selections below.	*Make specific DualPro and mounting hardware configuration selections below.

1 Switch-It DualPro and Headrest Mounting Hardware Packages

<input type="checkbox"/>	DPA-WH-CON	DPA - Long Pads, Whitmyer Compatible			
	DPA - Long Pads	<input type="checkbox"/>	DPA-WH-AXYS	<input type="checkbox"/>	DPA-WH-LINX
			with Axys Mount		with Linx Mount
				<input type="checkbox"/>	DPA-WH-LINX2
					with Linx2 Mount
<input type="checkbox"/>	DPP-WH-CON	DPP - Spot Pads, Whitmyer Compatible			
	DPP - Spot Pads	<input type="checkbox"/>	DPP-WH-AXYS	<input type="checkbox"/>	DPP-WH-LINX
			with Axys Mount		with Linx Mount
				<input type="checkbox"/>	DPP-WH-LINX2
					with Linx2 Mount

2 Switch-It DualPro External Mode Switch Packages

Mode Switch Mount attached to DualPro Lateral Pad Rod

<input type="checkbox"/>	FSE-DP	Proximity Sensor, 12" Flexible Mount, Infinity Connector			
	MicroLite Switch, Infinity Connector	<input type="checkbox"/>	ML-RGD-DP	12" Rigid Rod	<input type="checkbox"/>
					ML-FXE-DP
					12" Flexible Mount
	1" Single Switch, Infinity Connector	<input type="checkbox"/>	PIKO-RGD-DP	12" Rigid Rod	<input type="checkbox"/>
					PIKO-FXE-DP
					12" Flexible Mount
<input type="checkbox"/>	ES-RGD-DP	Egg Switch, 12" Rigid Rod, Infinity Connector			

Mode Switch Mount on separate Swing Away

<input type="checkbox"/>	ML-RGD-DP-SWA-WH	MicroLite Switch, 12" Rigid Rod, Swing-Away Connector	
<input type="checkbox"/>	PIKO-RGD-DP-SWA-WH	1" Single Switch, 12" Rigid Rod, Swing-Away Connector	
<input type="checkbox"/>	ES-RGD-DP-SWA-WH	Egg Switch, 12" Rigid Rod, Swing-Away Connector	

Vigo Head Drive

<input type="checkbox"/>	256757	Vigo Head Drive	
--------------------------	--------	-----------------	--

Cubo Mounting Packages (optional)

<input type="checkbox"/>	169966	Cubo Mounting for Omni2	\$ 153.00	<input type="checkbox"/>	256762	Cubo Mounting for Omni1	
--------------------------	--------	-------------------------	-----------	--------------------------	--------	-------------------------	--

Head Controls

Switch-It Switched Head Array Systems and Headrest Mounting Hardware Packages

3-Switched Head Array Systems - Short-Straight Wings (-SS)

<input type="checkbox"/>	3-Switch HDR, SS, 2x 1/8" Mono, Whitmyer Connector	<input type="checkbox"/> HDR3-SS-M2-WH-AXYS Axy Mount	<input type="checkbox"/> HDR3-SS-M2-WH-LINX Linx Mount	<input type="checkbox"/> HDR3-SS-M2-WH-LINX2 Linx2 Mount
--------------------------	--	--	---	---

3-Switched Head Array Systems - Short-Curved Wings (-SC)

<input type="checkbox"/>	3-Switch HDR, SC, 2x 1/8" Mono, Whitmyer Connector	<input type="checkbox"/> HDR3-SC-M2-WH-AXYS Axy Mount	<input type="checkbox"/> HDR3-SC-M2-WH-LINX Linx Mount	<input type="checkbox"/> HDR3-SC-M2-WH-LINX2 Linx2 Mount
--------------------------	--	--	---	---

3-Switched Head Array Systems - Long-Straight Wings (-LS)

<input type="checkbox"/>	3-Switch HDR, LS, 2x 1/8" Mono, Whitmyer Connector	<input type="checkbox"/> HDR3-LS-M2-WH-AXYS Axy Mount	<input type="checkbox"/> HDR3-LS-M2-WH-LINX Linx Mount	<input type="checkbox"/> HDR3-LS-M2-WH-LINX2 Linx2 Mount
--------------------------	--	--	---	---

3-Switched Head Array Systems - Long-Curved Wings (-LC)

<input type="checkbox"/>	3-Switch HDR, LC, 2x 1/8" Mono, Whitmyer Connector	<input type="checkbox"/> HDR3-LC-M2-WH-AXYS Axy Mount	<input type="checkbox"/> HDR3-LC-M2-WH-LINX Linx Mount	<input type="checkbox"/> HDR3-LC-M2-WH-LINX2 Linx2 Mount
--------------------------	--	--	---	---

4-Switched Head Array Systems - Long-Straight Wings (-LS) Only

<input type="checkbox"/>	4-Switch HDR, LS, Mode Outer Left, Whitmyer Connector	<input type="checkbox"/> HDR4-LS-M1L-WH-AXYS Axy Mount	<input type="checkbox"/> HDR4-LS-M1L-WH-LINX Linx Mount	<input type="checkbox"/> HDR4-LS-M1L-WH-LINX2 Linx2 Mount
<input type="checkbox"/>	4-Switch HDR, LS, Mode Outer Right, Whitmyer Connector	<input type="checkbox"/> HDR4-LS-M1R-WH-AXYS Axy Mount	<input type="checkbox"/> HDR4-LS-M1R-WH-LINX Linx Mount	<input type="checkbox"/> HDR4-LS-M1R-WH-LINX2 Linx2 Mount

Sip N' Puff Combo Systems

<input type="checkbox"/>	Combo S&P HDR, LC, Whitmyer Connector	<input type="checkbox"/> HDR-S&P-LC-WH-AXYS Axy Mount	<input type="checkbox"/> HDR-S&P-LC-WH-LINX Linx Mount	<input type="checkbox"/> HDR-S&P-LC-WH-LINX2 Linx2 Mount
--------------------------	---------------------------------------	--	---	---

Switch-It Switched Head Array Systems and Headrest Mounting Hardware Packages

Switch-It Switched Head Array Systems Mode Switch Packages

<input type="checkbox"/>	ML-SWA-WH	Micro Light Switch, 12" Rigid HDR S/A Mount, Whitmyer Connector
<input type="checkbox"/>	ML-FXE-SWA-WH	Micro Light Switch, 12" Flexible HDR S/A Mount, Whitmyer Connector
<input type="checkbox"/>	PIKO-SWA-WH	1" Single Switch, 12" Rigid HDR S/A Mount, Whitmyer Connector
<input type="checkbox"/>	PIKO-FXE-SWA-WH	1" Single Switch, 12" Flexible HDR S/A Mount, Whitmyer Connector
<input type="checkbox"/>	ES-SWA-WH	Egg Switch, 12" Rigid HDR S/A Mount, Whitmyer Connector

Switch-It RIM Proportional Head Control

<input type="checkbox"/>	RIM-WOM	RIM Proportional Head Control without Mounting Hardware Make specific RIM and mounting hardware configuration selections below
--------------------------	---------	---

Switch-It RIM Control Mounting Hardware Packages

<input type="checkbox"/>	RIM - Whitmyer Connector, Quick Release	<input type="checkbox"/> RIM-WH-AXYS Axy Mount	<input type="checkbox"/> RIM-WH-LINX Linx Mount	<input type="checkbox"/> RIM-WH-LINX2 Linx2 Mount
--------------------------	---	---	--	--

RIM Proportional Head Control Mode Switch and Mode Switch Mounting Options

<input type="checkbox"/>	ML-RGD-RIM-SWA-WH	Micro Light Switch, 16" Rigid Rod, Whitmyer Swing-Away
<input type="checkbox"/>	PIKO-RGD-RIM-SWA-WH	1" Single Switch, 16" Rigid Rod, Whitmyer Swing-Away
<input type="checkbox"/>	ES-RGD-RIM-SWA-WH	Egg Switch, 16" Rigid Rod, Whitmyer Swing-Away

Other Drive Controls

Switch-It Joystick in a Can - Arm and Foot Control

<input type="checkbox"/>	JSC	Joystick Only	
1 Switch-It Joystick in a Can - Mounting Options			
<input type="checkbox"/>	MA-SA-JSC	Armrest Mount, Joystick in a Can Connector	
<input type="checkbox"/>	MA-MID-JSC	Midline Mount (BodyPoint), Joystick in a Can Connector	
<input type="checkbox"/>	MA-AC-JSC	Arm Control Mounting Hardware, Large	
<input type="checkbox"/>	MA-FC-MED-JSC	Foot Control Mounting Hardware, Medium	
<input type="checkbox"/>	MA-FC-LG-JSC	Foot Control Mounting Hardware, Large	

2 Switch-It Joystick in a Can - Mode Switch Options

See Switches section for unmounted switch selections. Headrest-mounted mode switch options also available

Switch-It Cool Cube - Preset Packages

Cool Cube - Proximity/Proportional Switched System

<input type="checkbox"/>	CCB-4PS	Cool Cube with DB9 Connector for Proximity/Proportional Driving - 4 Switch Includes: 4 Pro Spot Proportional Switches (1 red, 1 blue, 1 green, and 1 yellow). Requires Omni for operation.	
<input type="checkbox"/>	CCB-5PS	Cool Cube with DB9 Connector for Proximity/Proportional Driving - 5 Switch Includes: 5 Pro Spot Proportional Switches (1 red, 1 blue, 1 green, 1 yellow, and 1 orange). Requires Omni for operation.	

If on/off and/or mode functionality is desired to be controlled via the Cool Cube, a single switch or switches must be selected from the list below.

Switch-It Cool Cube - Personalized Packages

Cool Cube - Proximity/Proportional Driving System

<input type="checkbox"/>	CCB	Cool Cube with DB9 Connector for Proximity/Proportional Driving - Personalized Pro Spot Proportional Switches are required if proportional driving is desired and must be selected below. If non-proportional driving is desired, single switches can be selected in-lieu of the Pro Spot switches from the full list below. Requires Omni for operation	
--------------------------	-----	---	--

If on/off and/or mode functionality is desired to be controlled via the Cool Cube, a single switch or switches must be selected from the list below.

Pro Spot - Proportional Switches

<input type="checkbox"/>	PROSP-R	Pro Spot Proportional Switch - Red	QTY	
<input type="checkbox"/>	PROSP-G	Pro Spot Proportional Switch - Green	QTY	
<input type="checkbox"/>	PROSP-B	Pro Spot Proportional Switch - Blue	QTY	
<input type="checkbox"/>	PROSP-Y	Pro Spot Proportional Switch - Yellow	QTY	
<input type="checkbox"/>	PROSP-O	Pro Spot Proportional Switch - Orange	QTY	

Switches

Single Switches (1/8" monoplug outlet unless otherwise stated)

2.5" Single Switch

<input type="checkbox"/>	244505	2.5" Single Switch - Black	QTY		
<input type="checkbox"/>	244025	2.5" Single Switch - Green	QTY		
<input type="checkbox"/>	244504	2.5" Single Switch - Red	QTY		
<input type="checkbox"/>	244506	2.5" Single Switch - Blue	QTY		
<input type="checkbox"/>	244507	2.5" Single Switch - Yellow	QTY		

1" Single Switch Button

<input type="checkbox"/>	244027	1" Single Switch - Black	QTY		
<input type="checkbox"/>	247620	1" Single Switch - Green	QTY		
<input type="checkbox"/>	247619	1" Single Switch - Red	QTY		
<input type="checkbox"/>	247617	1" Single Switch - Blue	QTY		
<input type="checkbox"/>	247621	1" Single Switch - Yellow	QTY		
<input type="checkbox"/>	247618	1" Single Switch - Orange	QTY		

Egg Switch

<input type="checkbox"/>	244026	Egg Switch - Black	QTY		
<input type="checkbox"/>	244557	Egg Switch - Green	QTY		
<input type="checkbox"/>	244558	Egg Switch - Red	QTY		
<input type="checkbox"/>	244559	Egg Switch - Blue	QTY		
<input type="checkbox"/>	244560	Egg Switch - Yellow	QTY		

Spec Switch

<input type="checkbox"/>	244030	Spec Switch - Black	QTY		
<input type="checkbox"/>	244499	Spec Switch - Green	QTY		
<input type="checkbox"/>	244501	Spec Switch - Red	QTY		
<input type="checkbox"/>	244502	Spec Switch - Blue	QTY		
<input type="checkbox"/>	244503	Spec Switch - Yellow	QTY		

Mini Button - Single Switch

<input type="checkbox"/>	MB-1-GRN	Mini Button Switch-GREEN	QTY		
<input type="checkbox"/>	MB-1-RED	Mini Button Switch-RED	QTY		
<input type="checkbox"/>	MB-1-BLU	Mini Button Switch-BLUE	QTY		

Flex Switch

Works with VG, VG-EZ, DP-A, DP-P only.

<input type="checkbox"/>	FS	Flex Switch, Proximity Sensor on 6" gooseneck, stereo plug output	QTY		
<input type="checkbox"/>	FSE	Flex Switch, Proximity Sensor on 12" gooseneck, stereo plug output	QTY		

Twister Switch

<input type="checkbox"/>	6003579	Twister Switch - Black	QTY		
<input type="checkbox"/>	6003580	Twister Switch - Red	QTY		
<input type="checkbox"/>	6003581	Twister Switch - Green	QTY		

Micro Light Switch

<input type="checkbox"/>	244028	Micro Light	QTY		
--------------------------	--------	-------------	-----	--	--

Uni-Wob Switch

<input type="checkbox"/>	LLS	Uni-Wob Long Lever Switch	QTY		
--------------------------	-----	---------------------------	-----	--	--

Switches

Multi-Function Switches

Dual Switch Selections

<input type="checkbox"/>	MB-2-M	Two Mini Button Switches-RED & GREEN, dual 1/8" mono jack output	QTY		
<input type="checkbox"/>	MB-2-S	Two Mini Button Switches-RED & GREEN, single 1/8" stereo jack output	QTY		

Quad Switch Selections

<input type="checkbox"/>	MB-4-IL	Mini 4 Button - Inline Enclosure, DB9 connector	QTY		
--------------------------	---------	---	-----	--	--

P&G Seat Controller

<input type="checkbox"/>	158709	Ctrl+5 with Voice Activation-Button	QTY		
<input type="checkbox"/>	158708	Ctrl+5 with Voice Activation-Toggle	QTY		

Cable Assemblies

Switch Output Cable Assembly

<input type="checkbox"/>	139025	Switch Output Cable Assm - DB9 to 4 Male Mono	QTY		
--------------------------	--------	---	-----	--	--

DB9 Split Cable

<input type="checkbox"/>	244569	DB9 Split Cable - 1 Female X 2 Male	QTY		
--------------------------	--------	-------------------------------------	-----	--	--

Stereo Cable Extension 90°

<input type="checkbox"/>	244361	Stereo Cable Extension 90° Male to Female	QTY		
--------------------------	--------	---	-----	--	--

DB9 Extension Cable

<input type="checkbox"/>	139026	DB9 Extension Cable	QTY		
--------------------------	--------	---------------------	-----	--	--

Mono Cable Male to Female

<input type="checkbox"/>	244022	Mono Cable Male - to Female - 6"	QTY		
--------------------------	--------	----------------------------------	-----	--	--

Stereo Cable Extension 90°

<input type="checkbox"/>	243795	Stereo Cable Extension Male to Female	QTY		
--------------------------	--------	---------------------------------------	-----	--	--

Switch-It Safety Devices

Switch-It Remote Stop Switch

RSS-XLR version requires Omni or Enhanced Display power cycle after remote button is pushed and chair is inhibited.

<input type="checkbox"/>	RSS-XLR	Remote Stop*, battery charger port connector	QTY		
<input type="checkbox"/>	RSS-MR	Remote Stop, ring connectors to battery terminals	QTY		
<input type="checkbox"/>	RSS-OBC	Remote Stop, OBC connector	QTY		

OptiStop

<input type="checkbox"/>	OS	Fiber Optic Stop, DB9 connector	QTY		
<input type="checkbox"/>	OS-SP	Fiber Optic Stop w/ Sip 'N Puff Mounting Hardware, DB9 connector	QTY		

Joysticks

1 Switch-It MicroGuide Proportional Joystick



2 SwitchIt VersaGuide, VersaGuide EZ



3 VersaGuide with Flex Switches



Head Controls

1 Switch-It DualPro Proportional Head Array



2 Proximity Sensor, 12" Flex Mount



3 1" Single Switch, 12" rigid rod



4 Egg Switch, 12" Rigid Rod, Swing Away Connector



5 Vigo Head Drive




6 Cubo, Omni2 Mount



7 Switch-It Switched Head Array



8 Sip 'N Puff Combo Systems



9 Switch-It RIM Proportional Head Control



Other Drive Controls

1 Switch-It Joystick in a Can



2 Foot Control Hardware



3 Arm Control Hardware



Other Drive Controls

4 Switch-It Cool Cube



Switches

1 2.5 Single Switch

2 1" Single Switch

3 Egg Switch

4 Mini Button - Single Switch

5 Flex Switch



6 Twister Switch

7 Micro Light Switch

8 Uni-Wob Switch

8 Mini Two Button Switches



9 Mini Four Button Switches

10 Ctrl+5 with Voice Activation - Buttons

11 Ctrl+5 with Voice Activation - Toggles






Cable Assemblies

1 Switch Output Cable Assembly

2 DB9 Extension Cable

3 DB9 Extension Cable



4 Mono Cable Ext

5 Stereo Cable Extension 90°

6 Stereo Cable Extension

Switch-It Remote Stop Switch

