



ORDER FORM

# Q700M


with SEDEO<sup>®</sup> LITE<sup>™</sup>

EIPW25

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		Province	Postal Code
Mark For		Phone	

## Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



## Notes

## Chair Specifications

Weights Starting at: (with Batteries)	NPC: 294 lbs. SPC: 318 lbs.
Power Seating	50° of CG Tilt 12" seat elevator
Seat Type	SEDEO® LITE™
Seat Size Range	Width 15" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	16.5" to 19"

## Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs. 350 lbs. HD
Base Size	Width 25", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 5 mph HD, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Colour	9 shroud colours and 7 colour accents
Maximum Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC267		<b>Quickie Q700 M NPC No Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,993.00
250PC268		<b>Quickie Q700 M SPC Single Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,993.00
Heavy Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC269		<b>Quickie Q700 M NPC HD No Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$9,969.00
250PC270		<b>Quickie Q700 M SPC HD Single Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$9,969.00
Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW25 (LITE™ Seat)</i>			\$
250PC271		<b>Quickie Q700 M NPD No Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$10,781.00
250PC272		<b>Quickie Q700 M SPD Single Power Option for SEDEO® LITE™ Seat</b> Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$10,781.00

BASE			
OPTION		DESCRIPTION	RETAIL
User Weight			\$
250WG23		Up to 300 lbs.	NC
250WG15		301 to 350 lbs.	NC
Motor Package			\$
250PF98		4-pole 8 kph / 5 mph	NC
250PF99		4-pole 10 kph / 6 mph	NC
250PF100		4-pole 12.5 kph / 8 mph	NC
250PF101		4-pole 12.5 kph / 8 mph	\$1,268.00
Tracking			\$
250MI01		SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,696.00
Drive Wheels			\$
250DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			\$
250CW55		Caster Wheels 6", Solid, Black	NC
Shroud Colour			\$
250C119		Matte Black	NC
250C125		Pearl White	NC
250C165		Candy Apple Red	NC
250C166		Candy Blue	NC
250C167		Candy Purple	NC
250C176		Green Apple	NC
250C178		Pearl Pink	NC
250C182		Silver	NC
Colour Accents			\$
<i>Note: Colour Accents and Shroud Colours may not match.</i>			
250EZ20		Apple Green	NC
250EZ22		Black	NC
250EZ21		Blue	NC
250EZ23		Orange	NC
250EZ24		Red	NC
250EZ19		White	NC
250EZ25		Silver	NC

BASE			
OPTION		DESCRIPTION	RETAIL
<b>ISO Transit Tie Downs</b>			\$
250T050		ISO Transit Package	NC
<b>LED Lights</b>			\$
250LI08		LED Fender Lights	\$646.00
250LI02		Omit LED Lights	NC
<b>Seat Heights</b>			\$

Note: Available seat-to-floor heights are dependent on 2 options: **1- Battery Size** and **2- Powered Seating**.  
Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

	Seat Height	Power Seating	No Power Seating	
250HT103	16.5"	Group 34	Group 34	NC
250HT55	17"	Group 34	GP24 or GP34	NC
250HT104	17.5"	Group 34	GP24 or GP34	NC
250HT56	18"	GP24 or GP34	GP24 or GP34	NC
250HT105	18.5"	GP24 or GP34	GP24 or GP34	NC
250HT57	19"	GP24 or GP34	GP24 or GP34	NC

### Batteries and Chargers

#### Group 24 Battery Options (Recommended for max daily range)

250AG25		GP24 Flat Top Pair w/Strap MK Gel	\$955.00
250AG26		GP24 Flat Top Pair w/Strap AGM	\$915.00

#### Group 34 Battery Options (Required for lower seat-to-floor heights)

250AG27		GP34 Flat Top Pair w/Straps MK Gel	\$955.00
250AG29		GP34 Flat Top Pair w/Straps AGM	\$863.00

#### Omit Batteries

250AG24		Battery Harnesses Only (Omit Batteries)	NC
---------	--	---	----

#### Chargers

250BJ28		10A Dual Mode (Compact)	NC
---------	--	-------------------------	----

**POWERED SEATING**

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

**Power Seating Options** \$**Seat Elevator**

250PI08	C-Me® Seat Elevator, Powered, 300 mm / 12" (Not available with Group 4) C-Me speed up to 3 mph.	\$3,477.00
---------	--	------------

**Seat Tilt**

250PA35	Posterior 50° Tilt	\$8,663.00
250PA36	Posterior 50° HD Tilt	\$8,923.00

**Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)** \$

250FA10	0° (Adds 0" to STFH)	NC
250FA3	3° (Adds 1" to STFH)	NC

**SEDEO® LITE™ Backrest Frame (Angle Adjustable 86° - 106°)** \$**Back Frame Width**

250BF21	Width Adjustable, 385 - 485 mm / 15 - 19"	NC
250BF22	Width Adjustable, 460 - 510 mm / 18 - 20"	NC

**Back Cane Height**

250PH30	Height Adjustable, Short w/Trailing Push Handles, 20" to 22"	NC
250PH31	Height Adjustable, Tall w/Trailing Push Handles, 23" to 26"	NC

**Backrest Frame Option** \$

250BF23	Flip Down, Quick Release	\$179.00
---------	--------------------------	----------

**Back Cushion Type** \$

250BT148	SEDEO® LITE™ Sling Back with Tension Adjustable Straps (Choose Seating in Section A)	NC
250BT135	JAY® J3® (Choose Seating in Section B)	NC
250BT133	OMIT (For personalized mounting. No headrest mount)	NC

**A - SEDEO® LITE™ Sling Back Dimensions** \$**Installed Back Height - (From Pan to Top of Back)**

250JR71	18"	NC
250JR72	19"	NC
250JR73	20"	NC
250JR74	21"	NC
250JR75	22"	NC
250JR76	23"	NC

**POWERED SEATING**

OPTION	DESCRIPTION	RETAIL
<b>B - JAY® J3® Back</b>		
<b>B - JAY® J3® Back - Contour</b>		<b>\$</b>
250JP01	J3® PA Back for 16"-20" Wide	\$912.00
250JP11	J3® PL Back for 16"-20" Wide	\$986.00
250JP13	J3® PD Back for 16"-20" Wide <i>(18" and 21" shell height only)</i>	\$1,086.00
<b>B - JAY® J3® Back Mounting Hardware</b>		<b>\$</b>
250JP20	TH	NC
<b>B - JAY® J3® Back Width</b>		<b>\$</b>
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>		
250JQ02	16"	NC
250JQ03	18"	NC
250JQ04	20"	NC
<b>B - JAY® J3® Back Shell Options and Height Ranges</b>		<b>\$</b>
<b>JAY® J3® Shell Height</b>		
250JR14	16.5" tall shell <i>(Upper Thoracic Short)</i>	NC
250JR15	18" tall shell <i>(Upper Thoracic Med)</i>	NC
250JR16	19.5" tall shell <i>(Upper Thoracic Tall)</i>	NC
250JR17	21" tall shell <i>(Shoulder Height)</i>	NC
250JR27	22.5" tall shell <i>(Shoulder Height)</i>	NC
250JR28	24" tall shell <i>(Shoulder Height)</i>	NC
<b>JAY® J3® Installed Back Height (From Pan to Top of Back)</b>		
250JR29	21"	NC
250JR19	22"	NC
250JR20	23"	NC
250JR21	24"	NC
250JR22	25"	NC
250JR23	26"	NC
250JR24	27"	NC
250JR25	28"	NC
250JR30	29"	NC
250JR31	30"	NC

**POWERED SEATING**

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

**B - JAY® Anterior Thoracic Support**

\$

*Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)*

250K10	Medium (ATS8302)	\$322.00
250K11	Large (ATS8303)	\$322.00

**B - JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery, & Mount)**

\$

**JAY® Lateral Thoracic Support Pads**

250JW35	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Right	\$75.00
250JW38	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Left	\$75.00
250JW36	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$75.00
250JW39	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$75.00
250JW37	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$75.00
250JW40	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$75.00

**JAY® Lateral Thoracic Support Pad Upholstery**

250NU01	Syncra™ Premium Leatherette	NC
250NU02	Syncra™ Stretch	NC

**JAY® Lateral Thoracic Swing-Away Mount**

250JW41	Lateral Thoracic Swing-Away Bracket - Right	\$334.00
250JW42	Lateral Thoracic Swing-Away Bracket - Left	\$334.00

**Headrests (Requires Pad, Upholstery, & Mount)**

\$

**Headrest Pad**

250HE36	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$227.00
250HE38	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$377.00
250HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$377.00
250HE44	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$377.00
250HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium ( <i>Adjustable from 7 to 16" Wide x 4" Tall</i> )	\$445.00
250HE48	OMIT Headrest Pad	NC

**Headrest Pad Upholstery**

250NU10	SEDEO® Upholstery	NC
250NU11	Whitmyer® Lycra	NC

**Headrest Mount**

250MH32	SEDEO® Multi-Position Headrest Mount with Quick Release ( <i>Only available with SEDEO® Pad</i> )	\$557.00
250MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$247.00
250MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$298.00
250MH35	Whitmyer® AXYS Headrest Mount with Quick Release	\$247.00
250MH17	OMIT Headrest Mount	NC

**POWERED SEATING**

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

**Seat Cushion Type** \$

250ST80	JAY® (See Section C)	NC
250ST75	Omit	NC

**B - Seat Size** \$**Seat Pan Width**

250R78	15"	NC
250R53	16"	NC
250R54	17"	NC
250R55	18"	NC
250R56	19"	NC
250R57	20"	NC

**Seat Pan Depth (Front of seat pan to front of back cushion – effective depth)**

250RF51	16"	NC
250RF45	17"	NC
250RF46	18"	NC
250RF47	19"	NC
250RF48	20"	NC
250RF49	21"	NC
250RF50	22"	NC

**C - JAY® Cushions****C - JAY® Seat Cushion** \$

250CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$657.00
250CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$764.00
250CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$1,045.00
250CU104	JAY® Fusion® with Cryo® Technology (Sizes 16" - 21" Wide)	\$1,152.00

## POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
<b>Knee and Thigh Swing-Away Support Pads</b>		<b>\$</b>
<b>Lateral Knee Support Pads and Hardware (must pick pad and mount)</b>		
250XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
250XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
250XX08	Lateral Swing-Away Bracket Knee - Right	\$325.00
250XX09	Lateral Swing-Away Bracket Knee - Left	\$325.00
<b>Lateral Thigh Support Pads and Hardware (must pick pad and mount)</b>		
250KF22	Lateral Thigh Support Pad Short, Right (5"x5")	\$146.00
250KF23	Lateral Thigh Support Pad Short, Left (5"x5")	\$146.00
250KF24	Lateral Thigh Support Pad Long, Right (5"x7")	\$146.00
250KF25	Lateral Thigh Support Pad Long, Left (5"x7")	\$146.00
250KF26	Lateral Thigh Removable Mount - Right	\$325.00
250KF27	Lateral Thigh Removable Mount - Left	\$325.00
<b>Armrest Structure</b>		<b>\$</b>
250EG35	Reclining, Flip Back (Height Adjustable 8" to 13")	NC
250EG36	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 13" x .625")	\$351.00
250EG34	Single Post, SEDEO LITE, Removable (N/A with thigh support), (Height Adjustable 10" to 15")	NC
250EG50	Dual Post HA Removable	NC
<b>Arm Pads</b>		<b>\$</b>
<b>SEDEO® Armrest Pad</b>		
250AB62	300 mm / 12", Integral Foam, Right	\$113.00
250AE62	300 mm / 12", Integral Foam, Left	\$113.00
250AB63	400 mm / 16", Integral Foam, Right	\$113.00
250AE63	400 mm / 16", Integral Foam, Left	\$113.00
250AB61	SEDEO LITE Pad, Right	NC
250AE61	SEDEO LITE Pad, Left	NC
<b>Comfort Arm Pad (Foam)</b>		
250AB46	2" Wide Right	\$113.00
250AE46	2" Wide Left	\$113.00
<b>Comfort Arm Pad (Gel)</b>		
250AB48	2" Wide Right	\$213.00
250AE48	2" Wide Left	\$213.00
<b>Dual Post HA Armrest Pads</b>		
250AB74	Armpad 10" classic dsk foam R	NC
250AB75	Armpad 14" Classic FII Foam R	NC
250AE74	Armpad 10" Classic Dsk Foam L	NC
250AE75	Armpad 14" Classic FII Foam L	NC
250AB76	Armpad 10" WtrfII Dsk Foam R	NC
250AB77	Armpad 14" WtrfII FII Foam R	NC
250AE76	Armpad 10" WtrfII Dsk Foam L	NC
250AE77	Armpad 14" WtrfII FII Foam L	NC

**Ottobock® Arm Troughs**

250AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
250AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
250AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
250AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
2500L01	Ottobock Hand Pad, Med (7"), Left	\$70.00
2500R01	Ottobock Hand Pad, Med (7"), Right	\$70.00
2500L02	Ottobock Hand Pad, Large (9"), Left	\$70.00
2500R02	Ottobock Hand Pad, Large (9"), Right	\$70.00
250AE78	Ottobock One-Piece Arm Trough (18"), Left	\$198.00
250AB78	Ottobock One-Piece Arm Trough (18"), Right	\$198.00

**POSITIONING SUPPORTS & ACCESSORIES**

OPTION	DESCRIPTION	RETAIL
<b>Elbow Support</b>		<b>\$</b>
2503A02	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$228.00
2503A03	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$228.00
2503A06	Elbow Support Pad, RH 3" x 3"	\$167.00
2503A07	Elbow Support Pad, LH 3" x 3"	\$167.00
2503A04	Swing-Away Elbow Support Hardware, RH	\$279.00
2503A05	Swing-Away Elbow Support Hardware, LH	\$279.00

**Fixed Centermount Legrest**

<b>Fixed Centermount Legrest</b>		<b>\$</b>
250H183	SEDEO® Fixed Centermount Legrest <i>(3" to 17.5" Leg Length)</i>	\$351.00
250H184	Fixed Centermount Legrest, Short <i>(Includes Calf Pads &amp; Single Footplate) (11.5" to 14.5" Leg Length)</i>	NC
250H185	Fixed Centermount Legrest, Long <i>(Includes Calf Pads &amp; Single Footplate) (13.5" to 18.5" Leg Length)</i>	NC
250H136	Angle Adjustable Centermount Legrest <i>(11.5" to 18.5" Leg Length)</i>	\$605.00

**Power Centermount Legrest**

<b>Power Centermount Legrest</b>		<b>\$</b>
250H175	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Mechanical Articulation) (5" to 17.5" Leg Length)</i>	\$2,487.00
250H109	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i>	\$3,435.00

**Power Centermount Legrest Accessories**

<b>Power Centermount Legrest Accessories</b>		<b>\$</b>
<b>SEDEO Centermount Footplates</b>		
250F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$284.00
250F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$284.00

**Footplate Accessories (Available with Power Centermounts)**

250DI68	Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i>	\$315.00
250DI69	Footplate Connector Link <i>(Attaches footplates together)</i>	\$96.00
250DI57	Gel Footpads <i>(Pair, LH and RH)</i>	\$401.00

**Angle Adjustable Centermount Footplates**

250F45		Single / Dual Footplate	\$284.00
250F44		Single Platform Footplate	\$284.00

**SEDEO Centermount Calf Pads**

250DI34		Calf Pads, Pair, Standard Integral Foam	NC
250DI35		Calf Pads, Pair, Padded	\$408.00

**Angle Adjustable Centermount Calf Pads**

250DI49		Calf Pads, Pair, Standard Integral Foam	NC
250DI50		Extra Large Single Calf Pad	\$197.00

**POSITIONING SUPPORTS & ACCESSORIES**

OPTION		DESCRIPTION	RETAIL
--------	--	-------------	--------

**SEDEO® PRO™ Swing-Away Legrest (Pair)**

			\$
250H177		Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length)	NC
250H178		Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$605.00
250H179		Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$3,435.00
250DI76		Omit Legrests	NC

**Swing-Away Legrest Footplates**

250F369		Angle Adjustable Flip-Up Footplates	\$284.00
---------	--	-------------------------------------	----------

**Swing-Away Legrest Calf Supports**

250DI42		Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$512.00
---------	--	---	----------

**Swing-Away Legrest Options**

250DI44		Heel Loops	NC
250KG4		Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$282.00

**SEDEO® Lower Extremity Swing-Away Support**

250VD01		Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$512.00
250VD02		Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$512.00

**Pelvic Positioning Belts**

250A837		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
250A839		Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
250A840		Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$103.00
250A845		Bodypoint – Evoflex Pelvic Stabilizer Belt	\$216.00
250A846		Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	\$271.00

**Accessories**

			\$
--	--	--	----

**Ventilator Option**

250VT12		Vent Tray Fixed (Pulmonetic LTV Models)	\$609.00
---------	--	---	----------

**O2 Tank Option**

250A718		O2 Tank Holder Mounted to Backrest (D and E Styles)	\$391.00
---------	--	---	----------

**Wheelchair Tray Table**

2507A09		Tray, Swing-Away	\$381.00
---------	--	------------------	----------

**Backpack**

250A1		Black Utility	\$112.00
-------	--	---------------	----------

**Medial Thigh Support**

2508A80	Medial Thigh Support	\$271.00
2508A81	Pommel Swing-Away Hardware	\$390.00

**Misc. Accessories**

2508A59	Transfer Bar <i>(Pair)</i>	\$493.00
2508A52	Ram Mount Phone	\$77.00
2508A53	Ram Mount Tablet	\$308.00
2508A54	Ram Mount Cup Holder	\$77.00
2508A58	Rear View Mirror	\$77.00
250A611	Crutch Holder	\$174.00
2508A73	Bag for Armrest	\$77.00

**JOYSTICK / ELECTRONICS**

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

**R-net Joystick *(Remote)***

\$

250J121	Standard LED Joystick <i>(w/o Screen Display)</i>	NC
250J116	Colour Joystick <i>(2.3" Colour LCD Screen Display, CJSM1)</i>	NC
250J118	High Resolution Colour Joystick,(CJSM2+BT) <i>3.5" Colour LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,162.00
250J119	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC

**Joystick Mounting**

\$

250JM34	Fixed on Armrest	NC
250JM35	Parallel Swing-Away	\$447.00
250JM36	Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>	\$447.00

**Joystick Mount (Location)**

\$

250HM16	Right Hand	NC
250HM17	Left Hand	NC

**Display Language Selection R-net**

\$

250L1	English	NC
250L3	French	NC
250L7	Spanish	NC

**Joystick Handle (Optional)**

\$

250JH13	Ball Handle	\$110.00
250JH15	T-Handle	\$110.00
250JH24	Bodypoint - Rubber Dome Joystick Handle	\$216.00
250JH25	Bodypoint - U-Shaped Joystick Handle <i>(4 inch)</i>	\$216.00
250JH26	Bodypoint - Grooved Mushroom Joystick Handle	\$216.00

## Electronic Programmers \$

250PR13	R-net PC Wireless Programmer	\$660.00
---------	------------------------------	----------

## R-net Power Modules \$

250ER11	Power Module, R-net, 120A, Non-Expandable	NC
250ER10	Power Module, R-net, 120A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$102.00
<input type="checkbox"/> 250EX09	Expandable 120A Controller - R-net	NC
<input type="checkbox"/> 250EX10	Harness for 120A Expandable Controller	NC
250ER18	Power Module, R-net, 90A, Non-Expandable	NC
250ER17	Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$102.00
<input type="checkbox"/> 250EX15	Expandable 90A Controller - R-net	NC
<input type="checkbox"/> 250EX16	Harness for 90A Expandable Controller	NC

## JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

## Seating Control \$

### Powered Seating Control

250E040	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
250E037	Thru-Drive Control for 1 Actuator	\$1,443.00
250E038	Thru-Drive Control for 2+ Actuators	\$2,903.00

### Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*

250E033	Dual Toggle for 1 Actuator	NC
250E034	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
250E041	Ctrl+5 Switch Box for 2 Actuators or More	NC

### Additional Power Seating Control *(Available with Thru-Drive only)*

250QG11	Ctrl+5 Switch Box with Voice Activation	\$1,078.00
---------	---	------------

### Dual Toggle Mounting Location

250TN1	Toggle Mount, Right	NC
250TN2	Toggle Mount, Left	NC

### Ctrl+5 Mounting Location *(Will swing away if mounted in same location as joystick)*

250EY19	Right Armrest	NC
250EY20	Left Armrest	NC
250EY21	Back Cane, Right	NC
250EY24	Back Cane, Left	NC

### Ctrl+5 Switch Type

250EY22	Push Buttons	NC
250EY23	Toggle Switches	\$356.00

## Attendant Controls \$

2500A05	Dual Control Module for Attendant	\$1,225.00
2500A06	R-net Colour LCD (CJSM1), Attendant Control Only	\$1,225.00

Attendant Mounting Location			€\$
250DN5		Right Back Cane	NC
250DN6		Left Back Cane	NC

### R-net Programming €\$

#### Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

250PM12		Include Assignable Buttons	NC
---------	--	----------------------------	----

#### Power Centermount Footplate Extension Option

250IE03		Include 7.5" Powered Footplate Extension for Transfer Assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.	\$592.00
---------	--	--	----------



Q700 M with SEDEO® LITE™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Signature

