

ORDER FORM




Q300M with SEDEO[®] LITE[™]

EIPW40

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		Province	Postal Code
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

Weights starting at: (with Batteries)	NPC: 250 lbs SPC: 274 lbs
Power Seating	50° or 30° of CG Tilt / 12" seat elevator
Seat Type	SEDEO® LITE™
Seat Size Range	Width 15" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 20.5", Length 33"
Minimum Turning Radius	17" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A / LiNX 120A
Driving Range	Up to 17 miles
Drive Motors	4-pole 6 mph / 10 kph
Drive Wheels	12" w/flat-free inserts (black) 14" w/flat-free inserts (black)
Casters	7" front / 7" rear (black)
Charger	8A Dual Mode (Compact)
Battery Size	50Ah Group 22
Color	9 shroud colors and 7 color accents
Ground Clearance	3"

BASE

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Standard Duty Mid Wheel Drive Base *End Item Number EIPW40 (LITE™ Seat)* \$

400PC294	Quickie Q300 M NPC No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,072.00
400PC295	Quickie Q300 M SPC Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q300 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,072.00

User Weight \$

400WG23	Up to 300 lbs.	NC
---------	----------------	----

Motor Package \$

400PF99	4-pole 10 kph / 6 mph	NC
---------	-----------------------	----

Tracking \$

400MI01	R-net SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,696.00
400MI04	LiNX Traction Control System for LiNX <i>(GYRO Technology)</i>	\$925.00

Drive Wheels \$

400DW13	Drive Wheel, 12", Solid, Knobby, Black	NC
400DW58	Drive Wheel, 14", Solid, Knobby, Black	NC

Caster Wheels \$

400CW59	Caster Wheels 7", Solid, Black	NC
---------	--------------------------------	----

Shroud Color \$

400C164	Black	NC
400C165	Candy Apple Red	NC
400C166	Candy Blue	NC
400C167	Candy Purple	NC
400C176	Green Apple	NC
400C178	Pearl Pink	NC
400C182	Silver	NC
400C125	White	NC

Color Accents \$

Note: Color Accents and Shroud Colors may not match.

400EZ20	Apple Green	NC
400EZ22	Black	NC
400EZ21	Blue	NC
400EZ23	Orange	NC
400EZ24	Red	NC
400EZ19	White	NC
400EZ25	Silver	NC

BASE			
OPTION		DESCRIPTION	RETAIL
ISO Transit Tie Downs			\$
400T050		ISO Transit Package	NC
LED Lights			\$
400LI08		LED Fender Lights	\$646.00
400LI02		Omit LED Lights	NC
Seat Heights			\$
Seat Height <i>Measured from front of seat pan to floor with 0° frame angle. Not all seat heights are available with all seat depths and power seating. Please see online configurator for full availability.</i>			
400HT103		16.5"	NC
400HT55		17"	NC
400HT104		17.5"	NC
400HT56		18"	NC
400HT105		18.5"	NC
400HT57		19"	NC
Batteries and Chargers			\$
Group 22 Battery			
400AG35		GP22 / 50Ah Flat Top Pair	\$678.00
Omit Batteries			
400AG24		Battery Harnesses Only <i>(Omit Batteries)</i>	NC
Chargers			
400BJ30		10A Dual Mode (Compact)	\$162.00
400BJ19		8A Dual Mode (Compact)	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Power Seating Options \$**Seat Elevator**

400PI07	Seat Elevator, Powered, 300 mm / 12"	\$3,477.00
400PI08	C-Me® Seat Elevator, Powered, 300 mm / 12" <i>(C-Me speed up to 3 mph) (N/A with VR-2)</i>	\$3,477.00

Seat Tilt

400PA48	Posterior 50° Tilt	\$6,496.00
---------	--------------------	------------

Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height) \$

400FA10	0° <i>(Adds 0" to STFH)</i>	NC
400FA3	3° <i>(Adds 1" to STFH)</i>	NC

SEDEO® LITE™ Backrest Frame (Angle Adjustable 86° to 106°) \$**Back Frame Width**

400BF21	Width Adjustable, 385 to 485 mm / 15" to 19"	NC
400BF22	Width Adjustable, 460 to 510 mm / 18" to 20"	NC

Back Cane Height

400B263	Standard, 18" to 20"	NC
400B264	Tall, 21" to 23"	NC
400PH30	Height Adjustable, Short w/Trailing Push Handles, 20" to 22"	NC
400PH31	Height Adjustable, Tall w/Trailing Push Handles, 23" to 26"	NC

Backrest Frame Option \$

400BF23	Flip Down, Quick Release	\$179.00
---------	--------------------------	----------

Back Cushion Type \$

400BT148	SEDEO® LITE™ Sling Back with Tension Adjustable Straps <i>(Choose Seating in Section A)</i>	NC
400BT135	JAY® J3® <i>(Choose Seating in Section B)</i>	NC
400BT133	OMIT <i>(For personalized Mounting. No Headrest Mount)</i>	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
A - SEDEO® LITE™		
A - SEDEO® LITE™ Sling Back Dimensions		
Installed Back Height - (From Pan to Top of Back)		
400JR71	18"	NC
400JR72	19"	NC
400JR73	20"	NC
400JR74	21"	NC
400JR75	22"	NC
400JR76	23"	NC
B - JAY® J3® Back		
B - JAY® J3® Back - Contour		
400JP01	J3® PA Back for 16"-20" Wide	\$912.00
400JP11	J3® PL Back for 16"-20" Wide	\$986.00
400JP13	J3® PD Back for 16"-20" Wide <i>(18" and 21" Shell Height Only)</i>	\$1,068.00
B - JAY® J3® Back Mounting Hardware		
400JP20	TH	NC
B - JAY® J3® Back Width		
400JQ02	16"	NC
400JQ03	18"	NC
400JQ04	20"	NC
B - JAY® J3® Back Shell Options and Height Ranges		
JAY® J3® Shell Height		
400JR14	16.5" tall shell <i>(Upper Thoracic Short)</i>	NC
400JR15	18" tall shell <i>(Upper Thoracic Med)</i>	NC
400JR16	19.5" tall shell <i>(Upper Thoracic Tall)</i>	NC
400JR17	21" tall shell <i>(Shoulder Height)</i>	NC
400JR27	22.5" tall shell <i>(Shoulder Height)</i>	NC
400JR28	24" tall shell <i>(Shoulder Height)</i>	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
JAY® J3® Installed Back Height (From Pan to Top of Back)		
400JR29	21"	NC
400JR19	22"	NC
400JR20	23"	NC
400JR21	24"	NC
400JR22	25"	NC
400JR23	26"	NC
400JR24	27"	NC
400JR25	28"	NC
400JR30	29"	NC
400JR31	30"	NC

B - JAY® Anterior Thoracic Support

\$

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY® backrest)

400K10	Medium (ATS8302)	\$325.00
400K11	Large (ATS8303)	\$325.00

B - JAY® Lateral Thoracic Swing-Away Supports (Requires pad, upholstery, and mount)

\$

Lateral Thoracic Support Pads

400JW35	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Right	\$73.00
400JW38	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Left	\$73.00
400JW36	Lateral Thoracic (LxH) 100 mm x 152 mm / 4"x 6" Pad - Right	\$73.00
400JW39	Lateral Thoracic (LxH) 100 mm x 152 mm / 4"x 6" Pad - Left	\$73.00
400JW37	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$73.00
400JW40	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$73.00

Lateral Thoracic Support Pad Upholstery

400NU01	Premium Leatherette	NC
400NU02	Stretch	NC

Lateral Thoracic Swing-Away Mount

400JW41	Lateral Thoracic Swing-Away Bracket - Right	\$325.00
400JW42	Lateral Thoracic Swing-Away Bracket - Left	\$325.00

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Headrests (Requires Pad, Upholstery, & Mount)

\$

Headrest Pad

400HE36	SEDEO® Pad Small with Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with quick-release and SEDEO® upholstery</i>	\$227.00
400HE38	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$377.00
400HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$377.00
400HE44	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$377.00
400HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (<i>Adjustable 7" to 16" Wide x 4" Tall</i>)	\$445.00
400HE48	OMIT Headrest Pad	NC

Headrest Pad Upholstery

400NU10	SEDEO® Upholstery	NC
400NU11	Whitmyer® Lycra	NC

Headrest Mount

400MH32	SEDEO® Multi-Position Headrest Mount with Quick Release (<i>Only available with SEDEO® Pad</i>)	\$557.00
400MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$247.00
400MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$298.00
400MH35	Whitmyer® AXYS Headrest Mount with Quick Release	\$247.00
400MH17	OMIT Headrest Mount	NC

Seat Cushion Type

\$

400ST80	JAY® (See Section C)	NC
400ST75	Omit	NC

Seat Size

\$

Seat Pan Width

400R78	15"	NC
400R53	16"	NC
400R54	17"	NC
400R55	18"	NC
400R56	19"	NC
400R57	20"	NC

Seat Pan Depth (Front of seat pan to front of back cushion – effective depth)

400RF51	16"	NC
400RF45	17"	NC
400RF46	18"	NC
400RF47	19"	NC
400RF48	20"	NC
400RF49	21"	NC
400RF50	22"	NC

POWERED SEATING

OPTION	DESCRIPTION	RETAIL
C - JAY® Seat Cushion		
C - JAY® Seat Cushion		\$
400CU59	JAY® Ion® (Sizes 16" to 20" Wide)	\$505.00
400CU61	JAY® Union® (Sizes 16" to 20" Wide)	\$764.00
400CU64	JAY® Fusion® (Sizes 16" to 20" Wide)	\$1,045.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Knee and Thigh Swing-Away Support Pads		\$
Lateral Knee Support Pads and Hardware (must pick pad and mount)		
400XY05	Lateral Knee Support Pad, Right (4" round)	\$247.00
400XY06	Lateral Knee Support Pad, Left (4" round)	\$247.00
400XX08	Lateral Swing-Away Bracket Knee - Right	\$325.00
400XX09	Lateral Swing-Away Bracket Knee - Left	\$325.00
Lateral Thigh Support Pads and Hardware (must pick pad and mount)		
400KF22	Lateral Thigh Support Pad Short, Right (5"x5")	\$146.00
400KF23	Lateral Thigh Support Pad Short, Left (5"x5")	\$146.00
400KF24	Lateral Thigh Support Pad Long, Right (5"x7")	\$146.00
400KF25	Lateral Thigh Support Pad Long, Left (5"x7")	\$146.00
400KF26	Lateral Thigh Removable Mount - Right	\$325.00
400KF27	Lateral Thigh Removable Mount - Left	\$325.00

Armrest Structure

Armrest Structure		\$
400EG41	Reclining, Flip Back, Short (Best used with shorter seat depths) (Height Adjustable 9" to 15")	NC
400EG36	Reclining, Flip Back with 2-point post locking arm. (Height Adjustable 11" to 15" x 0.625")	\$351.00
400EG34	Single Post, SEDEO LITE, Removable (N/A with thigh support) (Height Adjustable 10" to 15")	NC
400EG50	Dual Post, Flip Back, Removable (Height Adjustable)	NC

Arm Pads

Arm Pads		\$
SEDEO® Armrest Pad		
400AB62	300 mm / 12", Integral Foam, Right	\$54.00
400AE62	300 mm / 12", Integral Foam, Left	\$54.00
400AB63	400 mm / 16", Integral Foam, Right	\$54.00
400AE63	400 mm / 16", Integral Foam, Left	\$54.00
400AB61	SEDEO LITE Pad, Right	NC
400AE61	SEDEO LITE Pad, Left	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Dual Post Armrest Pad		
400AB74	250 mm / 10", Classic Desk Foam, Right	NC
400AE74	250 mm / 10", Classic Desk Foam, Left	NC
400AB76	250 mm / 10", Waterfall Desk, Right	NC
400AE76	250 mm / 10", Waterfall Desk, Left	NC
400AB75	350 mm / 14", Classic Full Foam, Right	NC
400AE75	350 mm / 14", Classic Full Foam, Left	NC
400AB77	350 mm / 14", Waterfall Full, Right	NC
400AE77	350 mm / 14", Waterfall Full, Left	NC
Comfort Arm Pad (Foam)		
400AB46	2" Wide Right	\$113.00
400AE46	2" Wide Left	\$113.00
Comfort Arm Pad (Gel)		
400AB48	2" Wide Right	\$213.00
400AE48	2" Wide Left	\$213.00
Ottobock® Arm Troughs		
400AE06	Ottobock Arm Trough, Med (12"), Left	\$139.00
400AB06	Ottobock Arm Trough, Med (12"), Right	\$139.00
400AE07	Ottobock Arm Trough, Large (14"), Left	\$139.00
400AB07	Ottobock Arm Trough, Large (14"), Right	\$139.00
400OL01	Ottobock Hand Pad, Med (7"), Left	\$70.00
400OR01	Ottobock Hand Pad, Med (7"), Right	\$70.00
400OL02	Ottobock Hand Pad, Large (9"), Left	\$70.00
400OR02	Ottobock Hand Pad, Large (9"), Right	\$70.00
400AE78	Ottobock 1 Piece Arm Trough (18"), Left	\$198.00
400AB78	Ottobock 1 Piece Arm Trough (18"), Right	\$198.00
Elbow Support \$		
4003A02	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$203.00
4003A03	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$203.00
4003A06	Elbow Support Pad, RH 3" x 3"	\$167.00
4003A07	Elbow Support Pad, LH 3" x 3"	\$167.00
4003A04	Swing-Away Elbow Support Hardware, RH	\$279.00
4003A05	Swing-Away Elbow Support Hardware, LH	\$279.00
Fixed Centermount Legrest \$		
400H136	Angle Adjustable Centermount Legrest	\$633.00
400H184	Fixed Centermount Legrest, Short <i>(Includes Calf Pads) (11.5" to 14.5" Leg Length)</i>	NC
400H185	Fixed Centermount Legrest, Long <i>(Includes Calf Pads) (13.5" to 18.5" Leg Length)</i>	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Power Centermount Legrest		\$
400H109	SEDEO® Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation)</i> <i>(3" to 17.5" Leg Length)</i>	\$3,435.00
Centermount Legrest Accessories		\$
Power Centermount Footplates		
400F127	Dual, Height, Width and Angle Adjustable, Standard <i>(14"W x 10"D)</i>	\$284.00
400F128	Dual, Height, Width and Angle Adjustable, Large <i>(16"W x 12"D)</i>	\$284.00
Footplate Accessories (Available with Power Centermounts)		
400D168	Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i>	\$315.00
400D169	Footplate Connector Link <i>(Attaches footplates together)</i>	\$96.00
400D157	Gel Footpads <i>(Pair, LH and RH)</i>	\$401.00
Angle Adjustable / Fixed Centermount Footplates		
400F45	Single/Dual Footplate	NC
Power Centermount Calf Pads		
400D134	Calf Pads, Pair, Standard Integral Foam	NC
400D135	Calf Pads, Pair, Padded <i>(N/A w/ std. footplate)</i>	\$356.00
Angle Adjustable Centermount Calf Pads		
400D149	Calf Pads, Pair, Standard Integral Foam	NC
400D150	Extra Large Single Calf Pad	\$356.00
SEDEO® PRO™ Swing-Away Legrest (Pair) (Footplates, calf supports & options must be selected separately.)		\$
400H177	Legrests, Swing-Away, 70°, SEDEO® PRO™ <i>(14" to 21" Leg Length)</i>	NC
400H178	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ <i>(Includes Padded Calf Pads)</i> <i>(14.5" to 21.5" Leg Length)</i>	\$380.00
400H179	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ <i>(Includes Padded Calf Pads)</i> <i>(15.5" to 22.5" Leg Length)</i>	\$3,435.00
400HI76	Omit Legrests	NC
Swing-Away Legrest Footplates		
400F369	Angle Adjustable Flip-Up Footplates	\$284.00
Swing-Away Legrest Calf Supports		
400D142	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, depending on chair width, x 9.5" Tall)</i>	\$512.00
Swing-Away Legrest Options		
400D144	Heel Loops	NC
400KG4	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$282.00
SEDEO® Lower Extremity Swing-Away Support		
400VD01	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$512.00
400VD02	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$512.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Pelvic Positioning Belts		
400A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
400A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
400A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$103.00
400A845	Bodypoint – Evoflex Pelvic Stabilizer Belt	\$216.00
400A846	Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	\$271.00
Accessories		
Ventilator Option		
400VT12	Vent Tray Fixed (<i>Pulmonetic LTV Models</i>)	\$609.00
O2 Tank Option		
400A718	O2 Tank Holder Mounted to Backrest (<i>D and E Styles</i>)	\$391.00
Wheelchair Tray Table		
4007A09	Tray, Swing-Away	\$381.00
Backpack		
400A1	Black Utility	\$112.00
Medial Thigh Support		
4008A80	Medial Thigh Support	\$244.00
4008A81	Pommel Swing-Away Hardware	\$358.00
Misc. Accessories		
4008A59	Transfer Bar (<i>Pair</i>)	\$493.00
4008A52	Ram Mount Phone	\$77.00
4008A53	Ram Mount Tablet	\$308.00
4008A54	Ram Mount Cup Holder	\$77.00
4008A58	Rear View Mirror	\$77.00
400A611	Crutch Holder	\$115.00
4008A73	Bag for Armrest	\$115.00
4009A10	Utility Hooks	\$58.00

JOYSTICK / ELECTRONICS				
OPTION		TYPE	DESCRIPTION	RETAIL
Electronics Package				\$
400EL01		VR-2	VR-2 or R-net Electronics	NC
400EL02			LiNX Electronics <i>(SitSmarter not available on Tilt Only Configuration)</i>	NC
VR-2 Joystick				\$
400J127		VR-2	VR-2 Joystick	NC
R-net Joystick <i>(Remote)</i>				\$
400J121		R-net	Standard LED Joystick <i>(w/o Screen Display)</i>	NC
400J116		R-net	Color Joystick <i>(2.3" Color LCD Screen Display, CJSM1)</i>	NC
400J118		R-net	High Resolution Color Joystick, (CJSM2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,162.00
400J119		R-net	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
LiNX Joystick				\$
400J143		LiNX	REM110 - Drive Only Remote Joystick <i>(Adaptive Load Compensation)</i>	NC
400J144		LiNX	REM200 Series - LED Joystick <i>(Adaptive Load Compensation, LED icon-based menu, Two Rocker Switches)</i>	NC
400J145		LiNX	REM400 - Touch Screen Color Display Joystick <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation)</i>	\$1,339.00
400J146		LiNX	REM400 - Touch Screen Color Display Joystick with Toggle Switches <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, & Adaptive Load Compensation)</i>	\$1,695.00
400J147		LiNX	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
Joystick Mounting				\$
400JM34			Fixed on Armrest	NC
400JM35			Parallel Swing-Away <i>(Flip-Back - Swing In/Out, Single Post - Swing Out)</i>	\$447.00
400JM36			Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>	\$447.00
Joystick Mount (Location)				\$
400HM16			Right Hand Armrest	NC
400HM17			Left Hand Armrest	NC
Display Language Selection				\$
400L1			English	NC
400L3			French	NC
400L7			Spanish	NC

JOYSTICK / ELECTRONICS

OPTION		TYPE	DESCRIPTION	RETAIL
Joystick Handle (Optional)				\$
400JH13		Both	Ball Handle	\$110.00
400JH15		R-net	T-Handle	\$110.00
400JH24		Both	Bodypoint – Rubber Dome Joystick Handle	\$216.00
400JH25		Both	Bodypoint – U-Shaped Joystick Handle (3")	\$216.00
400JH26		Both	Bodypoint – Grooved Mushroom Joystick Handle	\$216.00
400JH30		LiNX	Chin Cup	\$179.00
400JH31		LiNX	Golf Ball	\$179.00
Electronic Programmers				\$
400PR13		R-net	R-net PC Wireless Programmer	\$660.00
400PR14		LiNX	LiNX Access Key (LAK)	\$356.00
Power Modules				\$
400ER12		VR-2	Power Module, VR-2 90A, Non-Expandable	NC
400ER22		R-net	Power Module, R-net, 90A, Non-Expandable <i>(Only available with Single Power)</i>	NC
400ER13		R-net	Power Module, R-net, 90A, Expandable <i>(Required for Specialty Controls or 2+ Actuators)</i>	\$102.00
400EX11		R-net	Expandable 90A Controller - R-net	NC
400EX12		R-net	Harness for 90A Expandable Controller	NC
400ER24		LiNX	Power Module, LiNX, 120A, Non-Expandable <i>(Only available with Single Power)</i>	NC
400ER25		LiNX	Power Module, LiNX, 120A, Expandable <i>(Required for Specialty Controls or 2+ Actuators)</i>	\$102.00
400EX19		LiNX	Expandable 120A Controller - LiNX	NC
400EX20		LiNX	Harness for 120A Expandable Controller - LiNX	NC
Seating Control				\$
Powered Seating Control				
400E040			Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
400E037			Thru-drive Control for 1 Actuator	\$1,443.00
400E038			Thru-drive Control for 2+ Actuators	\$2,903.00
Switch Control <i>(Available with Omit Thru-Drive only)</i>				
400E033		R-net	Dual Toggle for 1 Actuator	NC
400E034		Both	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
400E041		R-net	Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more	NC
400E042		LiNX	LiNX 10 Way Switch Box for Actuator Control for 1 Actuator or more	NC
Additional Power Seating Control <i>(Available with Thru-Drive only)</i>				
400QG11		R-net	Ctrl+5 Switch Box with Voice Activation	\$1,078.00
400QG16		LiNX	LiNX 10 Way Switch Box for Actuator Control	NC
Dual Toggle Mounting Location				
400TN1			Toggle Mount, Right	NC
400TN2			Toggle Mount, Left	NC

JOYSTICK / ELECTRONICS

OPTION	TYPE	DESCRIPTION	RETAIL
--------	------	-------------	--------

Power Seating Control Mounting Location *(Will swing away if mounted on the same side as the joystick)*

400EY19		Right Armrest	NC
400EY20		Left Armrest	NC
400EY21		Back Cane, Right	NC
400EY24		Back Cane, Left	NC

Ctrl+5 Switch Type

400EY22		Push Buttons	NC
400EY23		Toggle Switches	\$356.00

Attendant Controls

\$

4000A05	R-net	Dual Control Module for Attendant	\$1,225.00
4000A06	R-net	R-net Color LCD (CJSM1), Control Attendant Only	\$1,225.00
4000A11	LiNX	LiNX Attendant Control Unit (ACU200)	\$1,225.00

Attendant Mounting Location

\$

400DN5		Right Back Cane	NC
400DN6		Left Back Cane	NC

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

400PM12		Include Assignable Buttons	NC
---------	--	----------------------------	----

Power Centermount Footplate Extension Option

400IE03		Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation</i>	\$592.00
---------	--	--	----------

Memory Seating Positions

Memory Seating Positions - SitSmarter - allows for positions to be saved and then can be accessed with a single command, moving all related actuators at once to the preset location. Positions can be locked for clinical goals and have others unlocked for easily adjustable positions for comfort and functionality that can be adjusted by the user. Angles can be shown on screen to provide live feedback.

400PM19	LiNX	Memory Seating Positions – SitSmarter <i>(SitSmarter not compatible with Tilt Only configurations)</i>	NC
---------	------	--	----



Q300 M with SEDEO® LITE™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Signature

