



Quickie QRI
May 2025
EIR5

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ Prov: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ Prov: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

| ✓ | Option # | Definition | Price | ✓ | Option # | Definition | Price |
|-----------------------------------|----------|------------|-------|---|----------|------------|-------|
| 1. FRAME STYLE & SETUP | | | | | | | |

Base Model

| | | | |
|--------------------------|---------|--------------------------------------|------------|
| <input type="checkbox"/> | 202WM81 | QRI (265 lbs) Classic Frame Style | \$3,013.00 |
|--------------------------|---------|--------------------------------------|------------|

Transit

| | | | |
|--------------------------|--------|---------|----|
| <input type="checkbox"/> | 202TO1 | Transit | NC |
|--------------------------|--------|---------|----|

Frame Angle (see Figure 1)

| | | | |
|--------------------------|--------|-----|----|
| <input type="checkbox"/> | 202FA1 | 75° | NC |
| <input type="checkbox"/> | 202FA8 | 80° | NC |
| <input type="checkbox"/> | 202FA2 | 85° | NC |

Frame Inset Per Side (See Figure 2)

| | | | |
|--------------------------|--------|----------------------|----|
| <input type="checkbox"/> | 202FF1 | 0" | NC |
| <input type="checkbox"/> | 202FF2 | 1" (N/A w/ 12" wide) | NC |

| Frame Inset | Footrest Width |
|-------------|------------------------------|
| 0 | Subtract 0 from Frame Width |
| 1" | Subtract 2" from Frame Width |

Figure 1





Figure 2



Seat Width

Measure from outside of seat tube to outside of opposite seat tube.

| | | | |
|--------------------------|--------|-----|----------|
| <input type="checkbox"/> | 202S3 | 12" | NC |
| <input type="checkbox"/> | 202S4 | 13" | NC |
| <input type="checkbox"/> | 202S5 | 14" | NC |
| <input type="checkbox"/> | 202S6 | 15" | NC |
| <input type="checkbox"/> | 202S7 | 16" | NC |
| <input type="checkbox"/> | 202S8 | 17" | NC |
| <input type="checkbox"/> | 202S9 | 18" | NC |
| <input type="checkbox"/> | 202S10 | 19" | NC |
| <input type="checkbox"/> | 202S48 | 20" | \$308.00 |

Seat Depth

Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

| | | | |
|--------------------------|---------|-----|----------|
| <input type="checkbox"/> | 202SD5 | 12" | NC |
| <input type="checkbox"/> | 202SD6 | 13" | NC |
| <input type="checkbox"/> | 202SD7 | 14" | NC |
| <input type="checkbox"/> | 202SD8 | 15" | NC |
| <input type="checkbox"/> | 202SD9 | 16" | NC |
| <input type="checkbox"/> | 202SD10 | 17" | NC |
| <input type="checkbox"/> | 202SD11 | 18" | NC |
| <input type="checkbox"/> | 202SD12 | 19" | NC |
| <input type="checkbox"/> | 202SD13 | 20" | \$308.00 |


| ✓ | Option # | Definition | Price | ✓ | Option # | Definition | Price |
|------------------------------------|----------|------------|-------|---|----------|------------|-------|
| 2. FRONT SEAT HEIGHT SET UP | | | | | | | |

| Front Seat Height | | | | | | | |
|--------------------------|---------|-------|----|--------------------------|---------|-------|----|
| See Figure 3 | | | | | | | |
| <input type="checkbox"/> | 202SH12 | 15.0" | NC | <input type="checkbox"/> | 145SH20 | 18.5" | NC |
| <input type="checkbox"/> | 202SH13 | 15.5" | NC | <input type="checkbox"/> | 145SH21 | 19.0" | NC |
| <input type="checkbox"/> | 202SH14 | 16.0" | NC | <input type="checkbox"/> | 145SH23 | 19.5" | NC |
| <input type="checkbox"/> | 202SH15 | 16.5" | NC | <input type="checkbox"/> | 145SH24 | 20.0" | NC |
| <input type="checkbox"/> | 202SH16 | 17.0" | NC | | | | |
| <input type="checkbox"/> | 202SH17 | 17.5" | NC | | | | |
| <input type="checkbox"/> | 202SH18 | 18.0" | NC | | | | |

| Caster Wheels | | | |
|--------------------------|---------|------------------------------|----------|
| <input type="checkbox"/> | 202CW3 | 4" Low Profile Polyurethane | NC |
| <input type="checkbox"/> | 202CW44 | 4" x 1.25" Semi Pneumatic | NC |
| <input type="checkbox"/> | 202CW64 | 4" x 1.5" Aluminum Soft Roll | \$110.00 |
| <input type="checkbox"/> | 202CW5 | 5" Low Profile Polyurethane | NC |
| <input type="checkbox"/> | 202CW65 | 5"x 1" Aluminum Polyurethane | \$126.00 |
| <input type="checkbox"/> | 202CW7 | 5"x 1.5" Semi Pneumatic | \$110.00 |
| <input type="checkbox"/> | 202CW66 | 5"x 1.5" Aluminum Soft Roll | \$110.00 |
| <input type="checkbox"/> | 202CW34 | 6" Polyurethane | NC |
| <input type="checkbox"/> | 202CW35 | 6"x 1.5" Semi Pneumatic | \$110.00 |
| <input type="checkbox"/> | 202CW68 | 6"x 1.5" Aluminum Soft Roll | \$110.00 |

| Frame Depth | Max Seat Squeeze |
|-------------|------------------|
| 12"/13" | 2.5" |
| 14"/15" | 3.0" |
| 16"/17" | 3.5" |
| 18"/19" | 4.0" |
| 20" | 4.5" |

Figure 3



| Caster Fork | | | |
|-----------------------------|---------|-------------------------------------|----------|
| Frog Legs N/A w/ 6" casters | | | |
| <input type="checkbox"/> | 202CF26 | Standard | NC |
| <input type="checkbox"/> | 202CF53 | Frog Legs 0 weight to 100 lbs | \$574.00 |
| <input type="checkbox"/> | 202CF54 | Frog Legs 101 lbs weight to 190 lbs | \$574.00 |
| <input type="checkbox"/> | 202CF55 | Frog Legs 191 lbs weight to 265 lbs | \$574.00 |

3. FOOTPLATE OPTIONS

| Footrest | | | |
|--|---------|-----------------------------------|----------|
| Leg strap not included, see Footplate Option section | | | |
| <input type="checkbox"/> | 202F308 | Aluminum Tubular | NC |
| <input type="checkbox"/> | 202F309 | Aluminum Tubular w/ Snap on Cover | NC |
| <input type="checkbox"/> | 202F310 | Aluminum Platform (Ang. Adj.) | \$152.00 |

| Footrest Options | | | |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 202DI19 | Leg Strap | \$61.00 |
| <input type="checkbox"/> | 202A129 | Extra Leg Strap | \$83.00 |
| <input type="checkbox"/> | 202DI06 | Positioning Plates for Platform Footrest | \$78.00 |
| <input type="checkbox"/> | 202A500 | Gel Cover Hanger/Frame | \$162.00 |

| Adaptive Platform Spacer | | | |
|---------------------------|--------|-----------------------------|----------|
| N/A with Tubular Footrest | | | |
| <input type="checkbox"/> | 202AD1 | Adaptive Platform 2" Spacer | \$182.00 |
| <input type="checkbox"/> | 202AD2 | Adaptive Platform 4" Spacer | \$182.00 |
| <input type="checkbox"/> | 202AD3 | Adaptive Platform 6" Spacer | \$182.00 |

Tubular and Platform Footrest Caster/Fork Availability

Lower Leg Length Range

| Frame Angle | Footrest Type | Range | FSH Range 15"-17" | | FSH Range 17.5"-20" | |
|-------------|---------------------------|-------|----------------------|------|------------------------|------|
| | | Inset | 0" | 1" | 0" | 1" |
| 75° | Platform | Min | 12.3 | 12.3 | 14.3 | 14.3 |
| | | Max | 16.0 | 16.0 | 18.0 | 18.0 |
| | Tubular | Min | 12.3 | 12.3 | 14.3 | 14.3 |
| | | Max | 15.5 | 15.5 | 17.5 | 17.5 |
| | Long Platform and Tubular | Min | 12.3 | 13.9 | 14.3 | 14.3 |
| | | Max | 17.5 | 17.5 | 19.5 | 19.5 |
| 80° | Platform | Min | 12.2 | 12.1 | 14.2 | 14 |
| | | Max | 16.1 | 16.1 | 18.1 | 18.1 |
| | Tubular | Min | 12.2 | 12.2 | 14.2 | 14.2 |
| | | Max | 15.6 | 15.6 | 17.6 | 17.6 |
| | Long Platform and Tubular | Min | 12.2 | 14.2 | 14.2 | 14.2 |
| | | Max | 17.8 | 17.8 | 19.8 | 19.8 |
| 85° | Platform | Min | 12.2 | 12.5 | 14.2 | 14 |
| | | Max | 16.1 | 16.1 | 18.1 | 18.1 |
| | Tubular | Min | 12.2 | 12.5 | 14.2 | 14.2 |
| | | Max | 15.5 | 15.5 | 17.5 | 17.5 |
| | Long Platform and Tubular | Min | 12.2 | 14.5 | 14.2 | 14.5 |
| | | Max | 17.7 | 17.7 | 19.7 | 19.7 |

Quote or order acknowledgement will display standard or long footrest used with configuration of chair. Subtract appropriate length if using platform spacers.

| | Caster | Fork | Inset | |
|-----------------|--------|---------------|-----------|-----------|
| | | | 0" | 1" |
| 75° Frame Angle | 4 | Standard Fork | | |
| | | Frog Legs | | |
| | 5 | Standard Fork | | |
| 80° Frame Angle | 4 | Standard Fork | | |
| | | Frog Legs | | |
| | 5 | Standard Fork | | |
| | | Frog Legs | N/A | |
| | 6 | Standard Fork | N/A w/1.5 | |
| | | Frog Legs | N/A | N/A |
| 85° Frame Angle | 4 | Standard Fork | | |
| | | Frog Legs | N/A | N/A w/1.5 |
| | 5 | Standard Fork | N/A w/1.5 | |
| | | Frog Legs | N/A | N/A |
| | 6 | Standard Fork | N/A | N/A |
| | | Frog Legs | N/A | N/A |

N/A w/ 1.5 = N/A with the 1.5" wide caster. Blank = available. No Go's may still apply that are not listed based on actual configuration of the chair.

| ✓ | Option # | Definition | Price | ✓ | Option # | Definition | Price |
|-----------------------------------|----------|------------|-------|---|----------|------------|-------|
| 4. REAR SEAT HEIGHT SET UP | | | | | | | |

| Rear Seat Height | | | |
|--|---------|-------|----|
| <i>14" and 14.5" only available w/ 22" rear wheels</i> | | | |
| <input type="checkbox"/> | 202RH13 | 14.0" | NC |
| <input type="checkbox"/> | 202RH14 | 14.5" | NC |
| <input type="checkbox"/> | 202RH15 | 15.0" | NC |
| <input type="checkbox"/> | 202RH16 | 15.5" | NC |
| <input type="checkbox"/> | 202RH17 | 16.0" | NC |
| <input type="checkbox"/> | 202RH18 | 16.5" | NC |
| <input type="checkbox"/> | 202RH19 | 17.0" | NC |
| <input type="checkbox"/> | 202RH26 | 17.5" | NC |
| <input type="checkbox"/> | 202RH20 | 18.0" | NC |
| <input type="checkbox"/> | 202RH27 | 18.5" | NC |
| <input type="checkbox"/> | 202RH21 | 19.0" | NC |
| <input type="checkbox"/> | 202RH37 | 19.5" | NC |
| <input type="checkbox"/> | 202RH22 | 20.0" | NC |

| Rear Wheel Size | | | |
|---|---------|---------------------------------------|------------|
| <i>Ultra Lightweight Wheels only available with aluminum anodized handrim</i> | | | |
| <i>Contact Built 4 Me for Spinergy color options</i> | | | |
| <input type="checkbox"/> | 202RW8 | 22" Mag Wheel | NC |
| <input type="checkbox"/> | 202RW1 | 24" Mag Wheel | NC |
| <input type="checkbox"/> | 202RW91 | 22" Lite Spoke Wheel | \$320.00 |
| <input type="checkbox"/> | 202RW92 | 24" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 202RW93 | 25" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 202RW94 | 26" Lite Spoke Wheel | \$126.00 |
| <input type="checkbox"/> | 202RW06 | 22" Ultra Lightweight Spoke Wheel 2.0 | \$357.00 |
| <input type="checkbox"/> | 202RW10 | 24" Ultra Lightweight Spoke Wheel 2.0 | \$357.00 |
| <input type="checkbox"/> | 202RW16 | 25" Ultra Lightweight Spoke Wheel 2.0 | \$357.00 |
| <input type="checkbox"/> | 202RW28 | 22" Spinergy SPOX Black | \$864.00 |
| <input type="checkbox"/> | 202RW95 | 24" Spinergy SPOX Black | \$864.00 |
| <input type="checkbox"/> | 202RW96 | 25" Spinergy SPOX Black | \$864.00 |
| <input type="checkbox"/> | 202RW97 | 26" Spinergy SPOX Black | \$864.00 |
| <input type="checkbox"/> | 202RW32 | 22" Spinergy LX Black | \$1,111.00 |
| <input type="checkbox"/> | 202RW98 | 24" Spinergy LX Black | \$1,111.00 |
| <input type="checkbox"/> | 202RW99 | 25" Spinergy LX Black | \$1,111.00 |
| <input type="checkbox"/> | 202RW40 | 26" Spinergy LX Black | \$1,111.00 |

| Rear Wheel Tire | | | |
|--------------------------------|---------|----------------------------------|----------|
| <i>See Availability Matrix</i> | | | |
| <input type="checkbox"/> | 202RT1 | Pneumatic | NC |
| <input type="checkbox"/> | 202RT2 | Pneumatic w/ Airless Insert | \$247.00 |
| <input type="checkbox"/> | 202RT3 | Full Poly | NC |
| <input type="checkbox"/> | 202RT6 | High-Pressure Clincher | \$105.00 |
| <input type="checkbox"/> | 202RT27 | Knobby V- Trak | \$121.00 |
| <input type="checkbox"/> | 202RT21 | Schwalbe Marathon Plus Evolution | \$157.00 |
| <input type="checkbox"/> | 202RT38 | SHOX Performance Poly | \$146.00 |

| Short Axle Plate/Wheel Size Matrix | | | | |
|---|------|------|------|------|
| Rear Seat to Floor Height | 22" | 24" | 25" | 26" |
| | 14.0 | | | |
| | 14.5 | 15.0 | 15.0 | 16.0 |
| | 16.0 | 15.5 | 15.5 | 17.5 |
| | 16.5 | 17.0 | 17.0 | 18.0 |
| | 17.0 | 17.5 | 17.5 | 18.5 |
| | 17.5 | 18.0 | 18.0 | 19.0 |

| Rear Tire & Wheel Matrix | | | | |
|--|-------|-------|-------------|-------------------|
| Tire | Wheel | Mag | Lite Spoke | Ultra Lightweight |
| Pneumatic | | 22,24 | 22,24,26 | |
| Pneumatic w/ Airless Insert | | 22,24 | 22,24 | |
| Full Poly | | 22,24 | 22,24 | |
| High-Pressure Clincher | | | 22,24,25,26 | 22,24,25 |
| Knobby V Trak by Kenda | | | | 24 |
| Schwalbe Marathon | | | 24,25,26 | 24,25,26 |
| SHOX Performance Poly | | 24 | 24,25 | 24,25 |
| 22" = 489 24" = 540 25" = 559 26" = 590 | | | | |

| Handrim | | | |
|--|---------|--------------------------------|----------|
| <i>Plastic not available on all wheels in all sizes. Natural Fits are N/A w/ 26"</i> | | | |
| <input type="checkbox"/> | 202HR1 | Aluminum Anodized | NC |
| <input type="checkbox"/> | 202HR4 | Plastic | \$146.00 |
| <input type="checkbox"/> | 202HR22 | Natural Fit - Standard Grip | \$374.00 |
| <input type="checkbox"/> | 202HR23 | Natural Fit - Super Grip | \$434.00 |
| <input type="checkbox"/> | 202HR30 | Natural Fit LT - Standard Grip | \$521.00 |
| <input type="checkbox"/> | 202HR31 | Natural Fit LT - Super Grip | \$592.00 |

| Center of Gravity Preset | | | |
|---|--------|------|----|
| <i>Measure from front of back post to center of rear axle</i> | | | |
| <input type="checkbox"/> | 202CC2 | 1.0" | NC |
| <input type="checkbox"/> | 202CC3 | 1.5" | NC |
| <input type="checkbox"/> | 202CC4 | 2.0" | NC |
| <input type="checkbox"/> | 202CC5 | 2.5" | NC |
| <input type="checkbox"/> | 202CC6 | 3.0" | NC |

| Camber | | | |
|--------------------------|--------|----|----|
| <input type="checkbox"/> | 202CI1 | 0° | NC |
| <input type="checkbox"/> | 202CI2 | 3° | NC |

| Axle | | | |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 202AX1 | Quick Release Steel | NC |
| <input type="checkbox"/> | 202AX13 | Quad Quick Release <i>Only available with Mag wheels</i> | \$182.00 |

| Axle Plates | | | |
|---|---------|-------------------------------|----|
| <i>See Short Axle Plate Matrix for RSH Availability</i> | | | |
| <input type="checkbox"/> | 202AL40 | Standard Short (COG 1" to 3") | NC |
| <input type="checkbox"/> | 202AL43 | Standard Long (COG 1" to 3") | NC |

| Wheel Lock | | | |
|---|---------|-----------------|----------|
| <i>Not all wheel locks available in all seat depths, wheel sizes and COG settings</i> | | | |
| <input type="checkbox"/> | 202WL1 | High Mount Push | NC |
| <input type="checkbox"/> | 202WL2 | High Mount Pull | NC |
| <input type="checkbox"/> | 202WL54 | Ergo Scissor | NC |
| <input type="checkbox"/> | 202WL30 | Compact | \$162.00 |
| <input type="checkbox"/> | 202WL10 | Omit Wheel Lock | NC |
| <input type="checkbox"/> | 202WL12 | Do Not Mount | NC |

| Extension Handle | | | |
|--|--------|--------------------------|---------|
| <i>Only available w/ High Mount Push or Pull Wheel Locks</i> | | | |
| <input type="checkbox"/> | 202EH7 | Extension handles (Pair) | \$73.00 |

5. ARMREST

| Armrest - Padded Swing Away | | | |
|---|---------|--|----------|
| <i>Preset to 10". Some limitations may exist with integral and folding push handles Range 9"-13".</i> | | | |
| <input type="checkbox"/> | 202AA18 | Padded Swing-Away - Standard | \$278.00 |
| Armrest - Height Adjustable | | | |
| <input type="checkbox"/> | 202AA11 | Single Post Ht. Adjust. Std (8.5"-13.5") | \$339.00 |
| <input type="checkbox"/> | 202AA29 | Single Post Ht. Adjust. Low (7"-11") | \$339.00 |

| Armpad - For Height Adjustable | | | |
|--------------------------------|---------|-----------------------------|----------|
| <input type="checkbox"/> | 202AM23 | Desk Length - Classic (10") | NC |
| <input type="checkbox"/> | 202AM24 | Full Length - Classic (14") | NC |
| Omit Armrest | | | |
| <input type="checkbox"/> | 202AA14 | Omit Armrest | NC |
| Armrest Accessories | | | |
| <input type="checkbox"/> | 2028A05 | Gel Cover Swing Away Arm | \$146.00 |

6. BACKREST

| Back Type | | | |
|----------------------------|---------|---|----------|
| <input type="checkbox"/> | 202BT82 | Folding Lock Down Angle Adj | NC |
| <input type="checkbox"/> | 202BT3 | Non-Folding Backrest | NC |
| Push Handles | | | |
| <input type="checkbox"/> | 202PH6 | Integral Push Handles | NC |
| <input type="checkbox"/> | 202PH5 | Fold-Down Push Handles <i>N/A w/ 10-14" Back Height or Standard Upholstery.</i> | \$168.00 |
| <input type="checkbox"/> | 202PH2 | Omit Push Handles | NC |
| Height Adjustable Backrest | | | |
| <input type="checkbox"/> | 202B110 | 10" - 14" Back Height | NC |
| <input type="checkbox"/> | 202B127 | 12" - 16" Back Height | NC |
| <input type="checkbox"/> | 202B111 | 14" - 18" Back Height | NC |
| <input type="checkbox"/> | 202B137 | 16" - 20" Back Height | NC |
| Back Upholstery | | | |
| <input type="checkbox"/> | 202BU1 | Standard Upholstery | NC |
| <input type="checkbox"/> | 202BU70 | Backrest Tension Adj Ballistic Nylon <i>*reference 186BU70 for billing</i> | \$190.00 |
| <input type="checkbox"/> | 202BU71 | Backrest Tension Adj 3DX Vented <i>*reference 186BU71 for billing</i> | \$471.00 |
| <input type="checkbox"/> | 202BU14 | Omit Upholstery (<i>Use w/ JAY Backs</i>) | NC |

| JAY J3 Back | | | |
|--------------------------|---------|-------------|--------|
| <input type="checkbox"/> | 202JB01 | JAY J3 Back | See OF |

7. SEATING

| Seat Sling | | | |
|--------------------------|--------|---------------------|----|
| <input type="checkbox"/> | 202SL1 | Standard Seat Sling | NC |
| <input type="checkbox"/> | 202SL6 | Omit Seat Sling | NC |

| JAY Cushions (Optional) | | | |
|---|---------|--|------------|
| <i>To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth</i> | | | |
| <input type="checkbox"/> | 202JC20 | JAY Bundled Cushion Cushion Name or Part # Width: _____ Depth: _____ | Price List |

| ✓ | Option # | Definition | Price | ✓ | Option # | Definition | Price |
|---|----------|------------|-------|---|----------|------------|-------|
|---|----------|------------|-------|---|----------|------------|-------|

8. Key Options & Accessories

| Anti Tips | | | |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 202AT1 | Anti Tips (Pair) | \$121.00 |
| <input type="checkbox"/> | 202AT3 | Omit Rear Anti Tip | NC |
| Side Guards | | | |
| <input type="checkbox"/> | 202SG33 | Plastic - Standard (8.5" High) | \$188.00 |
| <input type="checkbox"/> | 202SG29 | Plastic - Low (6.5" High) | \$188.00 |
| <input type="checkbox"/> | 202SG35 | Aluminum N/A w/ height adj arms, push/pull wheel locks only available in black | \$422.00 |
| <input type="checkbox"/> | 202SG36 | Aluminum w/ Fender N/A w/ height adj arms, push/pull wheel locks only available in black | \$495.00 |

| Pack Accessories | | | |
|--------------------------|-------|-----------------------------|----------|
| <input type="checkbox"/> | 202A1 | Quickie Backpack Black | \$112.00 |
| <input type="checkbox"/> | 202A2 | Quickie Kids Backpack Black | \$105.00 |
| <input type="checkbox"/> | 202A8 | Seat Pouch Black | \$83.00 |

| Positioning Belt | | | |
|--------------------------|---------|-------------------------------------|---------|
| <input type="checkbox"/> | 202A141 | 1.5" Auto Buckle Universal w/clamps | \$97.00 |
| <input type="checkbox"/> | 202A120 | Positioning Belt Slide Clamps | \$61.00 |

| Spoke Guards | | | |
|--------------------------|---------|---------------------------|----------|
| <input type="checkbox"/> | 202A302 | Quickie Clear Spoke Guard | \$157.00 |
| <input type="checkbox"/> | 202A23 | Quickie Black Spoke Guard | \$157.00 |

| Miscellaneous | | | |
|--------------------------|----------|--|----------|
| <input type="checkbox"/> | 202DI07 | Neoprene Impact Guards | \$61.00 |
| <input type="checkbox"/> | 202A2100 | Crutch Holder Only available with Integral or Fold Down push handles | \$115.00 |
| <input type="checkbox"/> | 202A85 | Bolt-On Push Handles | \$139.00 |
| <input type="checkbox"/> | 202A83 | Tool Kit | \$43.00 |

9. COLOR

| Standard Frame Colors | | | |
|--------------------------|---------|--------------------|----|
| <input type="checkbox"/> | 202C211 | Celestial Sapphire | NC |
| <input type="checkbox"/> | 202C212 | Desert Sage | NC |
| <input type="checkbox"/> | 202C213 | Desert Tan | NC |
| <input type="checkbox"/> | 202C214 | Gloss Black | NC |
| <input type="checkbox"/> | 202C215 | Matte Black | NC |
| <input type="checkbox"/> | 202C216 | Metallic Shadow | NC |
| <input type="checkbox"/> | 202C217 | Midnight Blue | NC |
| <input type="checkbox"/> | 202C218 | Mystic Purple | NC |
| <input type="checkbox"/> | 202C219 | Petal Pink | NC |
| <input type="checkbox"/> | 202C220 | Pink Power | NC |
| <input type="checkbox"/> | 202C221 | Red Pepper | NC |
| <input type="checkbox"/> | 202C222 | Royal Blue | NC |
| <input type="checkbox"/> | 202C223 | Salute Green | NC |
| <input type="checkbox"/> | 202C224 | Scarlet Noir | NC |
| <input type="checkbox"/> | 202C225 | Silver Ore | NC |
| <input type="checkbox"/> | 202C226 | Solar Orange | NC |
| <input type="checkbox"/> | 202C227 | Starlight White | NC |
| <input type="checkbox"/> | 202C228 | Storm Gray | NC |
| <input type="checkbox"/> | 202C229 | Turbo Gray | NC |
| <input type="checkbox"/> | 202C230 | Velocity Yellow | NC |

| Personalized Frame Color | | | |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 202C112 | Desert Camouflage | \$368.00 |
| <input type="checkbox"/> | 202C43 | Mossy Oak Camouflage | \$368.00 |
| <input type="checkbox"/> | 202C42 | Stars & Stripes | \$368.00 |
| <input type="checkbox"/> | 202C108 | Zebra (Choose base color) N/A w/ Gloss Black and Matte Black. | \$368.00 |

Additional Options Through Built for Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8