



Quickie QRi

July 2025

EIR5

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ Prov: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ Prov: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME STYLE & SETUP

Base Model

☐	202WM81	QRi (265 lbs) Classic Frame Style	\$3,013.00
--------------------------	---------	--------------------------------------	------------

Transit

☐	202TO1	Transit	NC
--------------------------	--------	---------	----

Frame Angle (see Figure 1)

☐	202FA1	75°	NC
☐	202FA8	80°	NC
☐	202FA2	85°	NC

Frame Inset Per Side (See Figure 2)

☐	202FF1	0"	NC
☐	202FF2	1" (N/A w/ 12" wide)	NC

Frame Inset	Footrest Width
0	Subtract 0 from Frame Width
1"	Subtract 2" from Frame Width

Figure 1





Figure 2



Seat Width

Measure from outside of seat tube to outside of opposite seat tube.

☐	202S3	12"	NC
☐	202S4	13"	NC
☐	202S5	14"	NC
☐	202S6	15"	NC
☐	202S7	16"	NC
☐	202S8	17"	NC
☐	202S9	18"	NC
☐	202S10	19"	NC
☐	202S48	20"	\$308.00

Seat Depth

Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.


☐	202SD5	12"	NC
☐	202SD6	13"	NC
☐	202SD7	14"	NC
☐	202SD8	15"	NC
☐	202SD9	16"	NC
☐	202SD10	17"	NC
☐	202SD11	18"	NC
☐	202SD12	19"	NC
☐	202SD13	20"	\$308.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
2. FRONT SEAT HEIGHT SET UP							

Front Seat Height						
See Figure 3						
☐	202SH12	15.0"	NC	☐	145SH20	18.5"
☐	202SH13	15.5"	NC	☐	145SH21	19.0"
☐	202SH14	16.0"	NC	☐	145SH23	19.5"
☐	202SH15	16.5"	NC	☐	145SH24	20.0"
☐	202SH16	17.0"	NC			
☐	202SH17	17.5"	NC			
☐	202SH18	18.0"	NC			

Frame Depth	Max Seat Squeeze
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Figure 3



Caster Wheels				
☐	202CW3	4" Low Profile Polyurethane		NC
☐	202CW44	4" x 1.25" Semi Pneumatic		NC
☐	202CW64	4" x 1.5" Aluminum Soft Roll		\$110.00
☐	202CW5	5" Low Profile Polyurethane		NC
☐	202CW65	5"x 1" Aluminum Polyurethane		\$126.00
☐	202CW7	5"x 1.5" Semi Pneumatic		\$110.00
☐	202CW66	5"x 1.5" Aluminum Soft Roll		\$110.00
☐	202CW34	6" Polyurethane		NC
☐	202CW35	6"x 1.5" Semi Pneumatic		\$110.00
☐	202CW68	6"x 1.5" Aluminum Soft Roll		\$110.00

Caster Fork				
Frog Legs N/A w/ 6" casters				
☐	202CF26	Standard		NC
☐	202CF53	Frog Legs 0 weight to 100 lbs		\$574.00
☐	202CF54	Frog Legs 101 lbs weight to 190 lbs		\$574.00
☐	202CF55	Frog Legs 191 lbs weight to 265 lbs		\$574.00

3. FOOTPLATE OPTIONS

Footrest				
Leg strap not included, see Footplate Option section				
☐	202F308	Aluminum Tubular		NC
☐	202F309	Aluminum Tubular w/ Snap on Cover		NC
☐	202F310	Aluminum Platform (Ang. Adj.)		\$152.00

Footrest Options				
☐	202DI19	Leg Strap		\$61.00
☐	202A129	Extra Leg Strap		\$83.00
☐	202DI06	Positioning Plates for Platform Footrest		\$78.00
☐	202A500	Gel Cover Hanger/Frame		\$162.00

Adaptive Platform Spacer				
N/A with Tubular Footrest				
☐	202AD1	Adaptive Platform 2" Spacer		\$182.00
☐	202AD2	Adaptive Platform 4" Spacer		\$182.00
☐	202AD3	Adaptive Platform 6" Spacer		\$182.00

Tubular and Platform Footrest Caster/Fork Availability

Lower Leg Length Range

Frame Angle	Footrest Type	Range	FSH Range 15"-17"		FSH Range 17.5"-20"	
			Inset	0"	1"	0"
75°	Platform	Min		12.3	12.3	14.3
		Max		16.0	16.0	18.0
	Tubular	Min		12.3	12.3	14.3
		Max		15.5	15.5	17.5
	Long Platform and Tubular	Min		12.3	13.9	14.3
		Max		17.5	17.5	19.5
80°	Platform	Min		12.2	12.1	14.2
		Max		16.1	16.1	18.1
	Tubular	Min		12.2	12.2	14.2
		Max		15.6	15.6	17.6
	Long Platform and Tubular	Min		12.2	14.2	14.2
		Max		17.8	17.8	19.8
85°	Platform	Min		12.2	12.5	14.2
		Max		16.1	16.1	18.1
	Tubular	Min		12.2	12.5	14.2
		Max		15.5	15.5	17.5
	Long Platform and Tubular	Min		12.2	14.5	14.2
		Max		17.7	17.7	19.7

Quote or order acknowledgement will display standard or long footrest used with configuration of chair. Subtract appropriate length if using platform spacers.

		Caster	Fork	Inset	
				0"	1"
75° Frame Angle	4		Standard Fork		
			Frog Legs		
	5		Standard Fork		
			Frog Legs		
	6		Standard Fork		
			Frog Legs	N/A	N/A
80° Frame Angle	4		Standard Fork		
			Frog Legs		
	5		Standard Fork		
			Frog Legs	N/A	
	6		Standard Fork	N/A w/1.5	
			Frog Legs	N/A	N/A
85° Frame Angle	4		Standard Fork		
			Frog Legs	N/A	N/A w/1.5
	5		Standard Fork	N/A w/1.5	
			Frog Legs	N/A	N/A
	6		Standard Fork	N/A	N/A
			Frog Legs	N/A	N/A

N/A w/ 1.5 = N/A with the 1.5" wide caster. Blank = available. No Go's may still apply that are not listed based on actual configuration of the chair.

✓	Option #	Definition	Price	✓	Option #	Definition	Price
4. REAR SEAT HEIGHT SET UP							

Rear Seat Height			
<i>14" and 14.5" only available w/ 22" rear wheels</i>			
☐	202RH13	14.0"	NC
☐	202RH14	14.5"	NC
☐	202RH15	15.0"	NC
☐	202RH16	15.5"	NC
☐	202RH17	16.0"	NC
☐	202RH18	16.5"	NC
☐	202RH19	17.0"	NC
☐	202RH26	17.5"	NC
☐	202RH20	18.0"	NC
☐	202RH27	18.5"	NC
☐	202RH21	19.0"	NC
☐	202RH37	19.5"	NC
☐	202RH22	20.0"	NC

Rear Wheel Size			
<i>Ultra Lightweight Wheels only available with aluminum anodized handrim</i>			
<i>Contact Built 4 Me for Spinergy color options</i>			
☐	202RW8	22" Mag Wheel	NC
☐	202RW1	24" Mag Wheel	NC
☐	202RW91	22" Lite Spoke Wheel	\$320.00
☐	202RW92	24" Lite Spoke Wheel	NC
☐	202RW93	25" Lite Spoke Wheel	NC
☐	202RW94	26" Lite Spoke Wheel	\$126.00
☐	202RW06	22" Ultra Lightweight Spoke Wheel 2.0	\$357.00
☐	202RW10	24" Ultra Lightweight Spoke Wheel 2.0	\$357.00
☐	202RW16	25" Ultra Lightweight Spoke Wheel 2.0	\$357.00
☐	202RW28	22" Spinergy SPOX Black	\$864.00
☐	202RW95	24" Spinergy SPOX Black	\$864.00
☐	202RW96	25" Spinergy SPOX Black	\$864.00
☐	202RW97	26" Spinergy SPOX Black	\$864.00
☐	202RW32	22" Spinergy LX Black	\$1,111.00
☐	202RW98	24" Spinergy LX Black	\$1,111.00
☐	202RW99	25" Spinergy LX Black	\$1,111.00
☐	202RW40	26" Spinergy LX Black	\$1,111.00

Rear Wheel Tire			
<i>See Availability Matrix</i>			
☐	202RT1	Pneumatic	NC
☐	202RT2	Pneumatic w/ Airless Insert	\$247.00
☐	202RT3	Full Poly	NC
☐	202RT6	High-Pressure Clincher	\$105.00
☐	202RT27	Knobby V- Trak	\$121.00
☐	202RT21	Schwalbe Marathon Plus Evolution	\$157.00
☐	202RT38	SHOX Performance Poly	\$146.00

Short Axle Plate/Wheel Size Matrix				
Rear Seat to Floor Height	22"	24"	25"	26"
	14.0			
	14.5	15.0	15.0	16.0
	16.0	15.5	15.5	17.5
	16.5	17.0	17.0	18.0
	17.0	17.5	17.5	18.5
	17.5	18.0	18.0	19.0

Rear Tire & Wheel Matrix				
Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight
Pneumatic		22,24	22,24,26	
Pneumatic w/ Airless Insert		22,24	22,24	
Full Poly		22,24	22,24	
High-Pressure Clincher			22,24,25,26	22,24,25
Knobby V Trak by Kenda				24
Schwalbe Marathon			24,25,26	24,25,26
SHOX Performance Poly		24	24,25	24,25
22" = 489 24" = 540 25" = 559 26" = 590				

Handrim			
<i>Plastic not available on all wheels in all sizes. Natural Fits are N/A w/ 26"</i>			
☐	202HR1	Aluminum Anodized	NC
☐	202HR4	Plastic	\$146.00
☐	202HR22	Natural Fit - Standard Grip	\$374.00
☐	202HR23	Natural Fit - Super Grip	\$434.00

Center of Gravity Preset			
<i>Measure from front of back post to center of rear axle</i>			
☐	202CC2	1.0"	NC
☐	202CC3	1.5"	NC
☐	202CC4	2.0"	NC
☐	202CC5	2.5"	NC
☐	202CC6	3.0"	NC

Camber			
☐	202CI1	0°	NC
☐	202CI2	3°	NC

Axle			
☐	202AX1	Quick Release Steel	NC
☐	202AX13	Quad Quick Release <i>Only available with Mag wheels</i>	\$182.00

Axle Plates			
<i>See Short Axle Plate Matrix for RSH Availability</i>			
☐	202AL40	Standard Short (COG 1" to 3")	NC
☐	202AL43	Standard Long (COG 1" to 3")	NC

Wheel Lock			
<i>Not all wheel locks available in all seat depths, wheel sizes and COG settings</i>			
☐	202WL1	High Mount Push	NC
☐	202WL2	High Mount Pull	NC
☐	202WL54	Ergo Scissor	NC
☐	202WL30	Compact	\$162.00
☐	202WL10	Omit Wheel Lock	NC
☐	202WL12	Do Not Mount	NC

Extension Handle			
<i>Only available w/ High Mount Push or Pull Wheel Locks</i>			
☐	202EH7	Extension handles (Pair)	\$73.00

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

5. ARMREST

Armrest - Padded Swing Away			
<i>Preset to 10". Some limitations may exist with integral and folding push handles Range 9"-13".</i>			
☐	202AA18	Padded Swing-Away - Standard	\$278.00

Armrest - Height Adjustable			
☐	202AA11	Single Post Ht. Adjust. Std (8.5"-13.5")	\$339.00
☐	202AA29	Single Post Ht. Adjust. Low (7"-11")	\$339.00

Armpad - For Height Adjustable			
☐	202AM23	Desk Length - Classic (10")	NC
☐	202AM24	Full Length - Classic (14")	NC

Omit Armrest			
☐	202AA14	Omit Armrest	NC

Armrest Accessories			
☐	2028A05	Gel Cover Swing Away Arm	\$146.00

6. BACKREST

Back Type			
☐	202BT82	Folding Lock Down Angle Adj	NC
☐	202BT3	Non-Folding Backrest	NC

Push Handles			
☐	202PH6	Integral Push Handles	NC
☐	202PH5	Fold-Down Push Handles N/A w/ 10-14" Back Height or Standard Upholstery.	\$168.00
☐	202PH2	Omit Push Handles	NC

Height Adjustable Backrest			
☐	202B110	10" - 14" Back Height	NC
☐	202B127	12" - 16" Back Height	NC
☐	202B111	14" - 18" Back Height	NC
☐	202B137	16" - 20" Back Height	NC

Back Upholstery			
☐	202BU1	Standard Upholstery	NC
☐	202BU70	Backrest Tension Adj Ballistic Nylon <i>*reference 186BU70 for billing</i>	\$190.00
☐	202BU71	Backrest Tension Adj 3DX Vented <i>*reference 186BU71 for billing</i>	\$471.00
☐	202BU14	Omit Upholstery (Use w/ JAY Backs)	NC

JAY J3 Back			
☐	202JB01	JAY J3 Back	See OF

7. SEATING

Seat Sling			
☐	202SL1	Standard Seat Sling	NC
☐	202SL6	Omit Seat Sling	NC

JAY Cushions (Optional)			
<i>To order any JAY cushions on-chair, please see the JAY price list.</i>			
<i>Cushion name or part number must be filled out with desired width and depth</i>			
☐	202JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: Depth:	

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

8. Key Options & Accessories

Anti Tips			
☐	202AT1	Anti Tips (Pair)	\$121.00
☐	202AT3	Omit Rear Anti Tip	NC
Side Guards			
☐	202SG33	Plastic - Standard (8.5" High)	\$188.00
☐	202SG29	Plastic - Low (6.5" High)	\$188.00
☐	202SG35	Aluminum <i>N/A w/ height adj arms, push/pull wheel locks only available in black</i>	\$422.00
☐	202SG36	Aluminum w/ Fender <i>N/A w/ height adj arms, push/pull wheel locks only available in black</i>	\$495.00

Pack Accessories			
☐	202A1	Quickie Backpack Black	\$112.00
☐	202A2	Quickie Kids Backpack Black	\$105.00
☐	202A8	Seat Pouch Black	\$83.00

Positioning Belt			
☐	202A141	1.5" Auto Buckle Universal w/clamps	\$97.00
☐	202A120	Positioning Belt Slide Clamps	\$61.00

Spoke Guards			
☐	202A302	Quickie Clear Spoke Guard	\$157.00
☐	202A23	Quickie Black Spoke Guard	\$157.00

Miscellaneous			
☐	202DI07	Neoprene Impact Guards	\$61.00
☐	202A2100	Crutch Holder <i>Only available with Integral or Fold Down push handles</i>	\$115.00
☐	202A85	Bolt-On Push Handles	\$139.00
☐	202A83	Tool Kit	\$43.00

9. COLOR

Standard Frame Colors			
☐	202C211	Celestial Sapphire	NC
☐	202C212	Desert Sage	NC
☐	202C213	Desert Tan	NC
☐	202C214	Gloss Black	NC
☐	202C215	Matte Black	NC
☐	202C216	Metallic Shadow	NC
☐	202C217	Midnight Blue	NC
☐	202C218	Mystic Purple	NC
☐	202C219	Petal Pink	NC
☐	202C220	Pink Power	NC
☐	202C221	Red Pepper	NC
☐	202C222	Royal Blue	NC
☐	202C223	Salute Green	NC
☐	202C224	Scarlet Noir	NC
☐	202C225	Silver Ore	NC
☐	202C226	Solar Orange	NC
☐	202C227	Starlight White	NC
☐	202C228	Storm Gray	NC
☐	202C229	Turbo Gray	NC
☐	202C230	Velocity Yellow	NC

Personalized Frame Color			
☐	202C112	Desert Camouflage	\$368.00
☐	202C43	Mossy Oak Camouflage	\$368.00
☐	202C42	Stars & Stripes	\$368.00
☐	202C108	Zebra (<i>Choose base color</i>) Std w/ Starlight White. <i>N/A w/ Gloss Black and Matte Black.</i>	\$368.00

Additional Options Through Built for Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com



Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8