



Quickie QXi

October 2026

EIF6

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ Prov: _____ PC: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ Prov: _____ PC: _____

Ship To Phone: () _____ Fax: () _____

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

1. FRAME

Base Model

Weight capacity 300 lbs

☐	196WM78	Quickie QXi (Adjustable Axle)	\$2,717.00
--------------------------	---------	-------------------------------	------------

Frame Style

See Seat Height Matrix on page 2

☐	196FS9	Standard Frame	NC
☐	196FS30	Hemi Frame	NC

Transit

☐	196TO38	Transit System WC19	NC
--------------------------	---------	---------------------	----

2. SIZE / SEATING

Frame Width

☐	196W5	14"	NC
☐	196W6	15"	NC
☐	196W7	16"	NC
☐	196W8	17"	NC
☐	196W9	18"	NC
☐	196W10	19"	\$308.00
☐	196W11	20"	\$308.00

Frame Width - Built 4 Me Options

Not available with EXO standard or Tension adjustable upholstery. Not available with Ballistic Nylon or 3DX Tension adjustable upholstery.

☐	196W59	B4ME - 21" Width	\$308.00
☐	196W60	B4ME - 22" Width	\$308.00

Seat Sling

☐	196SL1	Standard Seat Sling	NC
☐	196SL6	Omit Seat Sling	NC

Seat Depth

Frame Depths equals Seat Depths.

☐	196SD7	14"	NC
☐	196SD8	15"	NC
☐	196SD9	16"	NC
☐	196SD10	17"	NC
☐	196SD11	18"	NC
☐	196SD71	19"	\$308.00
☐	196SD13	20"	\$308.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	196JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____ Depth: _____	

Personalized Seating

☐	196QC01	Surefit	See OF
☐	196QC02	Configurfit	See OF
☐	196QC03	JayFit	See OF



Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair



Quickie QXi featuring the Jay ION™ cushion

✓	Option #	Definition	Price	✓	Option #	Definition	Price
3. CASTERS / SEAT HEIGHT							

Front Seat Height							
See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger							
12.5" & 13" only available w/ Alum Angle Adj Footplate							
☐	196SH34	12.5"	NC	☐	196SH14	16.0"	NC
☐	196SH35	13.0"	NC	☐	196SH15	16.5"	NC
☐	196SH40	13.5"	NC	☐	196SH16	17.0"	NC
☐	196SH41	14.0"	NC	☐	196SH17	17.5"	NC
☐	196SH42	14.5"	NC	☐	196SH18	18.0"	NC
☐	196SH12	15.0"	NC	☐	196SH20	18.5"	NC
☐	196SH13	15.5"	NC	☐	196SH21	19.0"	NC
				☐	196SH23	19.5"	NC

Caster Wheels			
☐	196CW3	4" Low Profile Polyurethane	NC
☐	196CW44	4" x 1.25" Semi Pneumatic	NC
☐	196CW5	5" Low Profile Polyurethane	NC
☐	196CW7	5" x 1.5" Semi Pneumatic	NC
☐	196CW34	6" Polyurethane	NC
☐	196CW35	6" x 1.5" Semi Pneumatic	NC
☐	196CW32	8" Polyurethane	NC
☐	196CW39	8" x 1.5" Semi Pneumatic	\$146.00

Caster Fork			
☐	196CF71	3.5" (standard stem bolt only)	NC
☐	196CF12	5.25"	NC
☐	196CF17	6"	NC
☐	196CF23	7"	NC

Caster Stem Bolt			
☐	196CS1	Standard	NC
☐	196CS5	3/4"	NC
☐	196CS9	1 1/2"	NC

Front Seat Height Matrix

Front Seat Height Matrix - Standard Stem Bolt									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0	
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0	
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0	
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0	
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0	
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0	
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	19.5	17.0	18.5	
Seat squeeze and seat depth can affect FSH availability									
For 3/4" and 1 1/2" Stem bolt matrix see page 5									

4. REAR WHEELS / REAR SEAT HEIGHTS

Rear Seat Height			
<i>See Rear Seat Height Matrix</i>			
☐	196RH9	12"	NC
☐	196RH10	12.5"	NC
☐	196RH11	13"	NC
☐	196RH12	13.5"	NC
☐	196RH13	14"	NC
☐	196RH14	14.5"	NC
☐	196RH15	15"	NC
☐	196RH16	15.5"	NC
☐	196RH17	16"	NC
☐	196RH18	16.5"	NC
☐	196RH19	17"	NC
☐	196RH26	17.5"	NC
☐	196RH20	18"	NC
☐	196RH27	18.5"	NC
☐	196RH21	19"	NC

Rear Wheel Size			
<i>See Rear Tire & Wheel Matrix for availability</i>			
☐	196RW3	20" Mag Wheel	NC
☐	196RW8	22" Mag Wheel	NC
☐	196RW1	24" Mag Wheel	NC
☐	196RW103	20" Mag Wheel 5 Spoke	NC
☐	196RW104	22" Mag Wheel 5 Spoke	NC
☐	196RW102	24" Mag Wheel 5 Spoke	NC
☐	196RW90	20" Lite Spoke Wheel	NC
☐	196RW91	22" Lite Spoke Wheel	NC
☐	196RW92	24" Lite Spoke Wheel	NC
☐	196RW93	25" Lite Spoke Wheel	NC
☐	196RW94	26" Lite Spoke Wheel	\$126.00
☐	196RW25	Omit	NC

Rear Seat Height Matrix		
Rear Wheel Size	Std Frame	Hemi Frame
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"
24" Rear Wheel	15.0" - 18.5"	13.5" - 16.5"
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.0"
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"

Rear Wheel Tire			
<i>See Rear Tire & Wheel Matrix for availability</i>			
☐	196RT1	Pneumatic	NC
☐	196RT2	Pneumatic w/ Airless Insert	\$247.00
☐	196RT3	Full Poly	NC
☐	196RT39	High-Pressure Clincher Black	\$105.00
☐	196RT21	Schwalbe Marathon Plus Evolution	\$157.00

Handrim			
<i>Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic N/A w/ 25, 26"</i>			
<i>See Appendix for Rear wheel / Handrim matrix</i>			
☐	196HR1	Aluminum Anodized	NC
☐	196HR4	Plastic	\$146.00
☐	196HR22	Natural Fit - Standard Grip	\$374.00
☐	196HR23	Natural Fit - Super Grip	\$434.00
☐	196HR39	Ellipse 3R	\$550.00
☐	196HR40	Carbolife - Gekko	\$575.00
☐	196HR41	Carbolife - Curve L Tetra Grip	\$625.00
☐	196HR28	Q-Grip Neoprene Coated	\$471.00

Seat Squeeze Matrix	
Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)

Rear Tire & Wheel Matrix

Tire	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Pneumatic	20,22,24	20,22,24,26	20,22,24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24
Full Poly	20,22,24	20,22,24	20,22,24
High-Pressure Clincher Blk	N/A	22,24,26	N/A
Schwalbe Marathon	N/A	24,25,26	N/A
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm			

One Arm Drive

Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim.

1AD option is only available on 14" - 19" frame widths.

☐ 196AL16	One Arm Drive Adj. Axle Plate	\$1,382.00
☐ 196HM1	☐ Right Hand Mount	NC
☐ 196HM2	☐ Left Hand Mount	NC

Axle

☐ 196AX1	Quick Release Steel	NC
☐ 196AX11	Quad Quick Release Steel N/A on Spoke	\$182.00
☐ 196AX45	Threaded Axle	\$150.00

Adjustable Axle Options

☐ 196AS4	Adjustable sleeve	NC
☐ 196AS5	Fixed Sleeve 0°	NC
☐ 196AS6	Fixed Sleeve 3°	NC

Center of Gravity

Measure from front of back post to center of rear axle

☐ 196CC2	1"	NC
☐ 196CC3	1.5"	NC
☐ 196CC4	2"	NC
☐ 196CC5	2.5"	NC
☐ 196CC6	3"	NC

5. WHEEL LOCKS

Wheel Lock

☐ 196WL1	High Mount Push	NC
☐ 196WL2	High Mount Pull	NC
☐ 196WL30	Compact	\$162.00
☐ 196WL24	Unilateral Wheel Lock Right	\$320.00
☐ 196WL23	Unilateral Wheel Lock Left	\$320.00
☐ 196WL10	Omit Wheel Lock	NC
☐ 196WL12	Do Not Mount	NC

Extension Handle

Ext. handles only available with high mount push or pull wheel locks

☐ 196EH7	Extension Handle (Pair)	\$73.00
---------------------------------	-------------------------	---------

6. HANGERS AND FOOTRESTS

Front Hangers

Footplate mounts with Extension Tube using clamping hardware

☐ 196H145	60° Swing In/Out Hanger	NC
☐ 196H139	70° Swing In/Out Hanger	NC
☐ 196H141	80° Swing In/Out Hanger	NC
☐ 196H159	Articulating Hanger N/A w/ Alum Angle Adj Front Mount footplates.	\$526.00
☐ 196H160	Elevating Hanger - Ext Mount N/A w/ Alum Angle Adj Front Mount footplates.	\$339.00
☐ 196H161	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$339.00
☐ 196H144	Omit Hanger & Footplate	NC

Articulating / Elevating Hanger Cover

☐ 196DI47	Neoprene Cover	\$325.00
☐ 196DI48	Gel Pad Cover	\$380.00

Extension Tubes

See Lower Leg Length Matrix for range

☐ 196E1	Short	NC
☐ 196E2	Medium	NC
☐ 196E5	Long	NC

Footplate

Heel loop not included, see Footplate Option section below

☐ 196F1	Composite Footplate	NC
☐ 196F271	Composite Angle Adj-(16"-22"w)	\$176.00
☐ 196F3	Aluminum Angle Adj. Footplate	\$176.00
☐ 196F5	Aluminum Only available w/Elevating Legrest	NC

Hanger & Footplate Options

☐ 196HL12	Adjustable Heel Loop	NC
☐ 1968A02	Gel Cover For Hangers (Pair)	\$162.00

7. BACKREST / SEATING

Back Type

☐ 196BT10	Standard	NC
☐ 196BT100	Angle Adjustable Back	\$357.00

Backrest Height

☐ 196B5	Med (13"-17")	NC
☐ 196B9	Tall (16"-20")	NC

Push Handles

Integral Push Handles, Folding Height Adjustable & Omit is 8° for Standard

Backrest, Straight for Angle Adj. Backrest

☐ 196PH6	Integral Push Handles	NC
☐ 196PH5	Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height	\$168.00
☐ 196PH17	Folding Push Handle with Ht. Adj. Back	\$310.00
☐ 196PH2	Omit Push Handles	NC

✓	Option #	Definition	Price	✓	Option #	Definition	Price
7. BACKREST / SEATING (CONT)							

Back Upholstery			
☐	196BU1	Standard Upholstery	NC
☐	196BU88	Standard Upholstery EXO	\$75.00
☐	196BU70	Backrest Tension Adj Ballistic Nylon <i>*Reference 186BU70 for billing</i>	\$190.00
☐	196BU71	Backrest Tension Adj 3DX Vented <i>*Reference 186BU71 for billing</i>	\$471.00
☐	196BU89	Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i>	\$449.00
☐	196BU14	Omit Back Upholstery	NC

JAY Backs			
☐	196JB01	JAY Back (see Order forms)	Price List

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away			
<i>Some limitations when used with the Angle Adjustable Back in the rear most setting</i>			
<i>Please circle a setting, will default to 11" if no selection is made</i>			
☐	196AA18	Padded Swing-Away STD (9" to 13")	\$278.00
Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"			

Armrest - Height Adjustable			
<i>Flip Back HA Classic N/A w / 14"-15" depth</i>			
☐	196AA104	Single Post Ht. Adj - Std (9.75"-14.75")	\$339.00
☐	196AA105	Single Post Ht Adj - Low (8.25"-12.25")	\$339.00
☐	196AA65	Flip Back Ht Adj LITE (10"-15")	\$381.00
☐	196AA103	Dual Post Flip Back Ht Adj	\$339.00
☐	196AA159	Flip Back Ht Adj Classic (11"-15")	\$381.00
☐	196AA14	Omit Armrest	NC

Armpad - For Height Adjustable			
☐	196AM23	Desk Length (10")	NC
☐	196AM24	Full Length - Full (14")	NC

Armrest Options			
☐	1968A01	Gel Cover For Swing Away Arms (pair)	\$162.00

Side Guards			
☐	196SG33	Plastic - Standard (8.5" High)	\$188.00
☐	196SG29	Plastic - Low (6.5" High)	\$188.00

9. KEY OPTIONS

Anti Tips			
☐	196AT1	Anti Tips (Pair)	\$121.00
☐	196AT86	Anti Tips (Pair) Low Profile	\$121.00

Positioning Belt			
☐	196A141	1.5" Auto Buckle w/ Clamps	\$97.00
☐	196A804	Positioning Belt Slide Clamps	\$77.00

JAY Positioning			
☐	196K1	JAY Lower Extremity S/A support ☐ Left ☐ Right Position	\$320.00
☐	196K2	Omit Hanger ☐ Left ☐ Right	NC
☐	196K7	Pad Size ☐ 6" x 8" ☐ 8" x 10"	NC

Pack Accessories			
☐	196A1	Backpack - Black/Charcoal Adult	\$112.00
☐	1961A01	Seat Pouch Black	\$83.00

Spoke Guards			
☐	1964A01	Quickie Black Spoke Guard	\$157.00
☐	1964A02	Quickie Clear Spoke Guard	\$157.00

Miscellaneous			
☐	1968A03	Caddy	\$89.00
☐	196A2100	Crutch Holder	\$115.00
☐	196A999	Quickie Water Bottle & Holder	\$81.00

10. COLOR

Standard Frame Colors			
☐	196C211	Celestial Sapphire	NC
☐	196C212	Desert Sage	NC
☐	196C213	Desert Tan	NC
☐	196C214	Gloss Black	NC
☐	196C215	Matte Black	NC
☐	196C216	Metallic Shadow	NC
☐	196C217	Midnight Blue	NC
☐	196C218	Mystic Purple	NC
☐	196C219	Petal Pink	NC
☐	196C220	Pink Power	NC
☐	196C221	Red Pepper	NC
☐	196C222	Royal Blue	NC
☐	196C223	Salute Green	NC
☐	196C224	Scarlet Noir	NC
☐	196C225	Silver Ore	NC
☐	196C226	Solar Orange	NC
☐	196C227	Starlight White	NC
☐	196C228	Storm Gray	NC
☐	196C229	Turbo Green	NC
☐	196C230	Velocity Yellow	NC

Personalized Frame Color			
☐	196C112	Desert Camouflage	\$368.00
☐	196C43	Mossy Oak Camouflage	\$368.00
☐	196C42	Stars & Stripes	\$368.00
☐	196C108	Zebra (Choose base color) N/A w/ Gloss Black and Matte Black.	\$368.00

Additional Options Through Built for Me
 To inquire about any unique personalized Quickie products, contact: Phone: (800) 263-3390, Email: cscanada@sunmed.com



APPENDIX

Front Seat Height Matrix Cont.

Front Seat Height Matrix - 3/4" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5	
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5	
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5	
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0	
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5	
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5	
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5	
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5	

Seat squeeze and seat depth can affect FSH availability


Front Seat Height Matrix - 1 1/2" Stem									
Caster	3.5" FORK				5 1/4" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5	
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0	
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5	
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Caster	6" FORK				7" FORK				
	STD		HEMI		STD		HEMI		
	Min	Max	Min	Max	Min	Max	Min	Max	
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5	
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0	
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
6 x 1.5 Semi Pneu	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0	
8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5	
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5	

Seat squeeze and seat depth can affect FSH availability

Lower Leg Length Range Matrix

Extension Mount Hanger Leg Lengths						
Hanger	Ext Tube	Seat-to-Footrest Range				
		Composite	Comp Angle Adj	Alum Angle Adj	Alum	
60° Swing In/Out	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A	
	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A	
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A	
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A	
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A	
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A	
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A	
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A	
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A	

All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).



Min Allowable Front Seat Heights (for Extension Mount Hangers)					
Hanger	Ext Tube	Min Allowable FSH			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	14.0	13.5	12.5	N/A
	Medium	14.5	14.5	12.5	N/A
	Long	16.0	16.0	13.5	N/A
70° Swing In/Out	Short	14.0	14.0	13.0	N/A
	Medium	15.0	15.0	13.5	N/A
	Long	16.5	16.5	15.0	N/A
80° Swing In/Out	Short	14.0	14.0	13.5	N/A
	Medium	14.5	14.5	14.5	N/A
	Long	16.5	16.5	15.5	N/A

ALR & ELR Hanger Leg Lengths				
Hanger	Extension	Seat Sling to Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	Short	8	11	10
	Med	10	13	12
	Long	12	15	14

Handrim	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Aluminum Anodized	20,22,24	20,22,24,25,26	20,22,24
Plastic Coat	20,22,24	22,24	20,22,24
Natural Fit - Standard Grip	24	22,24,25,26	22,24
Natural Fit - No Thumb Grip	24	22,24,25,26	22,24
Natural Fit - Super Grip	24	22,24,25,26	22,24
Surge - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Surge LT - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Q Grip Neoprene Coated	N/A	24,25	N/A



Sunrise Medical Canada Inc.

1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8