



ORDER FORM

# Horizon Stander

# EILK02

Robust and durable 3-in-1 standing frame for both young children and adult users.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

## Horizon Stander

Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	RETAIL
<b>Frame</b>		<b>\$</b>
Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals. (Knee supports and Footplates separate)		
<input type="checkbox"/> 141FO5	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	\$2,453
<input type="checkbox"/> 141FO6	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	\$2,798
<b>Base</b>		<b>\$</b>
Selection Required. MUST match Frame size. Includes worldwide adapters.		
<input type="checkbox"/> 141NB110	Size 2 Horizon Base/Chassis	\$4,827
<input type="checkbox"/> 141NB111	Size 3 Horizon Base/Chassis	\$6,026
<b>Headrest</b>		<b>\$</b>
Head Support Assembly selection required if Head Lateral Support/Extension Bracket is selected.		
<input type="checkbox"/> 141HE3	Size 2 and 3 Head Support Assembly - Black	\$394
<input type="checkbox"/> 141HE4	Size 2 and 3 Head Lateral Supports - Black (Pair)	\$339
<input type="checkbox"/> 141JG12	Headrest 4" Horizontal Extension Bracket	\$527
<input type="checkbox"/> 141Y029	Size 1-3 Contoured Head Support	\$861
<b>Hip Laterals</b>		<b>\$</b>
Selection required if Hip Belt is not selected.		
<input type="checkbox"/> 141KC7	Size 2 Horizon Hip Laterals - Black (Pair)	\$504
<input type="checkbox"/> 141KC8	Size 3 Horizon Hip Laterals - Black (Pair)	\$524
<b>Hip Belts</b>		<b>\$</b>
Selection required if Hip Laterals are not selected. Multiple Selections can be made.		
<input type="checkbox"/> 141A856	Size 2 Supine Belt - Black	\$302
<input type="checkbox"/> 141A857	Size 3 Supine Belt - Black	\$332
<b>Knee Supports</b>		<b>\$</b>
Selection of Standard Knee Assembly Required.		
Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" adjustment.		
Long Knee Bracket has 7" - 9" adjustment.		
<input type="checkbox"/> 141XX16	Size 2 and 3 Standard Knee Assembly - Black (Pair)	\$602
<input type="checkbox"/> 141JG14	Medium Knee Adjustment Bracket (Pair)	\$107
<input type="checkbox"/> 141JG15	Long Knee Adjustment Bracket (Pair)	\$111
<b>Lower Leg Positioning</b>		<b>\$</b>
Selection of Footplate Required.		
<input type="checkbox"/> 1418A107	Small Sandals with Straps (Pair)	\$273
<input type="checkbox"/> 1418A108	Medium Sandals with Straps (Pair)	\$291
<input type="checkbox"/> 1418A109	Large Sandals with Straps (Pair)	\$315
<input type="checkbox"/> 141F53	Size 2 and 3 Flip Up Footplate Assembly (Pair)	\$906
<input type="checkbox"/> 141F65	Size 1-3 3D Footplate Assembly (Pair) - One size fits all	\$1,821
<b>Upper Extremity Support Tray</b>		<b>\$</b>
Must match Horizon Stander Size		
<input type="checkbox"/> 141A209	Size 2 and 3 Supine Tray	\$434
<input type="checkbox"/> 141A210	Size 2 and 3 Prone Tray	\$783

Technical Specifications - Imperial		
Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

Technical Specifications - Metric		
Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm

Sunrise Medical Canada Inc.  
 1000 Creditstone Road, Unit 2  
 Concord, ON L4K 4P8, Canada



MK-129848-CAN Rev. 10 © 04/2025

