



ORDER FORM

Horizon Stander

EILK02

Robust and durable 3-in-1 standing frame for both young children and adult users.

| Customer | | | |
|--------------------------------|--------------------------------|------|--|
| Account | | | |
| Account # | | Date | |
| <input type="checkbox"/> Quote | <input type="checkbox"/> Order | PO # | |
| Buyer | | | |
| Mark For | | | |

| Shipping Address | | | |
|------------------|--|-------|--|
| Name | | | |
| Address | | | |
| City | | Prov | |
| PC | | Phone | |
| Email / Fax | | | |

Horizon Stander

Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 220 lbs

| ✓ OPTION | DESCRIPTION | RETAIL |
|---|--|---------|
| Frame | | \$ |
| Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals. (Knee supports and Footplates separate) | | |
| <input type="checkbox"/> 141FO5 | Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey | \$2,453 |
| <input type="checkbox"/> 141FO6 | Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey | \$2,798 |
| <input type="checkbox"/> 141FO7 | Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Vinyl | \$2,453 |
| <input type="checkbox"/> 141FO8 | Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Vinyl | \$2,798 |
| Base | | \$ |
| Selection Required. MUST match Frame size. Includes worldwide adapters. | | |
| <input type="checkbox"/> 141NB110 | Size 2 Horizon Base/Chassis | \$4,827 |
| <input type="checkbox"/> 141NB111 | Size 3 Horizon Base/Chassis | \$6,026 |
| Headrest | | \$ |
| Head Support Assembly selection required if Head Lateral Support/Extension Bracket is selected. | | |
| <input type="checkbox"/> 141HE3 | Size 2 and 3 Head Support Assembly - Black | \$394 |
| <input type="checkbox"/> 141HE4 | Size 2 and 3 Head Lateral Supports - Black (Pair) | \$339 |
| <input type="checkbox"/> 141HE63 | Size 2 and 3 Head Support Assembly - Vinyl | \$394 |
| <input type="checkbox"/> 141HE64 | Size 2 and 3 Head Lateral Supports - Vinyl (Pair) | \$339 |
| <input type="checkbox"/> 141JG12 | Headrest 4" Horizontal Extension Bracket | \$527 |
| <input type="checkbox"/> 141Y029 | Size 1-3 Contoured Head Support | \$861 |
| Hip Laterals | | \$ |
| Selection required if Hip Belt is not selected. | | |
| <input type="checkbox"/> 141KC7 | Size 2 Horizon Hip Laterals - Black (Pair) | \$504 |
| <input type="checkbox"/> 141KC8 | Size 3 Horizon Hip Laterals - Black (Pair) | \$524 |
| <input type="checkbox"/> 141KC29 | Size 2 Horizon Hip Laterals - Vinyl (Pair) | \$504 |
| <input type="checkbox"/> 141KC30 | Size 3 Horizon Hip Laterals - Vinyl (Pair) | \$524 |
| Hip Belts | | \$ |
| Selection required if Hip Laterals are not selected. Multiple Selections can be made. | | |
| <input type="checkbox"/> 141A856 | Size 2 Supine Belt - Black | \$302 |
| <input type="checkbox"/> 141A857 | Size 3 Supine Belt - Black | \$332 |
| <input type="checkbox"/> 141A883 | Size 2 Supine Belt - Vinyl | \$302 |
| <input type="checkbox"/> 141A884 | Size 2 Supine Belt - Vinyl | \$332 |
| Knee Supports | | \$ |
| Selection of Standard Knee Assembly Required. | | |
| Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" adjustment. | | |
| Long Knee Bracket has 7" - 9" adjustment. | | |
| <input type="checkbox"/> 141XX16 | Size 2 and 3 Standard Knee Assembly - Black (Pair) | \$602 |
| <input type="checkbox"/> 141XX17 | Size 2 and 3 Standard Knee Assembly - Vinyl (Pair) | \$602 |
| <input type="checkbox"/> 141JG14 | Medium Knee Adjustment Bracket (Pair) | \$107 |
| <input type="checkbox"/> 141JG15 | Long Knee Adjustment Bracket (Pair) | \$111 |
| Lower Leg Positioning | | \$ |
| Selection of Footplate Required. | | |
| <input type="checkbox"/> 1418A107 | Small Sandals with Straps (Pair) | \$273 |
| <input type="checkbox"/> 1418A108 | Medium Sandals with Straps (Pair) | \$291 |
| <input type="checkbox"/> 1418A109 | Large Sandals with Straps (Pair) | \$315 |
| <input type="checkbox"/> 141F53 | Size 2 and 3 Flip Up Footplate Assembly (Pair) | \$906 |
| <input type="checkbox"/> 141F65 | Size 1-3 3D Footplate Assembly (Pair) - One size fits all | \$1,821 |
| Upper Extremity Support Tray | | \$ |
| Must match Horizon Stander Size | | |
| <input type="checkbox"/> 141A209 | Size 2 and 3 Supine Tray | \$434 |
| <input type="checkbox"/> 141A210 | Size 2 and 3 Prone Tray | \$783 |

Technical Specifications - Imperial

| Frame Sizing Specification | Size 2 | Size 3 |
|---|-----------------------|-----------------------|
| Age Range | 9 -- 18 | 14 -- Adult |
| Horizontal Position L x W x H | 50.4" x 25.2" x 31.7" | 50.4" x 25.2" x 31.7" |
| Weight of Product | 110.6 lbs | 130 lbs |
| Maximum User Weight | 176 lbs | 220 lbs |
| User Standing Height (Min / Max) | 55.1" x 70.9" | 61" / 70.9" |
| Shoulder Height - Footplate to Top of Pad (Min / Max) | 39.4" / 57" | 39.4" / 57" |
| Hip Pad Height - Footplate to Top of Pad (Min / Max) | 29.5" / 37.4" | 29.5" / 45.3" |
| Hip Pad Width (Min / Max) | 8.3" / 15.3" | 12.2" / 18.5" |
| Knee Height - Footplate to Middle of Pad (Min / Max) | 15.7" / 21.6" | 15.7" / 25.6" |
| Width Between Chest Laterals (Min / Max) | 9.4" / 14.9" | 13.4" / 18.1" |
| Footplate Angle Adjustable + / - | 10° | 10° |
| Tray Sizes | 20.5" x 23.2" | 20.5" x 23.2" |
| Tray Angle Adjustment - Prone Version | 0° to 40° | 0° to 40° |
| Tray Angle Adjustment - Supine Version | -20° to +30° | -20° to +30° |
| Frame Clearance - Bottom of Ski to Ground | 7" | 7" |

Technical Specifications - Metric

| Frame Sizing Specification | Size 2 | Size 3 |
|---|---------------------------|---------------------------|
| Age Range | 9 -- 18 | 14 -- Adult |
| Horizontal Position L x W x H | 1280 mm x 640 mm x 805 mm | 1280 mm x 640 mm x 805 mm |
| Weight of Product | 50 kgs | 95 kgs |
| Maximum User Weight | 80 kgs | 100 kgs |
| User Standing Height (Min / Max) | 1400 mm / 1800 mm | 1550 mm / 1800 mm |
| Shoulder Height - Footplate to Top of Pad (Min / Max) | 1000 mm / 1450 mm | 1000 mm / 1450 mm |
| Hip Pad Height - Footplate to Top of Pad (Min / Max) | 750 mm / 950 mm | 750 mm / 1150 mm |
| Hip Pad Width (Min / Max) | 210 mm / 390 mm | 310 mm / 470 mm |
| Knee Height - Footplate to Middle of Pad (Min / Max) | 400 mm / 550 mm | 400 mm / 650 mm |
| Width Between Chest Laterals (Min / Max) | 240 mm / 380 mm | 340 mm / 460 mm |
| Footplate Angle Adjustable + / - | 10° | 10° |
| Tray Sizes | 520 mm x 590 mm | 520 mm x 590 mm |
| Tray Angle Adjustment - Prone Version | 0° to 40° | 0° to 40° |
| Tray Angle Adjustment - Supine Version | -20° to +30° | -20° to +30° |
| Frame Clearance - Bottom of Ski to Ground | 180 mm | 180 mm |

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, ON L4K 4P8, Canada



MK-129848-CAN Rev. 11 © 10/2025

