

# LECKEY<sup>®</sup> Pal

## Classroom Seat





# Classroom Seat




## Posture + Attention = Learning

Designed by the Leckey team of experienced therapists and engineers, Pal is an attractive, robust and easy to use classroom and nursery **seating system** for kids of **1-12 years** of age, with **mild to moderate postural needs**.

It is designed to give **improved stability**, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. Available in **4 sizes** and **4 colours** it will help kids up to 12 years of age to participate, learn and give their full attention in the classroom.

Available in 4 sizes and 4 colours





**Humpty Dumpty**  
Humpty Dumpty sat on a wall,  
Humpty Dumpty had a great fall.  
All the king's horses  
and all the king's men  
Couldn't put Humpty together a



**Little Miss Muffet**  
Little Miss Muffet  
Sat on a tuffet  
Eating her  
Along c  
Who s  
And f

Red

Brown

Orange

Yellow

Green

Blue

White

Purple

i

he

come

we

she

there

be



## Created for self assessment

The simplicity of Pal enables the therapist to assess the child's needs, instantly order their **Pal Classroom Seat** and quickly and easily set up the product when it arrives. Pal really is the simple, fast way to help kids join in with learning.



Pal Seat shown with  
footrest accessory



Pal Seat shown with  
castor leg accessory




Pal Seat with  
rocker leg accessory



Pal Seat with  
stabiliser leg  
accessory





## Pal Sizing Chart

	Age Range	Stature Range	User Weight	Seat Depth	Seat Width	Seat Height	Armrest Height	Backrest Height	Seat Weight	Table Height
<b>Size 1</b>	1 - 3 years	700 - 950mm	Max - 30kg	180 - 240mm	150 - 250mm	160 - 240mm	120-150mm	265mm	4.3kg	400mm
<b>Size 2</b>	2 - 5 years	900 - 1050mm	Max - 35kg	220 - 290mm	170 - 270mm	185 - 285mm	140-170mm	315mm	6kg	460mm
<b>Size 3</b>	4 - 8 years	1000 - 1300mm	Max - 40kg	250 - 330mm	230 - 330mm	250 - 350mm	180-230mm	380mm	8.2kg	530mm
<b>Size 4</b>	7-12 years	1200 - 1600mm	Max - 50kg	300 - 400mm	280 - 380mm	300 - 400mm	190-240mm	450mm	10kg	n/a

**Backrest Angle** 95 - 105°

The Pal Seat comes as standard with depth adjustable seat base, angle adjustable back rest, 2 point pelvic strap, height adjustable legs, height adjustable armrests and wipe clean back rest and seat base cushions.

## Pal Product Codes

X = Size 1, 2, 3 or 4

AA = Colour (01 Green, 02 Orange, 03 Blue, 04 Pink)



**Seat**  
135-X600-AA



**Ski with Castors (Pair)**  
135-X610



**Rocker (Pair)**  
135-X604



**Stabilisers (Pair)**  
135-X604



**Ramped Seat Cushion**  
135-X651-AA



**Footrest**  
135-X603



**25mm Side Pad (Pair)**  
135-X635-AA  
**50mm Side Pad (Pair)**  
135-X654-AA



**Depth Adjustable Activity Tray**  
135-X620  
**Tray Pad (Not shown)**  
135-X656



**Pommel**  
Size 1 - 135-1602-AA  
Size 2 & 3 - 135-3602-AA

24 hour postural  
care for babies,  
kids & adults.  
Sleeping, Sitting,  
Standing, Walking,  
Moving, Bathing,  
Toileting.