



Quickie M6™ - AADL W848

October 2024 EIF1

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ PR: _____ PC: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote ☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ PR: _____ PC: _____
Ship To Phone: () _____ Fax: () _____

✓ Option # Definition Price ✓ Option # Definition Price

1. FRAME

Base Model

☐ 115WM21 Quickie M6 \$4,740.00

2. SIZE / SEATING

Frame Width

| | | | |
|--------------------------|--------|-----|----|
| <input type="checkbox"/> | 115W13 | 22" | NC |
| <input type="checkbox"/> | 115W14 | 24" | NC |
| <input type="checkbox"/> | 115W15 | 26" | NC |
| <input type="checkbox"/> | 115W16 | 28" | NC |
| <input type="checkbox"/> | 115W17 | 30" | NC |

Frame Length

| | | | |
|--------------------------|---------|-------------------------|----|
| <input type="checkbox"/> | 115FL9 | Regular (18"-19" depth) | NC |
| <input type="checkbox"/> | 115FL18 | Long (18"-22" depth) | NC |

Seat Depth

| | | | |
|--------------------------|---------|-----------------------|----------|
| <input type="checkbox"/> | 115SD11 | 18" | NC |
| <input type="checkbox"/> | 115SD12 | 19" | NC |
| <input type="checkbox"/> | 115SD13 | 20" (Long frame only) | \$303.00 |
| <input type="checkbox"/> | 115SD14 | 21" (Long frame only) | \$303.00 |
| <input type="checkbox"/> | 115SD33 | 22" (Long frame only) | \$303.00 |

Seat Height

| | | | |
|--------------------------|---------|--|----------|
| <input type="checkbox"/> | 115HT65 | 17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster | \$552.00 |
| <input type="checkbox"/> | 115HT11 | 18.5" Only available w/ 5" casters | NC |
| <input type="checkbox"/> | 115HT14 | 20" | NC |

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out with desired width and depth

| | | | |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 115JC20 | JAY Bundled Cushion Cushion Name or Part # Width: _____ Depth: _____ | Price List |
|--------------------------|---------|--|------------|

3. CASTERS

Casters

| | | | |
|--------------------------|---------|------------------------------------|----------|
| <input type="checkbox"/> | 115CW6 | 5" Soft Roll | \$124.00 |
| <input type="checkbox"/> | 115CW36 | 8" x 2" Pneumatic | NC |
| <input type="checkbox"/> | 115CW37 | 8" x 2" Pneumatic w/Airless Insert | \$387.00 |

Seat Height Matrix

| Caster | Front/Rear |
|--------|-----------------------------|
| 5" | 17" (Different Axle & Fork) |
| 5" | 18.5" |
| 8" | 20" |

4. REAR WHEELS

Rear Wheel Type

| | | | |
|--------------------------|---------|--------------------|------------|
| <input type="checkbox"/> | 115RW2 | 24" Spoke Wheel | NC |
| <input type="checkbox"/> | 115RW29 | 24" Spinergy Wheel | \$1,850.00 |

Rear Wheel Tire

| | | | |
|--------------------------|---------|-------------------------------|----------|
| <input type="checkbox"/> | 115RT1 | Pneumatic Tire | NC |
| <input type="checkbox"/> | 115RT2 | Pneumatic Airless Insert Tire | \$243.00 |
| <input type="checkbox"/> | 115RT25 | Primo Orion Poly | \$267.00 |

Handrim

| | | | |
|--------------------------|--------|-------------------|----------|
| <input type="checkbox"/> | 115HR1 | Aluminum Anodized | NC |
| <input type="checkbox"/> | 115HR4 | Plastic Coated | \$144.00 |
| <input type="checkbox"/> | 115HR7 | Omit Handrim | NC |

Adjustable Axle Plate

| | | | |
|--------------------------|---------|------------------------------------|----|
| <input type="checkbox"/> | 115AL2 | Standard (for 18.5", 20" Seat Ht.) | NC |
| <input type="checkbox"/> | 115AL28 | Low (only for 17" Seat Ht.) | NC |

Axle

| | | | |
|--------------------------|--------|---------------|----|
| <input type="checkbox"/> | 115AX1 | Quick Release | NC |
|--------------------------|--------|---------------|----|

Rear Tire & Wheel Matrix

| Tire | Rear Wheel | |
|-----------------------------|-------------|----------------|
| | Spoke wheel | Spinergy wheel |
| Pneumatic | 24 | 24 |
| Pneumatic w/ Airless Insert | 24 | 24 |
| Primo Orion Poly | 24 | 24 |

| 5. WHEEL LOCKS | | | | | | | |
|--------------------------|---------|----------------------------|----|--------------------------|--------|---------|---------|
| Wheel Lock | | | | Extension Handle | | | |
| <input type="checkbox"/> | 115WL1 | High Mount Push Wheel Lock | NC | <input type="checkbox"/> | 115EH1 | 6" Pair | \$72.00 |
| <input type="checkbox"/> | 115WL12 | Do Not Mount Wheel Lock | NC | | | | |
| <input type="checkbox"/> | 115WL10 | Omit Wheel Lock | NC | | | | |

| 6. HANGERS AND FOOTREST | | | | | | | |
|--------------------------|---------|---|------------|--------------------------|-------|----------------------------------|----|
| Front Hangers | | | | Extension Tubes | | | |
| <input type="checkbox"/> | 115H1 | 60° Swing Away | NC | <input type="checkbox"/> | 115E1 | Short (13"-15.5") | NC |
| <input type="checkbox"/> | 115H146 | 70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i> | \$190.00 | <input type="checkbox"/> | 115E2 | Medium (15.5"-18") | NC |
| <input type="checkbox"/> | 115H132 | 80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i> | \$185.00 | <input type="checkbox"/> | 115E5 | Long (18"-20.5") NA w/ 5" caster | NC |
| <input type="checkbox"/> | 115H41 | ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i> | \$1,028.00 | | | | |

| Footplates | | | |
|---|--------|-----------------|----------|
| <i>Heel loop or leg strap not included, See Footplate Options Section</i> | | | |
| <input type="checkbox"/> | 115F1 | Composite | NC |
| <input type="checkbox"/> | 115F25 | Adult Angle Adj | \$166.00 |
| <input type="checkbox"/> | 115F5 | Aluminum | NC |
| Footplate Options | | | |
| <input type="checkbox"/> | 115HL1 | Heel Loop | NC |

| Lower Leg Length Range Matrix | | | | |
|--|-------------|-----------------|-----------------|-----------------|
| Hanger | Ext Tube | Footplate | | |
| | | Composite | Alum Angle Adj | Aluminum |
| 60° | Short | 11.50" - 13.25" | 11.75" - 14.00" | 11.25" - 13.50" |
| | Medium | 13.75" - 16.00" | 14.00" - 16.25" | 13.50" - 14.00" |
| | Long | 16.00" - 16.50" | 16.00" - 17.75" | NA |
| 70° | No Ext Tube | NA | 9.75" - 13.50" | NA |
| 80° | No Ext Tube | NA | 10.75" - 14.75" | NA |
| ELR | Short | NA | NA | 13.00" - 17.00" |
| | Medium | NA | NA | 14.25" - 17.00" |
| <i>All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle</i> | | | | |

| 7. BACKREST | | | | | | | |
|--------------------------|---------|--|----|--------------------------|--------|---------|----|
| Backrest | | | | Back Height Preset | | | |
| <input type="checkbox"/> | 115BT34 | 8" w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i> | NC | <input type="checkbox"/> | 115B91 | 16-3/4" | NC |
| | | | | <input type="checkbox"/> | 115B96 | 17-7/8" | NC |
| | | | | <input type="checkbox"/> | 115B39 | 19" | NC |
| | | | | <input type="checkbox"/> | 115B98 | 20-1/8" | NC |

| 8. ARMREST | | | |
|--------------------------|---------|---------------------------------|----------|
| Armrest | | | |
| <input type="checkbox"/> | 115AA33 | Dual Post Height Adj - Desk 10" | \$375.00 |
| <input type="checkbox"/> | 115AA34 | Dual Post Height Adj - Full 14" | \$375.00 |

| 9. KEY OPTIONS & ACCESSORIES | | | | | | | |
|------------------------------|--------|--------------------|----------|--------------------------|---------|-------------------------------------|----------|
| Rear Anti-Tip | | | | Positioning Belt | | | |
| <input type="checkbox"/> | 115AT1 | Rear Anti-Tip | \$119.00 | <input type="checkbox"/> | 115A115 | Belt 2" Aircraft Buckle Long | \$137.00 |
| <input type="checkbox"/> | 115AT3 | Omit Rear Anti-Tip | NC | <input type="checkbox"/> | 115A116 | Belt 2" Aircraft Buckle Padded Long | \$201.00 |
| Miscellaneous | | | | <input type="checkbox"/> | 115A142 | Belt 1.5" Auto Buckle | \$96.00 |
| <input type="checkbox"/> | 115A83 | Tool Kit | \$42.00 | | | | |

10. COLOR

| Standard Frame Colors | | | |
|--------------------------|---------|--------------------|----|
| <input type="checkbox"/> | 115C211 | Celestial Sapphire | NC |
| <input type="checkbox"/> | 115C212 | Desert Sage | NC |
| <input type="checkbox"/> | 115C213 | Desert Tan | NC |
| <input type="checkbox"/> | 115C214 | Gloss Black | NC |
| <input type="checkbox"/> | 115C215 | Matte Black | NC |
| <input type="checkbox"/> | 115C216 | Metallic Shadow | NC |
| <input type="checkbox"/> | 115C217 | Midnight Blue | NC |
| <input type="checkbox"/> | 115C218 | Mystic Purple | NC |
| <input type="checkbox"/> | 115C219 | Petal Pink | NC |
| <input type="checkbox"/> | 115C220 | Pink Power | NC |
| <input type="checkbox"/> | 115C221 | Red Pepper | NC |
| <input type="checkbox"/> | 115C222 | Royal Blue | NC |
| <input type="checkbox"/> | 115C223 | Salute Green | NC |
| <input type="checkbox"/> | 115C224 | Scarlet Noir | NC |
| <input type="checkbox"/> | 115C225 | Silver Ore | NC |
| <input type="checkbox"/> | 115C226 | Solar Orange | NC |
| <input type="checkbox"/> | 115C227 | Starlight White | NC |
| <input type="checkbox"/> | 115C228 | Storm Gray | NC |
| <input type="checkbox"/> | 115C229 | Turbo Green | NC |
| <input type="checkbox"/> | 115C230 | Velocity Yellow | NC |

