SUNRISE MEDICAL

QUICKIE®
Live without limits.

ORDER FORM




Q500M with SEDEO® PRO™ and JAY® Syncra™ Seating

EIPW30

Customer		Ship To		
Account #		Name		
Date		Address		
PO#		City		
Buyer		Province		Postal Code
Mark For		Phone		

Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

Chair Specifications

Weights starting at: (with Batteries)	NPC: 300 lbs SPC: 324 lbs MPC: 338 lbs
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction / 12" seat elevator
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 24", Length 33"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 13.5 miles
Drive Motors	4-pole 6 mph / 10 kph
Drive Wheels	13" w/flat-free inserts (black)
Casters	7" front / 7" rear (black)
Charger	8A Dual Mode (Compact)
Battery Size	Group 34 / Group 22
Colour	8 shroud colours and 6 colour accents
Ground Clearance	3"

BASE

OPTION		DESCRIPTION	RETAIL
Standard Duty Mid Wheel Drive Base *End Item Number EIPW30 (PRO™ Seat)*			$
☐ 300PC215		Quickie Q500 M NPC No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,398.00
☐ 300PC216		Quickie Q500 M SPC Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,398.00
☐ 300PC217		Quickie Q500 M MPC Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,398.00
User Weight			$
☐ 300WG23		Up to 300 lbs.	NC
Motor Package			$
☐ 300PF99		4-pole 10 kph / 6 mph	NC
Tracking			$
☐ 300MI01		SureTrac® Module & Motor Encoders *(Improves drive tracking)*	$1,696.00
Drive Wheels			$
☐ 300DW69		Drive Wheel, 13", Dual Density, Knobby, Black	NC
☐ 300DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			$
☐ 300CW59		Caster Wheels 7", Solid, Black	NC
Shroud Colour			$
☐ 300C119		Matte Black	NC
☐ 300C165		Candy Apple Red	NC
☐ 300C166		Candy Blue	NC
☐ 300C167		Candy Purple	NC
☐ 300C176		Green Apple	NC
☐ 300C178		Pearl Pink	NC
☐ 300C182		Silver	NC
Colour Accents *Note: Accent colours and shroud colours may not match.*			$
☐ 300EZ20		Apple Green	NC
☐ 300EZ22		Black	NC
☐ 300EZ21		Blue	NC
☐ 300EZ23		Orange	NC
☐ 300EZ24		Red	NC
☐ 300EZ25		Silver	NC
☐ 300EZ19		White	NC

BASE

	OPTION		DESCRIPTION	RETAIL
Service Decals				$
☐	300D10		Army	NC
☐	300D11		Air Force	NC
☐	300D13		Marines	NC
☐	300D14		Navy	NC
☐	300D15		Purple Heart	NC
ISO Transit Tie Downs				$
☐	300TO50		ISO Transit Package	NC
LED Lights				$
☐	300LI08		LED Fender Lights	$646.00
☐	300LI09		LED Fender Lights *(Available with multi-power only)*	NC
☐	300LI02		Omit LED Lights	NC
Seat Heights				$
			Seat Height *Measured from front of seat pan to floor with 0° frame angle.*	
☐	300HT103		16.5"	NC
☐	300HT55		17"	NC
☐	300HT104		17.5"	NC
☐	300HT56		18"	NC
☐	300HT105		18.5"	NC
☐	300HT57		19" *(Only available with Tilt or Seat Elevator.)*	NC
Batteries and Chargers				$
Group 34 Battery Options				
☐	300AG27		GP34 Flat Top Pair w/Straps MK Gel	$955.00
☐	300AG29		GP34 Flat Top Pair w/Straps AGM	$863.00
Group 22 Battery				
☐	300AG31		GP22 Flat Top Pair w/Straps	$678.00
Omit Batteries				
☐	300AG24		Battery Harnesses Only *(Omit Batteries)*	NC
Chargers				
☐	300BJ30		10A Dual Mode (Compact)	$162.00
☐	300BJ19		8A Dual Mode (Compact)	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
Power Seating Options				$
Seat Elevator				
☐	300PI07		Seat Elevator, Powered, 300 mm / 12"	$3,477.00
☐	300PI08		C-Me® Seat Elevator, Powered, 300 mm / 12" *(C-Me speed up to 3 mph) (N/A with VR-2)*	$3,477.00
Seat Tilt				
☐	300PA48		Posterior 50° Tilt	$8,663.00
Seat Recline				
☐	300PA37		150° Powered Recline w/ 4" Shear Reduction	$8,001.00
Seat Tilt/Recline				
☐	300PA41		Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	$13,275.00
Power Recline Backrest Preset Home Angle				
☐	300QK1		90°	NC
☐	300QK2		95° (Default Position)	NC
Pre Set Seat Frame Angle *(No Power Seating Options, Increases Seat to Floor Height)*				$
☐	300FA10		0° *(Adds 0" to STFH)*	NC
☐	300FA3		3° *(Adds 1" to STFH)*	NC
Backrest Frame Type				$
☐	300BF16		Tool Angle Adjustable Back, 90°-120°	NC
☐	300BF14		User Adjustable Manual Semi-Recline Back, 90°-130°, Right *(Adjustment Lever Mounted on Opposite Side of Joystick)*	$633.00
☐	300BF15		User Adjustable Manual Semi-Recline Back, 90°-130°, Left *(Adjustment Lever Mounted on Opposite Side of Joystick)*	$633.00
Backrest Frame Option				$
☐	300BF13		Flip Down, Quick Release	$179.00
Back Cushion Type				$
☐	300BT134		SEDEO® Back Cushion *(Choose Seating in Section A)*	NC
☐	300BT138		JAY® Syncra™ *(Choose Seating in Section B)*	NC
☐	300BT135		JAY® J3® *(Choose Seating in Section C)*	NC
☐	300BT137		SEDEO® PRO™ Aftermarket Back Adapter *(Height Adjustable Round Tubes for Traditional Mounting)*	$633.00
☐	300BT133		OMIT *(For personalized Mounting. No Headrest Mount)*	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
A - SEDEO® Backrest Dimensions				$
Back Width				
☐	300BN28		400mm / 16"	NC
☐	300BN29		460mm / 18"	NC
☐	300BN30		500mm / 20"	NC
☐	300BN31		560mm / 22"	NC

			Back Shell Height	Tool Angle Adjust.	150° Recline	Installed Back Height	
☐	300B255		460mm / 18"	22"-23"	23"-25"		NC
☐	300B256		560mm / 22"	24"-27"	26"-29"		NC
☐	300B257		660mm / 26"	28"-30"	30"		NC

Installed Back Height - *(From Pan to Top of Back)*

☐	300JR32		22"	NC
☐	300JR33		23"	NC
☐	300JR34		24"	NC
☐	300JR35		25"	NC
☐	300JR36		26"	NC
☐	300JR37		27"	NC
☐	300JR38		28"	NC
☐	300JR39		29"	NC
☐	300JR40		30"	NC

B - JAY® Syncra™ Seating

B - JAY® Syncra™ Back P. 21 $

☐	300JP14		JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	$1,339.00
☐	300JP15		JAY® Syncra™ Back 22" Wide with 3" Contour Depth	$1,339.00

B - JAY® Syncra™ Back Width $

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	300JQ07		16" Wide *(Adjustment range 16 to 17" Wide)**	NC
☐	300JQ08		18" Wide *(Adjustment range 18 to 19" Wide)**	NC
☐	300JQ09		20" Wide *(Adjustment range 20 to 21" Wide)**	NC
☐	300JQ10		22" Wide *(Adjustment range 22 to 23" Wide)*	NC

**Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.*

B - JAY® Syncra™ Back Shell Options and Height Ranges $

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.

JAY® Syncra™ Shell Height

☐	300JR42		16" Tall Shell *(22" to 24" Back Height)*	NC
☐	300JR43		18" Tall Shell *(24" to 26" Back Height)*	NC
☐	300JR44		20" Tall Shell *(26" to 28" Back Height)*	NC
☐	300JR45		22" Tall Shell *(28" to 30" Back Height)*	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
			JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)	
☐	300JR50		22" *(16" Shell Height)*	NC
☐	300JR51		23" *(16" Shell Height)*	NC
☐	300JR52		24" *(16" Shell Height)*	NC
☐	300JR53		24" *(18" Shell Height)*	NC
☐	300JR54		25" *(18" Shell Height)*	NC
☐	300JR55		26" *(18" Shell Height)*	NC
☐	300JR56		26" *(20" Shell Height)*	NC
☐	300JR57		27" *(20" Shell Height)*	NC
☐	300JR58		28" *(20" Shell Height)*	NC
☐	300JR67		28" *(22" Shell Height)*	NC
☐	300JR68		29" *(22" Shell Height)*	NC
☐	300JR69		30" *(22" Shell Height)*	NC
	B - JAY® Syncra™ Back Upholstery			$
☐	300BU96		Syncra™ Premium Leatherette	NC
☐	300BU97		Syncra™ Stretch	NC
	B - JAY® Syncra™ Fixed Lateral Support Accessory			$
			JAY® Fixed Lateral Support Accessory Hardware *(JAY® Fixed Lateral Support is height adjustable)*	
☐	300JW33		Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	$322.00
☐	300JW34		Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	$322.00
	C - JAY® J3® Back			
	C - JAY® J3® Back - Contour			$
☐	300JP01		J3® PA Back for 16"-20" Wide	$912.00
☐	300JP02		J3® PA Back for 22" Wide	$1941.00
☐	300JP11		J3® PL Back for 16"-20" Wide	$986.00
☐	300JP13		J3® PD Back for 16"-20" Wide *(18" and 21" Shell Height Only)*	$1,068.00
	C - JAY® J3® Back Width			$
			Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.	
☐	300JQ02		16"	NC
☐	300JQ03		18"	NC
☐	300JQ04		20"	NC
☐	300JQ05		22"	NC

POWERED SEATING

	OPTION		DESCRIPTION				RETAIL
C - JAY® J3® Back Shell Options and Height Ranges							$
JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.							
JAY® J3® Shell Height			**Tool Angle Adjust.**	**150° Recline**	**Installed Back Height**		
☐	300JR41		15" tall shell *(Upper Thoracic Med) for 22" Wide*	22"-23"	22"-23"		NC
☐	300JR14		16.5" tall shell *(Upper Thoracic Short)*	21"-23"	22"-24"		NC
☐	300JR15		18" tall shell *(Upper Thoracic Med)*	22"-24"	23"-25"		NC
☐	300JR16		19.5" tall shell *(Upper Thoracic Tall)*	24"-26"	24"-26"		NC
☐	300JR18		18" tall shell *(Shoulder Height) for 22" wide*	24"	24"		NC
☐	300JR17		21" tall shell *(Shoulder Height)*	24"-26"	25"-28"		NC
☐	300JR27		22.5" tall shell *(Shoulder Height)*	26"-28"	27"-30"		NC
☐	300JR28		24" tall shell *(Shoulder Height)*	27"-29"	28"-30"		NC
JAY® J3® Installed Back Height *(From Pan to Top of Back)*							
☐	300JR29		21"				NC
☐	300JR19		22"				NC
☐	300JR20		23"				NC
☐	300JR21		24"				NC
☐	300JR22		25"				NC
☐	300JR23		26"				NC
☐	300JR24		27"				NC
☐	300JR25		28"				NC
☐	300JR30		29"				NC
☐	300JR31		30"				NC
C - JAY® Anterior Thoracic Support							$
Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)							
☐	300K10		Medium *(ATS8302)*				$322.00
☐	300K11		Large *(ATS8303)*				$322.00

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)				$
JAY® Lateral Thoracic Support Pads				
☐	300JW35		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	$72.00
☐	300JW38		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	$72.00
☐	300JW36		Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	$72.00
☐	300JW39		Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	$72.00
☐	300JW37		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	$72.00
☐	300JW40		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	$72.00
JAY® Lateral Thoracic Support Pad Upholstery				
☐	300NU01		Syncra™ Premium Leatherette	NC
☐	300NU02		Syncra™ Stretch	NC
JAY® Lateral Thoracic Swing-Away Mount				
☐	300JW41		Lateral Thoracic Swing Away Bracket - Right	$322.00
☐	300JW42		Lateral Thoracic Swing Away Bracket - Left	$322.00
Headrests *(Requires Pad, Upholstery, & Mount)*				$
Headrest Pad				
☐	300HE36		SEDEO® Pad Small with Cover (LxH) 240 x 150 mm / 10" x 6" *Can only be used with Multi Position Headrest Mount with quick-release and SEDEO® upholstery*	$227.00
☐	300HE38		Whitmyer® Plush Headrest Pad 203 mm / 8"	$377.00
☐	300HE39		Whitmyer® Plush Headrest Pad 254 mm / 10"	$377.00
☐	300HE44		Whitmyer® Plush Headrest Pad 356 mm / 14"	$377.00
☐	300HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium *(Adjustable 7" to 16" Wide x 4" Tall)*	$445.00
☐	300HE48		OMIT Headrest Pad	NC
Headrest Pad Upholstery				
☐	300NU10		SEDEO® Upholstery	NC
☐	300NU04		Syncra™ Premium Leatherette	NC
☐	300NU11		Whitmyer® Lycra / Stretch	NC
Headrest Mount				
☐	300MH32		SEDEO® Multi-Position Headrest Mount with Quick Release *(Only available with SEDEO® Pad)*	$557.00
☐	300MH33		Whitmyer® Linx Headrest Mount with Quick Release	$247.00
☐	300MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	$298.00
☐	300MH35		Whitmyer® AXYS Headrest Mount with Quick Release	$247.00
☐	300MH17		OMIT Headrest Mount	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
Push Handles				$
☐	300PH27		Height Adjustable Push Handles	$61.00
Seat Frame Rails				$
☐	300SN14		Seat Rail Short, for Effective Seat Depth 14" - 18"	$806.00
☐	300SN15		Seat Rail Long, for Effective Seat Depth 17" - 22"	$806.00
Seat Cushion Type				$
☐	300ST81		JAY® Syncra™ (See Section D)	NC
☐	300ST80		JAY® (See Section E)	NC
☐	300ST75		Omit	NC
Seat Size				$
Seat Pan Width				
☐	300R53		16"	NC
☐	300R54		17"	NC
☐	300R55		18"	NC
☐	300R56		19"	NC
☐	300R57		20"	NC
☐	300R58		21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐	300R59		22"	NC
Seat Pan Depth ***(Front of seat pan to front of back cushion – effective depth)***				
☐	300RF52		14"	NC
☐	300RF53		15"	NC
☐	300RF51		16"	NC
☐	300RF45		17"	NC
☐	300RF46		18"	NC
☐	300RF47		19"	NC
☐	300RF48		20"	NC
☐	300RF49		21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐	300RF50		22"	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
D - JAY® Syncra™ Cushion	P. 21 (JAY Your Way modifications available through Customer Service (800) 263-3390)			$
☐	300CU89		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(16" to 20" wide)*	$1,017.00
☐	300CU90		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(22" wide)*	$1,125.00
JAY® Syncra™ Cushion Cover Upholstery				
☐	300YR73		Syncra™ Premium Leatherette	NC
☐	300YR74		Syncra™ Stretch	NC
E - JAY® Seat Cushion				$
☐	300CU59		JAY® Ion® *(Sizes 16" to 21" Wide)*	$505.00
☐	300CU60		JAY® Ion® *(Size 22" Wide)*	$698.00
☐	300CU61		JAY® Union® *(Sizes 16" to 21" Wide)*	$764.00
☐	300CU62		JAY® Union® *(Size 22" Wide)*	$806.00
☐	300CU64		JAY® Fusion® *(Sizes 16" to 21" Wide)*	$1,045.00
☐	300CU65		JAY® Fusion® *(Size 22" Wide)*	$1,172.00

POSITIONING SUPPORTS & ACCESSORIES

	OPTION		DESCRIPTION	RETAIL
Knee and Thigh Swing-Away Support Pads				$
Lateral Knee Support Pads and Hardware *(must pick pad and mount)*				
☐	300XY05		Lateral Knee Support Pad, Right *(4" round)*	$247.00
☐	300XY06		Lateral Knee Support Pad, Left *(4" round)*	$247.00
☐	300XX08		Lateral Swing-Away Bracket Knee - Right	$325.00
☐	300XX09		Lateral Swing-Away Bracket Knee - Left	$325.00
Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*				
☐	300KF22		Lateral Thigh Support Pad Short - Right *(4.8" H x 4.8 " L / 120mm x 120mm)*	$146.00
☐	300KF23		Lateral Thigh Support Pad Short - Left *(4.8" H x 4.8 " L / 120mm x 120mm)*	$146.00
☐	300KF24		Lateral Thigh Support Pad Long - Right *(4.8" H x 7" L / 120mm x 120mm)*	$146.00
☐	300KF25		Lateral Thigh Support Pad,Long - Left *(4.8" H x 7" L / 120mm x 120mm)*	$146.00
☐	300KF26		Lateral Thigh Removable Mount - Right	$325.00
☐	300KF27		Lateral Thigh Removable Mount - Left	$325.00
JAY® Lateral Thigh Supports *(Must pick pad, upholstery, and mount)*				$
JAY® Lateral Thigh Support Pads				
☐	300KF16		Lateral Thigh Support Pad Standard, Right *(8.25" x 3.25")*	$152.00
☐	300KF17		Lateral Thigh Support Pad Standard, Left *(8.25" x 3.25")*	$152.00
☐	300KF18		Lateral Thigh Support Pad Large, Right *(12" x 3.25")*	$15200
☐	300KF19		Lateral Thigh Support Pad Large, Left *(12" x 3.25")*	$152.00
JAY® Lateral Thigh Support Pad Upholstery				
☐	300NU07		Syncra™ Premium Leatherette	NC
☐	300NU08		Syncra™ Stretch	NC
JAY® Lateral Thigh Removable Mount				
☐	300KF20		Lateral Thigh Removable Mount - Right	$325.00
☐	300KF21		Lateral Thigh Removable Mount - Left	$325.00

POSITIONING SUPPORTS & ACCESSORIES

	OPTION		DESCRIPTION	RETAIL
Armrest Structure				**$**
☐	300EG41		Reclining, Flip Back, Short *(Best used with shorter seat depths) (Height adjustable 9" to 15")*	NC
☐	300EG36		Reclining, Flip Back with 2-point post locking arm. *(Height adjustable 11" to 15" x 0.625")*	$351.00
☐	300EG51		Reclining, Flip Back, Dual Link Armrest *(No Charge for Multi-Power) (Height adjustable 8" to 15" including length and lateral adjustment)*	NC
☐	300EG52		Reclining, Flip Back, Dual Link Armrest *(Height adjustable 8" to 15" including length and lateral adjustment)*	NC
☐	300EG32		Single Post, Flip Back *(N/A with thigh support) (Height adjustable 10" to 15")*	NC
☐	300EG50		Dual Post HA Removable	NC
Arm Pads				**$**
SEDEO® Armrest Pad				
☐	300AB43		300 mm / 12", Integral Foam, Right	NC
☐	300AE43		300 mm / 12", Integral Foam, Left	NC
☐	300AB44		400 mm / 16", Integral Foam, Right	NC
☐	300AE44		400 mm / 16", Integral Foam, Left	NC
Comfort Arm Pad *(Foam)*				
☐	300AB46		2" Wide Right	$113.00
☐	300AE46		2" Wide Left	$113.00
Comfort Arm Pad *(Gel)*				
☐	300AB48		2" Wide Right	$213.00
☐	300AE48		2" Wide Left	$213.00
Dual Post HA Armrest Pads				
☐	300AB74		ARMPAD 10" CLASSIC DSK FOAM R	NC
☐	300AB75		ARMPAD 14" CLASSIC FLL FOAM R	NC
☐	300AE74		ARMPAD 10" CLASSIC DSK FOAM L	NC
☐	300AE75		ARMPAD 14" CLASSIC FLL FOAM L	NC
☐	300AB76		ARMPAD 10" WTRFLL DSK FOAM R	NC
☐	300AB77		ARMPAD14" WTRFLL FLL FOAM R	NC
☐	300AE76		ARMPAD 10" WTRFLL DSK FOAM L	NC
☐	300AE77		ARMPAD 14" WTRFLL FLL FOAM L	NC
Ottobock® Arm Troughs				
☐	300AE06		Ottobock Arm Trough, Med (12") Left	$139.00
☐	300AB06		Ottobock Arm Trough, Med (12") Right	$139.00
☐	300AE07		Ottobock Arm Trough, Large (14") Left	$139.00
☐	300AB07		Ottobock Arm Trough, Large (14") Right	$139.00
☐	300OL01		Ottobock Hand Pad, Med (7") Left	$70.00
☐	300OR01		Ottobock Hand Pad, Med (7") Right	$70.00
☐	300OL02		Ottobock Hand Pad, Large (9") Left	$70.00
☐	300OR02		Ottobock Hand Pad, Large (9") Right	$70.00
☐	300AE78		Ottobock One Piece Arm Trough (18"), Left	$198.00

☐	300AB78		Ottobock One Piece Arm Trough (18"), Right	$198.00

Elbow Support $

☐	3003A02		Elbow Support Pad, RH 3" x 10" *(Includes side support)*	$228.00
☐	3003A03		Elbow Support Pad, LH 3" x 10" *(Includes side support)*	$228.00
☐	3003A06		Elbow Support Pad, RH 3" x 3"	$167.00
☐	3003A07		Elbow Support Pad, LH 3" x 3"	$167.00
☐	3003A04		Swing-Away Elbow Support Hardware, RH	$279.00
☐	3003A05		Swing-Away Elbow Support Hardware, LH	$279.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL

Side Guards *(Pair)* $

☐	300SG44		Side Guards Integral Foam *(11"L x 5.5"H) (N/A with Reclining Armrests)*	$356.00

Fixed Centermount Legrest $

☐	300H183		SEDEO® Seating Fixed Centermount Legrest *(3" to 17.5" Leg Length)*	$351.00
☐	300H184		Fixed Centermount Legrest, Short *(Includes Calf Pads and Single Footplate) (11.5" to 14.5" Leg Length)*	NC
☐	300H185		Fixed Centermount Legrest, Long *(Includes Calf Pads and Single Footplate) (13.5" to 18.5" Leg Length)*	NC

Power Centermount Legrest $

☐	300H175		Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Mechanical Articulation) (5" to 17.5" Leg Length)*	$2,487.00
☐	300H109		Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Powered Articulation) (3" to 17.5" Leg Length)*	$3,435.00

Centermount Legrest Accessories $

Footplates

☐	300F127		Dual, Height, Width and Angle Adjustable, Standard *(14"W x 10"D)*	$284.00
☐	300F128		Dual, Height, Width and Angle Adjustable, Large *(16"W x 12"D)*	$284.00

Footplate Accessories

☐	300DI68		Extra Large Centermount Footplate *(Single Footplate Overlay 16"W x 12"D)*	$315.00
☐	300DI69		Footplate Connector Link *(Attaches footplates together)*	$96.00
☐	300DI57		Gel Footplate *(Pair, LH and RH)*	$401.00

Calf Pad Options

☐	300DI34		Calf Pads, Pair, Standard Integral Foam	NC
☐	300DI35		Calf Pads, Pair, Padded *(N/A w/ Std footplate)*	$356.00

SEDEO® PRO™ Swing-Away Legrest *(Pair) (Footplates, calf supports & options must be selected separately.)* $

☐	300H177		Legrests, Swing-Away, 70°, SEDEO® PRO™ *(14" to 21" Leg Length)*	NC
☐	300H178		Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (14.5" to 21.5" Leg Length)*	$380.00
☐	300H179		Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)*	$3,435.00
☐	300DI76		Omit Legrests	NC

Swing-Away Legrest Footplates

☐	300F369	Angle Adjustable Flip-Up Footplates	$284.00

Swing-Away Legrest Calf Supports

☐	300DI42	Calf Support, Depth Adjustable, Padded *(6.5" to 8" Wide, depending on chair width, x 9.5" Tall)*	$512.00

Swing-Away Legrest Options

☐	300DI44	Heel Loops	NC
☐	300KG4	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	$282.00

SEDEO® Lower Extremity Swing-Away Support

☐	300VD01	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	$512.00
☐	300VD02	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	$512.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL

Pelvic Positioning Belts

☐	300A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐	300A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
☐	300A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	$103.00
☐	300A845	Bodypoint – Evoflex Pelvic Stabilizer Belt	$216.00
☐	300A846	Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	$271.00

Accessories $

Ventilator Option

☐	300VT12	Vent Tray Fixed *(Pulmonetic LTV Models)*	$609.00

O2 Tank Option

☐	300A718	O2 Tank Holder Mounted to Backrest *(D and E Styles)*	$391.00

Wheelchair Tray Table

☐	3007A09	Tray, Swing-Away	$381.00

Backpack

☐	300A1	Black Utility	$112.00

Medial Thigh Support

☐	3008A80	Medial Thigh Support	$271.00
☐	3008A81	Pommel Swing-Away Hardware	$390.00

Misc. Accessories

☐	3008A59		Transfer Bar *(Pair)*	$493.00
☐	3008A52		Ram Mount Phone	$77.00
☐	3008A53		Ram Mount Tablet	$308.00
☐	3008A54		Ram Mount Cup Holder	$77.00
☐	3008A58		Rear View Mirror	$77.00
☐	300A611		Crutch Holder	$135.00
☐	3008A73		Bag For Armrest	$77.00
☐	3008A75		Reflectors *(Backrest shroud mounted)*	$77.00
☐	3009A10		Utility Hooks	$58.00

JOYSTICK / ELECTRONICS

VR-2 Joystick			$
☐	300J127	VR-2 Joystick	NC

R-net Joystick *(Remote)*			$
☐	300J121	Standard LED Joystick *(w/o Screen Display)*	NC
☐	300J116	Colour Joystick *(2.3" Colour LCD Screen Display, CJSM1)*	NC
☐	300J118	High Resolution Colour Joystick, (CJSM2+BT) *3.5" Colour LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control*	$1,162.00
☐	300J119	Omit Remote Control *(Required with arm mounted specialty controls)*	NC

Joystick Mounting			$
☐	300JM34	Fixed on Armrest	NC
☐	300JM35	Parallel Swing-Away	$447.00
☐	300JM36	Link-It Mid Line Mount on Armrest *(Swing-Away)*	$447.00
☐	300JM49	Attendant *(Available with LED or CJSM1 only. Mounted to Push Handle)*	NC

Joystick Mount (Location)			$
☐	300HM16	Right Hand	NC
☐	300HM17	Left Hand	NC

Display Language Selection R-net			$
☐	300L1	English	NC
☐	300L3	French	NC
☐	300L7	Spanish	NC

Joystick Handle (Optional)			$
☐	300JH13	Ball Handle	$110.00
☐	300JH15	T-Handle	$110.00
☐	300JH24	Bodypoint – Rubber Dome Joystick Handle	$216.00
☐	300JH25	Bodypoint – U-Shaped Joystick Handle *(4 in)*	$216.00
☐	300JH26	Bodypoint – Grooved Mushroom Joystick Handle	$216.00

Electronic Programmers			$
☐	300PR13	R-net PC Wireless Programmer	$660.00

R-net Power Modules			$
☐	300ER12	Power Module, VR-2 90A, Non-Expandable	NC
☐	300ER22	Power Module, R-net, 90A, Non-Expandable *(Only available with single power)*	NC
☐	300ER13	Power Module, R-net, 90A, Expandable *(Required for Specialty Controls or 2+ Actuators)*	$102.00
☐	300EX11	Expandable 90A Controller - R-net	NC
☐	300EX12	Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

	OPTION	DESCRIPTION	RETAIL

Seating Control			$
Powered Seating Control			
☐	300EO40	Omit Thru-Drive *(Requires Switch Control)*	NC
☐	300EO37	Thru-drive Control for 1 Actuator	$1,443.00
☐	300EO38	Thru-drive Control for 2+ Actuators	$2,903.00
Switch Control *(Available with Omit Thru-Drive only) (1 Actuator for Tilt or Elevate only)*			
☐	300EO33	Dual Toggle for 1 Actuator	NC
☐	300EO34	2.5" Single Button Switch for 1 Actuator	NC
☐	300EO41	Ctrl+5 Switch Box for 2 Actuators or more	NC
Additional Power Seating Control *(Available with Thru-Drive only)*			
☐	300QG11	Ctrl+5 Switch Box with Voice Activation	$1,078.00
☐	300QG12	Ctrl+5 Switch Box with Voice Activation *(Multi-power only)*	NC
Dual Toggle Mounting Location			
☐	300TN1	Toggle Mount, Right	NC
☐	300TN2	Toggle Mount, Left	NC
Ctrl+5 Mounting Location *(Will swing away if mounted on the same side as the joystick)*			
☐	300EY19	Right Armrest	NC
☐	300EY20	Left Armrest	NC
☐	300EY21	Back Cane, Right	NC
☐	300EY24	Back Cane, Left	NC
Ctrl+5 Switch Type			
☐	300EY22	Push Buttons	NC
☐	300EY23	Toggle Switches	$356.00

Attendant Controls			$

☐ 3000A05		Dual Control Module for Attendant	$1,225.00
☐ 3000A06		R-net Colour LCD (CJSM1), Control Attendant Only	$1,225.00

Attendant Mounting Location			$
☐ 300DN5		Right Back Cane	NC
☐ 300DN6		Left Back Cane	NC

JOYSTICK / ELECTRONICS

OPTION		DESCRIPTION	RETAIL

R-net Programming $

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

☐ 300PM12		Include Assignable Buttons	NC

Power Centermount Footplate Extension Option

☐ 300IE03		Include 7.5" Powered Footplate Extension for Transfer Assist *(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation*	$592.00

R-net Omni2 Options *For Expandable Controls, choose all that apply.* $

☐ 300ON14		R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Choose Mount Below	$2,384.00
☐ 300ON15		R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.* $

☐ 300DS14		Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Right	$698.00
☐ 300DS15		Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Left	$698.00

SPECIALTY INPUT DEVICES

	OPTION		DESCRIPTION	RETAIL
Hand Control Packages	*Requires Standard Joystick Omit / Only one hand control device can be ordered.*			$
Proportional Input Devices				
☐	300EU14		Switch-It MicroGuide	$3,844.00
☐	300EU17		Switch-It VersaGuide	$2,557.00
☐	300EU18		Switch-It VersaGuide EZ	$3,023.00
Mounting Position *(Required)*				
☐	300EM04		Link-It to the Backrest *(Not Available w/ Omit Back Cushion Type)*	$1,338.00
☐	300EM05		Link-It to the Fixed Joystick Mount	$416.00
☐	300EM06		Link-It to the Parallel Swing-Away Joystick Mount	$416.00
☐	300EM07		Link-It to the Midline Swing-Away	$1,480.00
☐	300EM08		Omit Mount	NC
Mounting Platforms *(Optional)*	*Not available w/ Link-It to Backrest*			
☐	300MG16		SunDisk Hand Tray *(Only available w/ MicroSeries Joysticks)*	$416.00
☐	300MG15		Bullet Hand Tray	$416.00
Mounting Location *(Required)*				
☐	300IM09		Right	NC
☐	300IM10		Left	NC
Switch Options				
☐	–		See Switch Selections *(Maximum 2)*	–
Head Control Packages	*(Select headrest mount on next page)*			$
Only one head control device can be ordered.				
Proportional Input Devices				
☐	300EU25		Switch-It DualPro® Proportional Head Array - Facial Pads	$7,057.00
☐	300EU26		Switch-It DualPro® Proportional Head Array - Spot Pads	$7,057.00
☐	300EU51		Switch-It Vigo Head Drive	$6,293.00
Switched Input Devices				
☐	300EU28		Switch-It 3-Switch Head Array - Short Curved Wings	$6,154.00
☐	300EU30		Switch-It 3-Switch Head Array - Long Curved Wings	$6,154.00
Switch Options				
☐	–		See Switch Selections *(Maximum 2)*	–

SPECIALTY INPUT DEVICES

	OPTION		DESCRIPTION	RETAIL
Sip 'N Puff				$
☐	300ID73		Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware *(Requires headrest mount)*	$7,037.00
☐	300ID74		Switch-It OptiStop with Sip 'N Puff	$2,741.00
Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)*				$
☐	300IJ13		Whitmyer® Linx Headrest Mount with Quick Release	$247.00
☐	300IJ14		Whitmyer® Linx-2 Headrest Mount with Quick Release	$298.00
☐	300IJ15		Whitmyer® AXYS Headrest Mount with Quick Release	$247.00
Chin Control Packages				$
Proportional Input Devices				
☐	300EU34		Switch-It MicroGuide *(Switches mounted to the headrest only.)*	$4,529.00
☐	300EU36		Switch-It VersaGuide	$2,557.00
☐	300EU37		Switch-It VersaGuide EZ	$3,023.00
Mounting Position				
☐	300EM09		Link-It To Chin Harness, Small (150mm / 6")	$1,760.00
☐	300EM10		Link-It To Chin Harness, Large (220mm / 8")	$1,760.00
☐	300EM11		Link-It To Headrest, Swing-Away	$1,760.00
☐	300EM16		Omit Mount	NC
Input Device Mounting Location				
☐	300IM13		Right	NC
☐	300IM14		Left	NC
Switch Options				
☐	–		See Switch Selections	–
Safety Devices				$
☐	300SA06		Switch-It Remote Stop Switch - XLR Connector	$1,317.00
Switched Driving Packages				$
☐	300SW44		Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. *(Pro Spot proportional switches are required if proportional driving is desired.)* *(For non-proportional driving, single switches can be used.)*	$3,792.00

SPECIALTY INPUT DEVICES

OPTION		DESCRIPTION	RETAIL

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.
Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls
Switch ordering instructions:

1. Select switch type and quantity in Section "A" using colums below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections — **Switch / Quantity**

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	$
2.5" Single Switch - Black	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Green	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Red	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$230.00
Egg Switch - Black	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Blue	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Green	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Red	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Yellow	☐	☐	☐	☐	☐	☐	$131.00
MicroLite *(Includes colour labels)*	☐	☐	☐	☐	☐	☐	$195.00
Mini Button 1 - Blue *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
Mini Button 1 - Green *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
Mini Button 1 - Red *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
1" Single Switch - Black	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Green	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Orange	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Red	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$326.00
Electronic Switches							
Flex Switch (Headrest mount only)	☐	☐	☐	☐	☐	☐	$254.000
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Green	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Blue	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Yellow	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Orange	☐	☐	☐	☐	☐	☐	$562.00
Section B – Switch Mounting Position							
Link-It to the Backrest *(Up to 1 per side)*	☐	☐	☐	☐	☐	☐	$770.00
Link-It to the Headrest - Gooseneck & Rigid *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$770.00
Link-It to the Fixed Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$587.00
Link-It to the Parallel Swing-Away Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$587.00
Link-It to the Seat Rail *(Up to 2 per side)* *(Not avail w/knee & thigh supports)*	☐	☐	☐	☐	☐	☐	$467.00
Omit Mount *(includes Velcro)*	☐	☐	☐	☐	☐	☐	NC
Section C – Switch Mounting Location							
Right side of chair	☐	☐	☐	☐	☐	☐	NC
Left side of chair	☐	☐	☐	☐	☐	☐	NC

Joystick Integrated Switches *(Pair)*

☐	300J003		Proximity Switches on Flex Mount, Pair *(Used with VersaGuide/VersaGuide EZ)*	$1,178.00

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

Comments for special build

built4me

Q500 M with SEDEO® PRO™	$
Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Signature

Send E-Mail

SUNRISE MEDICAL

www.SunriseMedical.ca
800.263.3390
E-mail: cscanada@sunmed.com

Specifications are subject to change without notice.

© 2025 Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, Ontario L4K 4P8
10.2026 MK-129646-CAN Rev. 24