SUNRISE MEDICAL

QUICKIE®
Live without limits.



ORDER FORM


Q700M with SEDEO® PRO™, PRO™ ADVANCED and JAY® Syncra™ Seating

EIPW25

Customer		Ship To			
Account #		Name			
Date		Address			
PO#		City			
Buyer		Province		Postal Code	
Mark For		Phone			

Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

Chair Specifications

Typical Weight w/out Batteries	330 lbs.
Typical Weight with Batteries	410 lbs.
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction / 12" seat elevator
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 24.5", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Colour	9 shroud colours and 7 colour accents
Maximum Ground Clearance	3"

BASE

OPTION		DESCRIPTION	RETAIL
Standard Duty Mid Wheel Drive Base *End Item Number EIPW25 (PRO™ Seat)*			$
☐ 250PC197		Quickie Q700 M NPC No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,993.00
☐ 250PC198		Quickie Q700 M SPC Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,993.00
☐ 250PC203		Quickie Q700 M MPC Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,99300
Standard Duty Mid Wheel Drive Base *End Item Number EIPW25 (PRO™ Seat)*			$
☐ 250PC201		Quickie Q700 M NPD No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	$10,781.00
☐ 250PC202		Quickie Q700 M SPD Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	$10,781.00
☐ 250PC205		Quickie Q700 M MPD Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	$10,781.00

BASE

OPTION		DESCRIPTION	RETAIL
User Weight			$
☐ 250WG23		Up to 300 lbs.	NC
Motor Package			$
☐ 250PF99		4-pole 10 kph / 6 mph	NC
☐ 250PF100		4-pole 12.5 kph / 8 mph	NC
☐ 250PF101		4-pole 12.5 kph / 8 mph	$1,268.00
Tracking			$
☐ 250MI01		SureTrac® Module & Motor Encoders *(Improves drive tracking)*	$1,696.00
Drive Wheels			$
☐ 250DW58		Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			$
☐ 250CW55		Caster Wheels 6", Solid, Black	NC
Shroud Colour			$
☐ 250C119		Matte Black	NC
☐ 250C165		Candy Apple Red	NC
☐ 250C166		Candy Blue	NC
☐ 250C167		Candy Purple	NC
☐ 250C176		Green Apple	NC
☐ 250C178		Pearl Pink	NC
☐ 250C182		Silver	NC
☐ 250C125		Pearl White	NC
Colour Accents			$

Note: Colour Accents and Shroud Colours may not match.

☐ 250EZ20		Green	NC
☐ 250EZ22		Black	NC
☐ 250EZ21		Blue	NC
☐ 250EZ23		Orange	NC
☐ 250EZ24		Red	NC
☐ 250EZ25		Silver	NC
☐ 250EZ19		White	NC

BASE

	OPTION		DESCRIPTION	RETAIL
Service Decals				$
☐	250D10		Army	NC
☐	250D11		Air Force	NC
☐	250D13		Marines	NC
☐	250D14		Navy	NC
☐	250D15		Purple Heart	NC
ISO Transit Tie Downs				$
☐	250TO50		ISO Transit Package	NC
LED Lights				$
☐	250LI08		LED Fender Lights	$646.00
☐	250LI09		LED Fender Lights *(Available with Multi-Power Only)*	NC
☐	250LI02		Omit LED Lights	NC

Seat Heights

*Note: Available seat-to-floor heights are dependent on 2 options: **1**- Battery Size and **2**- Powered Seating.*
Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

					PRO ADVANCED Recline		
		Seat Height	Power Seating	No Power Seating	w/o Elevate	w/ Elevate	
☐	250HT103		16.5"	Group 34	Group 34	N/A	N/A	NC
☐	250HT55		17"	Group 34	GP24 or GP34	GP34	GP34	NC
☐	250HT104		17.5"	Group 34	GP24 or GP34	GP34	GP34	NC
☐	250HT56		18"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP34	NC
☐	250HT105		18.5"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC
☐	250HT57		19"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC

Batteries and Chargers

				$
Group 24 Battery Options *(Recommended for max daily range)*				
☐	250AG25		GP24 Flat Top Pair w/Strap MK Gel	$955.00
☐	250AG26		GP24 Flat Top Pair w/Strap AGM	$915.00
Group 34 Battery Options *(Required for lower seat-to-floor heights)*				
☐	250AG27		GP34 Flat Top Pair w/Straps MK Gel	$955.00
☐	250AG29		GP34 Flat Top Pair w/Straps AGM	$863.00
Omit Batteries				
☐	250AG24		Battery Harnesses Only *(Omit Batteries)*	NC
Chargers				
☐	250BJ28		10A Dual Mode (Compact)	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
Power Seating Options				$
Seat Elevator				
☐	250PI08		C-Me® Seat Elevator, Powered, 300 mm / 12" *(C-Me speed up to 3 mph)*	$3,477.00
☐	250PI07		Seat Elevator, Powered, 12" *(Required for Anterior Tilt)*	$3,477.00
Seat Tilt				
☐	250PA35		Posterior 50° Tilt	$8,663.00
Seat Recline				
☐	250PA37		150° Powered Recline w/ 4" Shear Reduction	$8,001.00
Seat Tilt/Recline				
☐	250PA41		Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	$13,275.00
☐	250PA62		Posterior 50° Tilt and PRO ADVANCED 165° Recline Powered, w/ 7" Shear Reduction	$13,785.00
Seat Anterior Tilt				
☐	250TG02		Anterior Tilt 0° to 30°	$3,329.00
Power Recline Backrest Preset Home Angle				
☐	250QK1		90°	NC
☐	250QK2		95° (Default Position)	NC
Pre Set Seat Frame Angle *(No Power Seating Options, Increases Seat to Floor Height)*				$
☐	250FA10		0° *(Adds 0" to STFH)*	NC
☐	250FA3		3° *(Adds 1" to STFH)*	NC
Backrest Frame Type *(N/A with Powered Recline)*				$
☐	250BF10		Tool Angle Adjustable Back, 90°-120°	NC
☐	250BF11		User Adjustable Manual Semi-Recline Back, 90°-130°, Right *(Adjustment lever mounted on opposite side of joystick)*	$824.00
☐	250BF12		User Adjustable Manual Semi-Recline Back, 90°-130°, Left *(Adjustment lever mounted on opposite side of joystick)*	$824.00
Backrest Frame Option *(N/A with PRO ADVANCED)*				$
☐	250BF13		Flip Down, Quick Release	$179.00
Back Cushion Type				$
☐	250BT134		SEDEO® Back Cushion *(Choose Seating in Section A) (N/A with PRO ADVANCED)*	NC
☐	250BT138		JAY® Syncra™ *(Choose Seating in Section B)*	NC
☐	250BT135		JAY® J3® *(Choose Seating in Section C)*	NC
☐	250BT137		SEDEO® Aftermarket Back Adapter *(Height adjustable round tubes for traditional mounting)*	$633.00
☐	250BT133		OMIT *(For personalized mounting. No headrest mount)*	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
A - SEDEO® Backrest Dimensions *(N/A with PRO ADVANCED)*				$
Back Width				
☐	250BN28		400mm / 16"	NC
☐	250BN29		460mm / 18"	NC
☐	250BN30		500mm / 20"	NC
☐	250BN31		560mm / 22"	NC

Back Shell Height

				Tool Angle Adjust.	150° Recline	Installed Back Height	
☐	250B255		460mm / 18"	22"-23"	23"-25"		NC
☐	250B256		560mm / 22"	24"-27"	26"-29"		NC
☐	250B257		660mm / 26"	28"-30"	30"		NC

Installed Back Height - *(From Pan to Top of Back)*

☐	250JR32		22"	NC
☐	250JR33		23"	NC
☐	250JR34		24"	NC
☐	250JR35		25"	NC
☐	250JR36		26"	NC
☐	250JR37		27"	NC
☐	250JR38		28"	NC
☐	250JR39		29"	NC
☐	250JR40		30"	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
B - JAY® Syncra™ Seating				
B - JAY® Syncra™ Back P. 23				$
☐	250JP14		JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	$1,339.00
☐	250JP15		JAY® Syncra™ Back 22" Wide with 3" Contour Depth	$1,339.00
B - JAY® Syncra™ Back Width				$
Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.				
☐	250JQ07		16" Wide *(Adjustment range 16 to 17" Wide)**	NC
☐	250JQ08		18" Wide *(Adjustment range 18 to 19" Wide)**	NC
☐	250JQ09		20" Wide *(Adjustment range 20 to 21" Wide)**	NC
☐	250JQ10		22" Wide *(Adjustment range 22 to 23" Wide)*	NC
**Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.*				
B - JAY® Syncra™ Back Shell Options and Height Ranges				$
JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.				
JAY® Syncra™ Shell Height				
☐	250JR42		16" Tall Shell *(22" to 24" Back Height)*	NC
☐	250JR43		18" Tall Shell *(24" to 26" Back Height)*	NC
☐	250JR44		20" Tall Shell *(26" to 28" Back Height)*	NC
☐	250JR45		22" Tall Shell *(28" to 30" Back Height)*	NC
JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)				
☐	250JR50		22" *(16" Shell Height)*	NC
☐	250JR51		23" *(16" Shell Height)*	NC
☐	250JR52		24" *(16" Shell Height)*	NC
☐	250JR53		24" *(18" Shell Height)*	NC
☐	250JR54		25" *(18" Shell Height)*	NC
☐	250JR55		26" *(18" Shell Height)*	NC
☐	250JR56		26" *(20" Shell Height)*	NC
☐	250JR57		27" *(20" Shell Height)*	NC
☐	250JR58		28" *(20" Shell Height)*	NC
☐	250JR67		28" *(22" Shell Height)*	NC
☐	250JR68		29" *(22" Shell Height)*	NC
☐	250JR69		30" *(22" Shell Height)*	NC
B - JAY® Syncra™ Back Upholstery				$
☐	250BU96		Syncra™ Premium Leatherette	NC
☐	250BU97		Syncra™ Stretch	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
C - JAY® J3® Back				
C - JAY® J3® Back - Contour				$
☐	250JP01		J3® PA Back for 16"-20" Wide	$912.00
☐	250JP02		J3® PA Back for 22" Wide	$1,941.00
☐	250JP11		J3® PL Back for 16"-20" Wide	$986.00
☐	250JP13		J3® PD Back for 16"-20" Wide *(18" and 21" shell height only)*	$1,068.00
C - JAY® J3® Back Width				$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	250JQ02		16"	NC
☐	250JQ03		18"	NC
☐	250JQ04		20"	NC
☐	250JQ05		22"	NC

C - JAY® J3® Back Shell Options and Height Ranges $

JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown, based on the back frame type and the amount of usable hardware adjustability.

JAY® J3® Shell Height

				Tool Angle Adjust.	150° Recline	PRO ADVANCED	Installed Back Height	
☐	250JR41		15" tall shell *(Upper Thoracic Med) for 22" Wide*	22"-23"	22"-23"	21"-24"		NC
☐	250JR14		16.5" tall shell *(Upper Thoracic Short)*	21"-23"	22"-24"	21"-24"		NC
☐	250JR15		18" tall shell *(Upper Thoracic Med)*	22"-24"	23"-25"	22"-26"		NC
☐	250JR16		19.5" tall shell *(Upper Thoracic Tall)*	24"-26"	24"-26"	23"-27"		NC
☐	250JR18		18" tall shell *(Shoulder Height) for 22" wide*	24"	24"	21"-25"		NC
☐	250JR17		21" tall shell *(Shoulder Height)*	24"-26"	25"-28"	25"-28"		NC
☐	250JR27		22.5" tall shell *(Shoulder Height)*	26"-28"	27"-30"	26"-29"		NC
☐	250JR28		24" tall shell *(Shoulder Height)*	27"-29"	28"-30"	28"-30"		NC

JAY® J3® Installed Back Height *(From Pan to Top of Back)*

☐	250JR29		21"	NC
☐	250JR19		22"	NC
☐	250JR20		23"	NC
☐	250JR21		24"	NC
☐	250JR22		25"	NC
☐	250JR23		26"	NC
☐	250JR24		27"	NC
☐	250JR25		28"	NC
☐	250JR30		29"	NC
☐	250JR31		30"	NC

POWERED SEATING

OPTION		DESCRIPTION	RETAIL
C - JAY® Anterior Thoracic Support			$
Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)			
☐ 250K10		Medium *(ATS8302)*	$322.00
☐ 250K11		Large *(ATS8303)*	$322.00
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)			$
JAY® Lateral Thoracic Support Pads			
☐ 250JW35		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	$75.00
☐ 250JW38		Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	$75.00
☐ 250JW36		Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	$75.00
☐ 250JW39		Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	$75.00
☐ 250JW37		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	$75.00
☐ 250JW40		Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	$75.00
JAY® Lateral Thoracic Support Pad Upholstery			
☐ 250NU01		Syncra™ Premium Leatherette	NC
☐ 250NU02		Syncra™ Stretch	NC
JAY® Lateral Thoracic Swing-Away Mount			
☐ 250JW41		Lateral Thoracic Swing Away Bracket - Right	$334.00
☐ 250JW42		Lateral Thoracic Swing Away Bracket - Left	$334.00
JAY® Syncra™ Fixed Lateral Support Accessory			$
JAY® Fixed Lateral Support Accessory Hardware *(JAY® Fixed Lateral Support is height adjustable)*			
☐ 250JW33		Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	$147.00
☐ 250JW34		Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	$147.00
Headrests (Requires Pad, Upholstery, & Mount)			$
Headrest Pad			
☐ 250HE36		SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" *Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery*	$227.00
☐ 250HE38		Whitmyer® Plush Headrest Pad 203 mm / 8"	$377.00
☐ 250HE39		Whitmyer® Plush Headrest Pad 254 mm / 10"	$377.00
☐ 250HE44		Whitmyer® Plush Headrest Pad 356 mm / 14"	$377.00
☐ 250HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium *(Adjustable from 7 to 16" Wide x 4" Tall)*	$438.00
☐ 250HE48		OMIT Headrest Pad	NC
Headrest Pad Upholstery			
☐ 250NU10		SEDEO® Upholstery	NC
☐ 250NU04		Syncra™ Premium Leatherette	NC
☐ 250NU05		Syncra™ Stretch	NC
☐ 250NU11		Whitmyer® Lycra	NC

POWERED SEATING

OPTION		DESCRIPTION	RETAIL
Headrest Mount			
☐ 250MH32		SEDEO® Multi-Position Headrest Mount with Quick Release *(Only available with SEDEO® Pad)*	$557.00
☐ 250MH33		Whitmyer® Linx Headrest Mount with Quick Release	$243.00
☐ 250MH34		Whitmyer® Linx-2 Headrest Mount with Quick Release	$298.00
☐ 250MH35		Whitmyer® AXYS Headrest Mount with Quick Release	$247.00
☐ 250MH17		OMIT Headrest Mount	NC
Push Handles			**$**
☐ 250PH27		Height Adjustable Push Handles	$61.00
Seat Frame Rails			**$**
☐ 250SN14		Seat Rail Short *(For effective seat depth 14" to 18")*	$806.00
☐ 250SN15		Seat Rail Long *(For effective seat depth 17" to 22")*	$806.00
Seat Cushion Type			**$**
☐ 250ST81		JAY® Syncra™ Cushion (See Section D)	NC
☐ 250ST80		JAY® (See Section E)	NC
☐ 250ST75		Omit	NC
C - Seat Size			**$**
Seat Pan Width			
☐ 250R53		16"	NC
☐ 250R54		17"	NC
☐ 250R55		18"	NC
☐ 250R56		19"	NC
☐ 250R57		20"	NC
☐ 250R58		21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐ 250R59		22"	NC
Seat Pan Depth *(Front of seat pan to front of back cushion – effective depth)*			
☐ 250RF52		14" *(N/A with PRO ADVANCED)*	NC
☐ 250RF53		15" *(N/A with PRO ADVANCED)*	NC
☐ 250RF51		16" *(N/A with PRO ADVANCED)*	NC
☐ 250RF45		17"	NC
☐ 250RF46		18"	NC
☐ 250RF47		19"	NC
☐ 250RF48		20"	NC
☐ 250RF49		21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐ 250RF50		22"	NC

POWERED SEATING

	OPTION		DESCRIPTION	RETAIL
D - JAY® Syncra™ Cushion	P. 23 (JAY Your Way modifications available through Customer Service (800-263-3390)			$
☐	250CU89		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(16" to 20" wide)*	$1,017.00
☐	250CU90		JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(22" wide)*	$1,125.00
JAY® Syncra™ Cushion Cover Upholstery				
☐	250YR73		Syncra™ Premium Leatherette	NC
☐	250YR74		Syncra™ Stretch	NC

E - JAY® Seat Cushion				$
☐	250CU59		JAY® Ion® *(Sizes 16" - 21" Wide)*	$657.00
☐	250CU60		JAY® Ion® *(Size 22" Wide)*	$698.00
☐	250CU61		JAY® Union® *(Sizes 16" - 21" Wide)*	$764.00
☐	250CU62		JAY® Union® *(Size 22" Wide)*	$806.00
☐	250CU64		JAY® Fusion® *(Sizes 16" - 21" Wide)*	$1,045.00
☐	250CU65		JAY® Fusion® *(Size 22" Wide)*	$1,188.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Knee and Thigh Swing-Away Support Pads			$
Lateral Knee Support Pads and Hardware *(must pick pad and mount)*			
☐ 250XY05		Lateral Knee Support Pad, Right *(4" round)*	$247.00
☐ 250XY06		Lateral Knee Support Pad, Left *(4" round)*	$247.00
☐ 250XX08		Lateral Swing-Away Bracket Knee - Right	$325.00
☐ 250XX09		Lateral Swing-Away Bracket Knee - Left	$325.00
Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*			
☐ 250KF22		Lateral Thigh Support Pad Short - Right *(5"x5")*	$146.00
☐ 250KF23		Lateral Thigh Support Pad Short - Left *(5"x5")*	$146.00
☐ 250KF24		Lateral Thigh Support Pad Long - Right *(5"x7")*	$146.00
☐ 250KF25		Lateral Thigh Support Pad Long - Left *(5"x7")*	$146.00
☐ 250KF26		Lateral Thigh Removable Mount - Right	$325.00
☐ 250KF27		Lateral Thigh Removable Mount - Left	$325.00
JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount)			$
JAY® Lateral Thigh Support Pads			
☐ 250KF16		Lateral Thigh Support Pad Standard, Right *(8.25" x 3.25")*	$146.00
☐ 250KF17		Lateral Thigh Support Pad Standard, Left *(8.25" x 3.25")*	$146.00
☐ 250KF18		Lateral Thigh Support Pad Large, Right *(12" x 3.25")*	$146.00
☐ 250KF19		Lateral Thigh Support Pad Large, Left *(12" x 3.25")*	$146.00
JAY® Lateral Thigh Support Pad Upholstery			
☐ 250NU07		Syncra™ Premium Leatherette	NC
☐ 250NU08		Syncra™ Stretch	NC
JAY® Lateral Thigh Removable Mount			
☐ 250KF20		Lateral Thigh Removable Mount - Right	$325.00
☐ 250KF21		Lateral Thigh Removable Mount - Left	$325.00
Armrest Structure			$
☐ 250EG41		Reclining, Flip Back, Short *(Best used with shorter seat depths) (Height adjustable 9" to 15")*	NC
☐ 250EG36		Reclining, Flip Back with 2-point post locking arm *(Height adjustable 11" to 15" x .625")*	$351.00
☐ 250EG51		Reclining, Flip Back, Dual Link Armrest *(No Charge for Multi-Power) (Required for PRO ADVANCED) (Height adjustable 8" to 15" including length and lateral adjustment)*	NC
☐ 250EG52		Reclining, Flip Back, Dual Link Armrest *(Height adjustable 8" to 15" including length and lateral adjustment)*	NC
☐ 250EG50		Dual Post HA Removable	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL
Arm Pads			**$**
SEDEO® Armrest Pad			
☐ 250AB43		300 mm / 12", Integral Foam, Right	NC
☐ 250AE43		300 mm / 12", Integral Foam, Left	NC
☐ 250AB44		400 mm / 16", Integral Foam, Right	NC
☐ 250AE44		400 mm / 16", Integral Foam, Left	NC
Comfort Arm Pad *(Foam)*			
☐ 250AB46		2" Wide Right	$113.00
☐ 250AE46		2" Wide Left	$113.00
Comfort Arm Pad *(Gel)*			
☐ 250AB48		2" Wide Right	$213.00
☐ 250AE48		2" Wide Left	$213.00
Dual Post HA Armrest Pads			
☐ 250AB74		ARMPAD 10" CLASSIC DSK FOAM R	NC
☐ 250AB75		ARMPAD 14" CLASSIC FLL FOAM R	NC
☐ 250AE74		ARMPAD 10" CLASSIC DSK FOAM L	NC
☐ 250AE75		ARMPAD 14" CLASSIC FLL FOAM L	NC
☐ 250AB76		ARMPAD 10" WTRFLL DSK FOAM R	NC
☐ 250AB77		ARMPAD14" WTRFLL FLL FOAM R	NC
☐ 250AE76		ARMPAD 10" WTRFLL DSK FOAM L	NC
☐ 250AE77		ARMPAD 14" WTRFLL FLL FOAM L	NC
Ottobock® Arm Troughs			
☐ 250AE06		Ottobock® Arm Trough, Med (12"), Left	$139.00
☐ 250AB06		Ottobock® Arm Trough, Med (12"), Right	$139.00
☐ 250AE07		Ottobock® Arm Trough, Large (14") Left,	$139.00
☐ 250AB07		Ottobock® Arm Trough, Large (14"), Right	$139.00
☐ 250OL01		Ottobock® Hand Pad, Med (7"), Left	$70.00
☐ 250OR01		Ottobock® Hand Pad, Med (7"), Right	$70.00
☐ 250OL02		Ottobock® Hand Pad, Large (9"), Left	$70.00
☐ 250OR02		Ottobock® Hand Pad, Large (9"), Right	$70.00
☐ 250AE78		Ottobock® One Piece Arm Trough, (18"), Left	$198.00
☐ 250AB78		Ottobock® One Piece Arm Trough, (18"), Right	$198.00

Elbow Support			$
☐	2503A02	Elbow Support Pad, RH 3" x 10" *(Includes side support)*	$228.00
☐	2503A03	Elbow Support Pad, LH 3" x 10" *(Includes side support)*	$228.00
☐	2503A06	Elbow Support Pad, RH 3" x 3"	$167.00
☐	2503A07	Elbow Support Pad, LH 3" x 3"	$167.00
☐	2503A04	Swing-Away Elbow Support Hardware, RH	$279.00
☐	2503A05	Swing-Away Elbow Support Hardware, LH	$279.00

Side Guards *(Pair)*			$
☐	250SG44	Side Guards Integral Foam *(11"L x 5.5"H) (N/A with Reclining Armrests)*	$356.00

Fixed Centermount Legrest			$
☐	250H183	SEDEO® Fixed Centermount Legrest *(3" to 17.5" Leg Length)*	$351.00
☐	250H184	Fixed Centermount Legrest, Short *(Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length)*	NC
☐	250H185	Fixed Centermount Legrest, Long *(Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length)*	NC

POSITIONING SUPPORTS & ACCESSORIES

OPTION		DESCRIPTION	RETAIL

Power Centermount Legrest			$
☐	250H175	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Mechanical Articulation) (5" to 17.5" Leg Length)*	$2,487.00
☐	250H109	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Powered Articulation) (3" to 17.5" Leg Length) (Required for PRO ADVANCED)*	$3,435.00

Power Centermount Legrest Accessories			$
Footplates			
☐	250F127	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	$284.00
☐	250F128	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	$284.00
Footplate Accessories *(Available with Power Centermounts)*			
☐	250DI68	Extra Large Centermount Footplate *(Single Footplate Overlay 16"W x 12"D)*	$315.00
☐	250DI69	Footplate Connector Link *(Attaches footplates together)*	$96.00
☐	250DI57	Gel Footpads *(Pair, LH and RH)*	$401.00
Calf Pad Options			
☐	250DI34	Calf Pads, Pair, Standard	NC
☐	250DI35	Calf Pads, Pair, Padded	$356.00
Knee Block *(Required for Anterior Tilt)*			
☐	250DI37	Knee Blocks: Dual, Height, Width and Angle Adjustable	$704.00
☐	250DI55	Knee Block Quick Release Removable Hardware	$179.00

SEDEO® PRO™ Swing-Away Legrest *(Pair)*			$
☐	250H177	Legrests, Swing-Away, 70°, SEDEO® PRO™ *(14" to 21" Leg Length)*	NC
☐	250H178	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)*	$605.00
☐	250H179	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)*	$3,435.00
☐	250DI76	Omit Legrests	NC

Swing-Away Legrest Footplates

☐	250F369	Angle Adjustable Flip-Up Footplates	$284.00

Swing-Away Legrest Calf Supports

☐	250DI42	Calf Support, Depth Adjustable, Padded *(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)*	$512.00

Swing-Away Legrest Options

☐	250DI44	Heel Loops	NC
☐	250KG4	Knee Pads, Integral Foam *(2.5" Wide x 4.5" Tall)*	$282.00

SEDEO® Lower Extremity Swing-Away Support

☐	250VD01	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	$512.00
☐	250VD02	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	$512.00

POSITIONING SUPPORTS & ACCESSORIES

	OPTION		DESCRIPTION	RETAIL

Pelvic Positioning Belts

☐	250A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐	250A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
☐	250A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	$103.00
☐	250A845	Bodypoint - Evoflex Pelvic Stabilizer Belt	$216.00
☐	250A846	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	$271.00

Chest Straps

☐	250ZR06	Chest Strap, 2-Point, Padded *(Required)*	NC
☐	250ZR07	Additional Chest Strap, 2-Point, Padded	$351.00

Chest Bars

☐	250ZR08	Chest Bar, Detachable, Round Pad	$704.00
☐	250ZR09	Chest Bar, Detachable, Flat Pad, 200 mm x 100 mm	$748.00

Accessories $

Ventilator Option

☐	250VT12	Vent Tray Fixed *(Pulmonetic LTV Models)*	$609.00

O2 Tank Option

☐	250A718	O2 Tank Holder Mounted to Backrest *(D and E Styles)*	$391.00

Wheelchair Tray Table

☐	2507A09	Tray, Swing-Away	$381.00

Backpack

☐	250A1	Black Utility	$112.00

Medial Thigh Support (N/A with Centermount Legrests)

☐	2508A80	Medial Thigh Support	$271.00
☐	2508A81	Pommel Swing-Away Hardware	$390.00

Misc. Accessories

☐	2508A59	Transfer Bar *(Pair)*	$493.00
☐	2508A52	Ram Mount Phone	$77.00
☐	2508A53	Ram Mount Tablet	$308.00
☐	2508A54	Ram Mount Cup Holder	$77.00
☐	2508A58	Rear View Mirror	$77.00
☐	250A611	Crutch Holder *(N/A w/Power Recline)*	$174.00
☐	2508A73	Bag for Armrest	$77.00
☐	2508A75	Reflectors *(Backrest shroud mounted)*	$115.00
☐	2509A10	Utility Hooks	$30.00

JOYSTICK / ELECTRONICS

	OPTION		DESCRIPTION	RETAIL
R-net Joystick *(Remote)*				**$**
☐	250J121		Standard LED Joystick *(w/o Screen Display)*	NC
☐	250J116		Colour Joystick (*2.3" Colour LCD Screen Display, CJSM1)*	NC
☐	250J118		High Resolution Colour Joystick,(CJSM2+BT) *3.5" Colour LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control*	$1,162.00
☐	250J119		Omit Remote Control *(Required with arm mounted specialty controls)*	NC
Joystick Mounting				**$**
☐	250JM34		Fixed on Armrest	NC
☐	250JM35		Parallel Swing-Away	$447.00
☐	250JM36		Link-It Mid Line Mount on Armrest *(Swing-Away)*	$447.00
☐	250JM49		Attendant *(Available with LED or CJSM1 only. Mounted to Push Handle)*	NC
Joystick Mount (Location)				**$**
☐	250HM16		Right Hand	NC
☐	250HM17		Left Hand	NC
Display Language Selection R-net				**$**
☐	250L1		English	NC
☐	250L3		French	NC
☐	250L7		Spanish	NC
Joystick Handle (Optional)				**$**
☐	250JH13		Ball Handle	$110.00
☐	250JH15		T-Handle	$110.00
☐	250JH24		Bodypoint – Rubber Dome Joystick Handle	$216.00
☐	250JH25		Bodypoint – U-Shaped Joystick Handle *(4 inch)*	$216.00
☐	250JH26		Bodypoint – Grooved Mushroom Joystick	$216.00
Electronic Programmers				**$**
☐	250PR13		R-net PC Wireless Programmer	$613.00

R-net Power Modules			$
☐	250ER11		Power Module, R-net, 120A, Non-Expandable	NC
☐	250ER10		Power Module, R-net, 120A, Expandable *Required for specialty controls or 3+ actuators*	$102.00
☐	250EX09		Expandable 120A Controller - R-net	NC
☐	250EX10		Harness for 120A Expandable Controller	NC
☐	250ER18		Power Module, R-net, 90A, Non-Expandable *(Available with No Power and Single Power)*	NC
☐	250ER17		Power Module, R-net, 90A, Expandable *Required for specialty controls or 3+ actuators*	$102.00
☐	250EX15		Expandable 90A Controller - R-net	NC
☐	250EX16		Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

OPTION		DESCRIPTION	RETAIL

Seating Control			$
Powered Seating Control			
☐	250EO40		Omit Thru-Drive *(Requires Switch Control)*	NC
☐	250EO37		Thru-Drive Control for 1 Actuator	$1,443.00
☐	250EO38		Thru-Drive Control for 2+ Actuators	$2,903.00
☐	250EO32		CxSM, Thru-Drive Control for 2+ Actuators *(Required for PRO ADVANCED)*	$2,903.00
Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*			
☐	250EO33		Dual Toggle for 1 Actuator	NC
☐	250EO34		2.5" Single Button Switch for 1 Actuator *(Button switches are not mounted)*	NC
☐	250EO41		Ctrl+5 Switch Box for 2 Actuators or More	NC
Additional Power Seating Control *(Available with Thru-Drive only)*			
☐	250QG11		Ctrl+5 Switch Box with Voice Activation	$1,078.00
☐	250QG12		Ctrl+5 Switch Box with Voice Activation *(Multi-power only)*	NC
Dual Toggle Mounting Location			
☐	250TN1		Toggle Mount, Right	NC
☐	250TN2		Toggle Mount, Left	NC
Ctrl+5 Mounting Location *(Will swing away if mounted in the same location as the joystick)*			
☐	250EY19		Right Armrest	NC
☐	250EY20		Left Armrest	NC
☐	250EY21		Back Cane, Right	NC
☐	250EY24		Back Cane, Left	NC
Ctrl+5 Switch Type			
☐	250EY22		Push Buttons	NC
☐	250EY23		Toggle Switches	$356.00

Attendant Controls			$
☐	2500A05		Dual Control Module for Attendant	$1,225.00
☐	2500A06		R-net Colour LCD (CJSM1), Attendant Control Only	$1,225.00

Attendant Mounting Location			$
☐	250DN5	Right Back Cane	NC
☐	250DN6	Left Back Cane	NC

JOYSTICK / ELECTRONICS

	OPTION		DESCRIPTION	RETAIL

R-net Programming $

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

☐	250PM12	Include Assignable Buttons	NC

Power Centermount Footplate Extension Option

☐	250IE03	Include 7.5" Powered Footplate Extension for Transfer Assist *(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.*	$592.00

R-net Omni2 Options *For Expandable Controls, choose all that apply.* $

☐	250ON14	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	$2,384.00
☐	250ON15	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.* $

☐	250DS14	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	$698.00
☐	250DS15	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	$698.00

PRO ADVANCED Programming $

Biometric Repositioning

☐	250PM03	Include Biometric Repositioning *(Requires Memory Position #1, #2, #3 below)*	NC

Memory Seating

"Memory Seating" feature allows for up to six custom seating positions that are accessed either through the joystick, Ctrl+5, joystick shortcut keys or programmable jacks on Omni2 and joystick. Choose desired preset settings below. Default postions are dealer adjustable.

☐	250PM06	Include #1 Memory Seating (Seated) *(Required)*	NC
☐	250PM07	Include #2 Memory Seating (Supine)	NC
☐	250PM08	Include #3 Memory Seating (Pressure Relief)	NC
☐	250PM09	Include #4 Memory Seating (Feet Above Heart)	NC
☐	250PM10	Include #5 Memory Seating (Transfer)	NC
☐	250PM11	Include #6 Memory Seating (C-Me)	NC

PRAR Pressure Relief Alert and Reminders

PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the R-net electronics that will prompt the user when it is time to perform a pressure relief. The R-net system will record how often the pressure relief is completed or ignored.

☐	250PM14	Include PRAR *(Minimum requirement Memory Position #1 and #3)*	NC

SPECIALTY INPUT DEVICES

	OPTION		DESCRIPTION	RETAIL
Hand Control Packages *Requires Standard Joystick Omit / Only one hand control device can be ordered.*				$
Proportional Input Devices				
☐	250EU14		Switch-It MicroGuide	$3,844.00
☐	250EU17		Switch-It VersaGuide	$2,557.00
☐	250EU18		Switch-It VersaGuide EZ	$3,023.00
Mounting Position *(Required)*				
☐	250EM04		Link-It to the Backrest *(Not Available w/ Omit Back Cushion Type)*	$1,338.00
☐	250EM05		Link-It to the Fixed Joystick Mount	$416.00
☐	250EM06		Link-It to the Parallel Swing-Away Joystick Mount	$416.00
☐	250EM07		Link-It to the Midline Swing-Away	$1,480.00
☐	250EM08		Omit Mount	NC
Mounting Platforms *(Optional)* *Not available w/ Link-It to Backrest*				
☐	250MG16		SunDisk Hand Tray *(Only available w/ MicroSeries Joysticks)*	$146.00
☐	250MG15		Bullet Hand Tray	$416.00
Mounting Location *(Required)*				
☐	250IM09		Right	NC
☐	250IM10		Left	NC
Switch Options				
☐	–		See Switch Selections *(Maximum 2)*	–
Head Control Packages *(Select headrest mount)*				$
Only one head control device can be ordered.				
Proportional Input Devices				
☐	250EU25		Switch-It DualPro® Proportional Head Array - Facial Pads	$7,057.00
☐	250EU26		Switch-It DualPro® Proportional Head Array - Spot Pads	$7,057.00
☐	250EU51		Switch-It Vigo Head Drive	$6,293.00
Switched Input Devices				
☐	250EU28		Switch-It 3-Switch Head Array - Short Curved Wings	$6,154.00
☐	250EU30		Switch-It 3-Switch Head Array - Long Curved Wings	$6,154.00
Switch Options				
☐	–		See Switch Selections *(Maximum 2)*	–
Sip 'N Puff				$
☐	250ID73		Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware *(Requires headrest mount)*	$7,037.00
☐	250ID74		Switch-It OptiStop with Sip 'N Puff	$2,741.00
Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)*				$
☐	250IJ13		Whitmyer® Linx Headrest Mount with Quick Release	$247.00
☐	250IJ14		Whitmyer® Linx-2 Headrest Mount with Quick Release	$298.00
☐	250IJ15		Whitmyer® AXYS Headrest Mount with Quick Release	$247.00

SPECIALTY INPUT DEVICES

OPTION		DESCRIPTION	RETAIL
Chin Control Packages			$
Proportional Input Devices			
☐ 250EU34		Switch-It MicroGuide *(Switches mounted to the headrest only.)*	$4,529.00
☐ 250EU36		Switch-It VersaGuide	$2,557.00
☐ 250EU37		Switch-It VersaGuide EZ	$3,023.00
Mounting Position			
☐ 250EM09		Link-It To Chin Harness, Small (150mm / 6")	$1,760.00
☐ 250EM10		Link-It To Chin Harness, Large (220mm / 8")	$1,760.00
☐ 250EM11		Link-It To Headrest (Swing-Away)	$1,760.00
☐ 250EM16		Omit Mount	NC
Input Device Mounting Location			
☐ 250IM13		Right	NC
☐ 250IM14		Left	NC
Switch Options			
☐ –		See Switch Selections	–
Safety Devices			$
☐ 250SA06		Switch-It Remote Stop Switch - XLR Connector	$1,317.00
Switched Driving Packages			$
☐ 250SW44		Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. *(Pro Spot proportional switches are required if proportional driving is desired.)* *(For non-proportional driving, single switches can be used.)*	$3,949.00

SPECIALTY INPUT DEVICES

OPTION		DESCRIPTION	RETAIL

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.
Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls
Switch ordering instructions:

1. Select switch type and quantity in Section "A" using colums below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections	**Switch / Quantity**						
Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	$
2.5" Single Switch - Black	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Green	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Red	☐	☐	☐	☐	☐	☐	$230.00
2.5" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$230.00
Egg Switch - Black	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Blue	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Green	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Red	☐	☐	☐	☐	☐	☐	$131.00
Egg Switch - Yellow	☐	☐	☐	☐	☐	☐	$131.00
MicroLite *(Includes colour labels)*	☐	☐	☐	☐	☐	☐	$195.00
Mini Button 1 - Blue *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
Mini Button 1 - Green *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
Mini Button 1 - Red *(Velcro)*	☐	☐	☐	☐	☐	☐	$125.00
1" Single Switch - Black	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Green	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Orange	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Red	☐	☐	☐	☐	☐	☐	$326.00
1" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$326.00
Electronic Switches							
Flex Switch (Headrest mount only)	☐	☐	☐	☐	☐	☐	$254.000
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Green	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Blue	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Yellow	☐	☐	☐	☐	☐	☐	$562.00
ProSpot Proportional Switch - Orange	☐	☐	☐	☐	☐	☐	$562.00
Section B – Switch Mounting Position							
Link-It to the Backrest *(Up to 1 per side)*	☐	☐	☐	☐	☐	☐	$770.00
Link-It to the Headrest - Gooseneck & Rigid *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$770.00
Link-It to the Fixed Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$587.00
Link-It to the Parallel Swing-Away Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$587.00
Link-It to the Seat Rail *(Up to 2 per side)* *(Not avail w/knee & thigh supports)*	☐	☐	☐	☐	☐	☐	$467.00
Omit Mount *(includes Velcro)*	☐	☐	☐	☐	☐	☐	NC
Section C – Switch Mounting Location							
Right side of chair	☐	☐	☐	☐	☐	☐	NC
Left side of chair	☐	☐	☐	☐	☐	☐	NC

Joystick Integrated Switches *(Pair)*

☐	250JO03	Proximity Switches on Flex Mount, Pair *(Used with VersaGuide/VersaGuide EZ)*	$1,178.00

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

Comments for special build

built4me

Q700 M with SEDEO® PRO™	$
Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Signature

Send E-Mail

SUNRISE MEDICAL

www.SunriseMedical.ca
800.263.3390
E-mail: cscanada@sunmed.com
Specifications are subject to change without notice

© 2025 Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2
Concord, Ontario L4K 4P8
10.2026 MK-129630 Rev. 28