MWD Crossover Compact Electric Wheelchair

EIPW60



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC		Phone	
Mark For				Email			


Please submit orders to: cscanada@sunmed.com

Seating Specifications			
Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options		
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" | Depth: 13-24"	STF Height Range	17", 18", 19"

Drive Base Specifications			
Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	350 lbs.	Drive Motors	4-pole 6mph
Base Size: Off-Road	Width 26" | Length 39" | Turning Radius 22.5"	Drive Wheels	14" Off-Road Drive Wheels | 14" Urban Drive Wheels
Base Size: Urban	Width 24" | Length 39" | Turning Radius 22.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24
		Ground Clearance	3"
		Color Accents	10 colors

Client Measurements

a)	Top of head to seat		in
b)	Top of shoulder to seat		in
c)	Underarm to seat		in
d)	Elbow to seat		in
e)	Top of leg to seat		in
f)	Shoulder width		in
g)	Chest width		in
i)	Hip width		in
j)	Actual sitting depth		in
m)	Foot length		in
n)	Knee to foot		in
o)	Seat width		in



Notes

✔ OPTION		DESCRIPTION	RETAIL
Seating Options			$
☐ 600SF4		**REHAB** Seat and Back - Back Canes and adjustable flat seat pan	NC
☐ 600SF5		**MPS** Seat and Back - Captain style seating	NC
☐ 600SF6		**MPS** Back and **REHAB** Seat - Captain style backrest and adjustable flat seat pan	NC
Base			
Magic 360 Base - Includes Transit and lights		End Item Number: EIPW60	$
☐ 600MW61		Magic Mobility - 360 No Power Seating Options	$ 27,547.00
☐ 600MW62		Magic Mobility - 360 Single Power Seating Options	$ 27,547.00
☐ 600MW63		Magic Mobility - 360 Multi Power Seating Options	$ 27,547.00
☐ 600MW64		Magic Mobility - 360 No Power Seating Options - MPS (Captain)	$ 27,547.00
☐ 600MW65		Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	$ 27,547.00
☐ 600MW66		Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	$ 27,547.00
Drive Wheels			$
☐ 600DW80		Off-Road Drive Wheels and Tires	NC
Optional Drive Tire Options			$
☐ 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	$ 53.00
Base Accessories			$
☐ 6008A91		Jack for Drive Tire	$ 280.00
Frame Color			$
☐ 600C185		Midnight Black	NC
☐ 600C186		Lucent Lime	NC
☐ 600C187		Sunset Red	NC
☐ 600C188		Horizon Blue	NC
☐ 600C189		Dusk Purple	NC
☐ 600C190		Vivid Pink	NC
☐ 600C191		Desert Orange	NC
☐ 600C192		Camouflage - Woods	NC
☐ 600C193		Camouflage - Desert	NC
☐ 600C171		Carbon Fiber	NC
☐ 600C194		American Flag	NC
Seat To Floor Height			$
Seat to Floor Heights start at 17" +/- 0.25"			
☐ 600HT44		Lowest Seat To Floor Height	NC
☐ 600HT45		Seat To Floor Height +1"	NC
☐ 600HT46		Seat To Floor Height +2"	NC
Batteries			$
☐ 600AG25		Group 24 Gel Battery (Pair) - Installed	$ 1,238.00
Battery Charger			$
☐ 600BJ28		10A Battery Charger	NC
☐ 600BJ38		12A Battery Charger	$ 716.00

Seating: Seat Functions

✔ OPTION		DESCRIPTION	RETAIL
No Power Seat Options			$
☐ 600ST87		No Power Seating System with Power Functionality	$ 9,608.00
Power Tilt			$
Utilized without Power Recline.			
☐ 600PA64		50° Power Tilt in Space	$ 10,543.00
Power Seat Elevate			$
Only available when Power Tilt or Power Tilt and Recline is selected.			
☐ 600PI16		12" Adjustable Seat Elevator	$ 9,185.00
Power Recline			$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.			
☐ 600PA65		150° Power Recline Only	$ 7,505.00
Power Tilt and Recline			$
☐ 600PA66		50° Tilt & 150° Power Recline	$ 16,650.00

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)			$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
☐ 600SI24		16" x 18" Vinyl Seat	NC
☐ 600SI25		18" x 18" Vinyl Seat	NC
☐ 600SI26		18" x 20" Vinyl Seat	NC
☐ 600SI27		20" x 18" Vinyl Seat	NC
☐ 600SI28		20" x 20" Vinyl Seat	NC
☐ 600SI29		22" x 20" Vinyl Seat	$ 636.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			$
REHAB Seat Pan			
Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.			
☐ 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	$ 640.00
☐ 600PN15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	$ 640.00

REHAB Seat Width

☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21" ☐ 22" ☐ 23"

REHAB Seat Depth

☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21" ☐ 22" ☐ 23" ☐ 24"

REHAB Seating Cushion Options			$
Note cushion size will be determined by Rehab Seat size selected above			
☐ 600CU59		JAY® Ion® (Sizes 16" - 21" Wide)	$ 677.00
☐ 600CU60		JAY® Ion® (Size 22" Wide)	$ 719.00
☐ 600CU61		JAY® Union® (Sizes 16" - 21" Wide)	$ 787.00
☐ 600CU62		JAY® Union® (Size 22" Wide)	$ 830.00
☐ 600CU64		JAY® Fusion® (Sizes 16" - 21" Wide	$ 1,076.00
☐ 600CU65		JAY® Fusion® (Size 22" Wide)	$ 1,224.00

Seating: Seat Sizes

✔ OPTION		DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			$
Lateral Thigh Support Pads			
☐ 600KF22		Lateral Thigh Support Pad Short - Right (5x5)	$ 77.00
☐ 600KF23		Lateral Thigh Support Pad Short - Left (5x5)	$ 77.00
Lateral Thigh Support Mounts			
☐ 600KF26		Lateral Thigh Removable Mount - Right	$ 335.00
☐ 600KF27		Lateral Thigh Removable Mount - Left	$ 335.00

Seating: Seatback Options

MPS Backrest Options - Comes standard with Manual Recline			$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	$ 466.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	$ 766.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	$ 422.00
MPS Back Accessories			$
☐ 600PH33		Push Handle - Needed for O2 holder	$ 217.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	$ 164.00
REHAB Seating Backrest Options			$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01		J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	$ 939.00
☐ 600JP02		J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	$ 1,999.00
☐ 600JP11		J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	$ 1,016.00
☐ 600JP13		J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	$ 1,100.00
REHAB Seating Backrest Options			$
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR41		15" tall shell (Upper Thoracic Med) for 22" wide	NC
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR18		18" tall shell (Shoulder Height) for 22" wide	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			$
☐ 600OD23		**REHAB** Flip Forward Backrest Option	$ 220.00

Seating: Seatback Options

✔	OPTION	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			$
JAY® Lateral Thoracic Support Pads			
☐	600JW35	Lateral Thoracic (L x H) 4" x 5" Pad - Right	$ 77.00
☐	600JW38	Lateral Thoracic (L x H) 4" x 5" Pad - Left	$ 77.00
☐	600JW37	Lateral Thoracic (L x H) 6" x 4" Pad - Right	$ 77.00
☐	600JW40	Lateral Thoracic (L x H) 6" x 4" Pad - Left	$ 77.00
JAY® Lateral Thoracic Swing-Away Mount			
☐	600JW41	Lateral Thoracic Swing Away Bracket - Right	$ 344.00
☐	600JW42	Lateral Thoracic Swing Away Bracket - Left	$ 344.00
JAY® Anterior Thoracic Support			$
☐	600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	$ 234.00
☐	600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	$ 234.00
☐	600K10	JAY® Anterior Thoracic Support - Medium ATS8302	$ 332.00
☐	600K11	JAY® Anterior Thoracic Support - Large ATS8303	$ 332.00
Pelvic Positioning Options - Select One			$
☐	600A848	Positioning Belt, Auto Style Buckle, 2" x 60"	NC
☐	600A849	Positioning Belt, Aircraft Buckle, Padded 2" x 50"	$ 106.00
☐	600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	$ 164.00
☐	600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	$ 280.00
Headrest Options - MPS Style			$
☐	600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	$ 569.00
☐	600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	$ 569.00
Headrest Options - REHAB Style			$
Headrest Pad			
☐	600HE39	Whitmyer® Plush Headrest Pad, 10"	$ 388.00
☐	600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	$ 458.00
Headrest Mount			
☐	600MH33	Whitmyer® Linx Headrest Mount with Quick Release	$ 254.00
☐	600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	$ 307.00
Armrest Style Options			$
☐	600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	$ 1,842.00
☐	600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	$ 2,615.00
☐	600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	$ 3,262.00
Armrest Pad Options			$
☐	600AB81	10" Waterfall Armrest Pad - Right	NC
☐	600AE81	10" Waterfall Armrest Pad - Left	NC
☐	600AB82	14" Waterfall Armrest Pad - Right	NC
☐	600AE82	14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✔ OPTION		DESCRIPTION	RETAIL
Center Mount Fixed Legrest			$
Center Mount Fixed Legrest			
☐ 600H210		Center Mount (14.5" - 17" Leg Length) with Flip-Up Footplate (11" x 10")	$ 628.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59		Fixed Center Mount Calf Pads	$ 718.00
☐ 600LS16		Foot Strap - 1.5" wide strap	$ 124.00
Center Mount Power Elevating Legrest			$
Center Mount Power Elevating Legrest			
☐ 600H213		Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	$ 5,750.00
Center Mount Power Elevating Legrest Extension			
☐ 600E42		Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
☐ 600E43		Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates			
☐ 600F117		Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
☐ 600F118		Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest			$
Swing-Away Legrest			
☐ 600H211		Swing-Away Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	$ 903.00
Swing-Away Power Elevating Legrest			$
Swing-Away Power Elevating Legrest			
☐ 600H214		Swing-Away Power Elevating Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	$ 6,277.00
S/A Power Legrest Positioning Options			
☐ 600HL15		Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	$ 176.00

Accessories			$
☐ 6008A93		Waterproof Wheelchair Cover	$ 229.00
☐ 6008A94		Waterproof Joystick Cover	$ 71.00
☐ 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	$ 53.00
☐ 600A221		Deluxe Backpack	$ 148.00
☐ 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	$ 296.00
☐ 600A719		Articulating Oxygen Holder for **MPS** Seating System	$ 1,537.00
☐ 600A612		Cane Holder for **MPS** Seating System - Includes push handles	$ 986.00
☐ 6008A96		Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	$ 514.00
☐ 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	$ 827.00
☐ 600A613		Cell Phone Holder - Adjustable to fit most cell phones	$ 458.00
☐ 600A614		Fishing Rod Holder	$ 458.00
Be Adaptive Products			
☐ 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	$ 2,309.00
☐ 6008A99		Be Adaptive Trigger Mechanism	$ 1,279.00
☐ 6008A100		Be Adaptive Limited Mobility Fishing Mount	$ 2,309.00

Electronics

✔ OPTION		DESCRIPTION	RETAIL
Joystick			$
☐ 600J121		Standard LED Joystick (w/o Screen Display	NC
☐ 600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
☐ 600J118		High Resolution Color Joystick (CJSM2+BT) 3.5" Color LCD Display	$ 780.00
Joystick Mount			$
☐ 600JM34		Fixed on Armrest - Depth adjustable 8.25"	NC
☐ 600JM45		Swing-Away Joystick Mount - Includes latch to lock in the forward position	$ 874.00
☐ 600JM46		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	$ 2,112.00
Joystick Mount - Location			$
☐ 600HM16		Right Hand Armrest	NC
☐ 600HM17		Left Hand Armrest	NC

Display Language Selection R-net								$
☐ 600L1	English	NC	☐ 600L3	French	NC	☐ 600L7	Spanish	NC

Joystick Handle Type			$
Comes standard with PG Cone-Shaped Joystick Handle			
☐ 600JH13		Ball Handle	$ 113.00
☐ 600JH26		Bodypoint Grooved Mushroom Joystick Handle	$ 247.00
☐ 600JH24		Bodypoint – Rubber Dome Joystick Handle	$ 264.00
☐ 600JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)	$ 264.00
Expandable Power Module			$
Must select both items			
☐ 600EX05		Expandable 120A Controller - R-net	$ 1,132.00
☐ 600EX06		Harness for 120A Expandable Controller	$ 791.00
Power Seat Function			$
☐ 600EO37		Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function	$ 1,486.00
☐ 600EO38		Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions	$ 2,990.00
Power Seat Control Option			$
☐ 600QG11		Ctrl+5 Switch Box with Free Assignable Functions	$ 1,048.00
☐ 600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)	NC
Ctrl+5 Switch Type			$
☐ 600EY22		Push Button	NC
☐ 600EY23		Toggle Switches	NC

Power Seat Control Option - Ctrl+5 Mounting location					$
☐ 600EY19	Right Armrest	NC	☐ 600EY20	Left Armrest	NC

Attendant Controls			$
☐ 6000A05		Dual Control Module for Attendant - Includes accessory mount and push handle	$ 2,534.00
Attendant Control Mounting Location			$

☐ 600DN5	Right Back Cane	NC	☐ 600DN6	Left Back Cane	NC

R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth			$
☐ 600ON14		R-net Omni2 Only Input Control Module	$ 3,719.00
☐ 600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices	NC
R-net Omni2 - Mounting Location			$

☐ 600DS23	Right Hand Armrest	NC	☐ 600DS24	Left Hand Armrest	NC

Head Control Packages			$
☐ 600EU25		Switch-It DualPro® Proportional Head Array - Facial Pads	$ 7,269.00
☐ 600EU26		Switch-It DualPro® Proportional Head Array - Spot Pads	$ 7,269.00