Quickie 5R-AADL W398

October 2026

EIR10

Mark For:

Date:
Dealer Acct #:
Dealer:
Dealer Contact:
Dealer Address:
Dealer City: Prov: PC:
Dealer Phone: ( ) Fax: ( )
Confirmation Email:
Confirm Via: o Fax o Email

Submitting for: o **Quote** o **Order**

PO#:

ADDITIONAL SHIPPING INFORMATION

Ship To:
Attention:
Address:
Address:
Ship To City: Prov: PC:
Ship To Phone: ( ) Fax: ( )

1. FRAME

Base Model

a	Option #		Definition	Price
o	129WM86		Quickie 5R (weight capacity 265 lbs)	$4,025.00

Frame Angle (see Figure 1)

a	Option #		Definition	Price
o	129FA1		75°	NC
o	129FA8		80°	NC
o	129FA2		85°	NC

Frame Inset Per Side (See Figure 2)

a	Option #		Definition	Price
o	129FF1		0"	NC
o	129FF2		1" ***(N/A w/ 12" wide)***	NC
o	129FF3		2" ***(N/A w/ 12", 13", or 14" wide)***	NC

Frame Inset	Footrest Width
0	Subtract 0 from Frame Width
1"	Subtract 2" from Frame Width
2"	Subtract 4" from Frame Width

Transit

Transit not available w/ XTR suspension & Xtender.

a	Option #		Definition	Price
o	129TO45		Transit System	NC

Frame Style (See Figure 3)

XTR Suspension N/A with Xtender

a	Option #		Definition	Price
o	129FS34		Classic Frame	NC
o	129FS35		Standard Frame (Performance)	NC
o	129FS36		XTR Suspension Frame ***Only available in Performance Frame Set Up***	$1,044.00
o	129SK1		Light Shock (up to 180 lbs.)	NC
o	129SK3		Heavy Shock (up to 265 lbs.)	NC

XTR available squeeze and RSH

Seat Depth	Squeeze (FSH - RSH)
14"	2.0" or 3.0"
15"	2.0" or 3.5"
16"	2.5" or 3.5"
17"	2.5" or 4.0"
18"	2.5" or 4.0"
19"	2.5" or 4.5"
20"	3.0" or 4.5"

XTR Tire - Available RSH	
Rear Tire	RSH
20" SHOX	13.5" thru 15.0"
20" all others	14.0" thru 15.5"
22" SHOX	14.5" thru 16.0"
22" Pneumatic	15.5" thru 16.5"
22" all others	15.0" thru 16.0"
24" SHOX	15.5" thru 16.5"
24" all others	16.0" thru 17.0"
25" SHOX	16.0" thru 17.0"
25" Marathon	16.0" thru 17.5"
26" all tires	17.0" thru 18.0"

Figure 1




Figure 2





Figure 3



Classic Frame - Traditional seat sling to front frame bend angle dimensions.

Standard Frame (Performance) - Seat Sling to front frame bend angle is 1.5" shorter than the classic position.
Front of back cane post to front edge of seat sling will remain the same for both Classic and Performance.


a	Option #	Definition	Price	a	Option #	Definition	Price

2. SIZE / SEATING

Depth Adjustable Backrest (Optional)

See diagram Depth Adjust Backrest extends the frame to the rear of the seat sling, select a factory setting. N/A on 19" and 20" Seat Depth or XTR frame.

o	129DP02		0 Setting ***Backrest at standard position, field adjustable to -1" or -2"***	$204.00
o	129DP03		-1 Setting ***Backrest -1" offset, field adjustable to flush or -2"***	$204.00
o	129DP04		-2 Setting ***Backrest -2" offset, field adjustable to flush or -1"***	$204.00



Seat Width

See Figure 2. Measure from outside of seat tube to outside of opposite seat tube. XTR and J3 Back N/A w/ 12" and 13" widths

o	129S3		12"	NC
o	129S4		13"	NC
o	129S5		14"	NC
o	129S6		15"	NC
o	129S7		16"	NC
o	129S8		17"	NC
o	129S9		18"	NC
o	129S10		19"	NC
o	129S48		20"	$308.00

Seat Depth

See Figure .1 Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered. XTR frame N/A with 12" and 13" depth

o	129SD5		12"	NC
o	129SD6		13"	NC
o	129SD7		14"	NC
o	129SD8		15"	NC
o	129SD9		16"	NC
o	129SD10		17"	NC
o	129SD11		18"	NC
o	129SD12		19"	NC
o	129SD13		20"	$641.00

Seat Sling

o	129SL1		Standard Seat Sling	NC
o	129SL6		Omit Seat Sling	NC

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out with desired width and depth

o	129JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part #		
		Width:	Depth:	

3. CASTERS / SEAT HEIGHT

Front Seat Height

See Figures 4 and 5

o	129SH14	16.0"	NC	o	129SH20	18.5"	NC
o	129SH15	16.5"	NC	o	129SH21	19.0"	NC
o	129SH16	17.0"	NC	o	129SH23	19.5"	NC
o	129SH17	17.5"	NC	o	129SH24	20.0"	NC
o	129SH18	18.0"	NC	o	129SH26	20.5"	NC
				o	129SH27	21.0"	NC

Front to Rear Seat Height Differences

Frame Depth	Max Seat Sqz
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Caster Wheels

o	129CW1		3" Micro	NC
o	129CW3		4" Low Profile Polyurethane	NC
o	129CW31		4" Micro Lighted	$110.00
o	129CW72		4" x 1" Aluminum Polyurethane	$73.00
o	129CW44		4" x 1.25" Semi Pneumatic	NC
o	129CW64		4" x 1.5" Aluminum Soft Roll	$110.00
o	129CW5		5" Low Profile Polyurethane	NC
o	129CW51		5" Micro Lighted	$110.00
o	129CW6		5" Soft Roll (Grey) ***N/A w/ Frog Legs***	$126.00
o	129CW65		5"x 1" Aluminum Polyurethane	$126.00
o	129CW7		5"x 1.5" Semi Pneumatic	$110.00
o	129CW66		5"x 1.5" Aluminum Soft Roll	$110.00
o	129CW34		6" Polyurethane	NC
o	129CW35		6"x 1.5" Semi Pneumatic	$110.00
o	129CW68		6"x 1.5" Aluminum Soft Roll	$110.00

Caster Fork

See Figures 5 and 6. Figure 6 in appendix pg 6.

o	129CF65		Multi Position Fork	NC
o	129CF69		Single Side Fork	$197.00

Caster Fork Suspension

o	129CF53		Frog Legs 0 weight to 100 lbs	$507.00
o	129CF54		Frog Legs 101 lbs weight to 190 lbs	$507.00
o	129CF55		Frog Legs 191 lbs weight to 265 lbs	$507.00

Figure 4




Figure 5

	Available FSH per Caster Fork		
Caster Size	Short (Multi, Single Sided)	Standard (Multi, Single Sided)	Frog Legs
3"	16"-20"	16.5"-21"	17"-21"
4"	16"-20.5"	16.5"-21"	17"-21"
5"	17,18,19,20,21"	17"-21"	18"-21"
6"x1"	N/A	18"-21"	19,20,21"
6"x1.5"	N/A	18"-21"	19,20,21"

4. REAR WHEELS / SEAT HEIGHT

Rear Seat Height

Please reference RSH-Axle plate-wheel diagram on bottom of the page.

a	Option #		Definition	Price
o	129RH11		13.0"	NC
o	129RH12		13.5"	NC
o	129RH13		14.0"	NC
o	129RH14		14.5"	NC
o	129RH15		15.0"	NC
o	129RH16		15.5"	NC
o	129RH17		16.0"	NC
o	129RH18		16.5"	NC
o	129RH19		17.0"	NC
o	129RH26		17.5"	NC
o	129RH20		18.0"	NC
o	129RH27		18.5"	NC
o	129RH21		19.0"	NC
o	129RH37		19.5"	NC
o	129RH22		20.0"	NC

Camber

0,3,6 available with std frame. 0,4,8 available with XTR frame. 3 N/A w/Suspension.

a	Option #		Definition	Price
o	129CI1		0°	NC
o	129CI2		3°	NC
o	129CI3		4°	NC
o	129CI4		6°	NC
o	129CI5		8°	NC

Axle Plate - Wheel - Available RSH

Axle Plate	20" Wheel	22" Wheel	24" Wheel	25" Wheel	26" Wheel
Std Short	13"-17"	13.5"-18"	14.5"-19"	15"-19.5"	15.5"-20"
Std Long	13"-18.5"	13.5"-19.5"	14.5"-20"	15"-20"	15.5"-20"
Amp Short	13.5"-17"	14.5"-18"	15.5"-19"	16"-19.5"	16.5"-20"
Amp Long	13.5"-19"	14.5"-20"	15.5"-20"	16"-20"	16.5"-20"
Susp	14"-17.5"	15"-18.5"	16"-19.5"	16.5"-20"	17"-20"
Xtender	N/A	17.5"-19.5"	18.5"-20"	N/A	N/A

*Note: other factors such as CG or Transit may limit RSH

Xtender Power Assist

RSH 17.5"-19.5" with 22" wheels, RSH 18.5"-20" with 24" wheels 0 camber, N/A w/ Amputee Axle. 22" & 24" wheels only. Not available with XTR or standard rear suspension. Anti-Tips no charge with Xtender. Dual Drive mode available w/B4ME

a	Option #		Definition	Price
o	129XT6		**NICKEL**	$9,996.00
o	129XT7		22" Xtender Power Assist JWX-2 **LITHIUM**	$11,604.00
o	129XT4		24" Xtender Power Assist JWX-2 **NICKEL**	$9,996.00
o	129XT5		24" Xtender Power Assist JWX-2 **LITHIUM**	$11,604.00

Please select one of the following handrim options

a	Option #		Definition	Price
o	129XH5		Plastic Coated	NC
o	129XH6		Stainless Steel	NC

Xtender - Other Options

a	Option #		Definition	Price
o	129XA10		Spare NiMH Power Pack - NICKEL	$1,925.00
o	129XA11		Spare Li-ion Power Pack - LITHIUM	$2,479.00

Xtender Wheel Options

a	Option #		Definition	Price
o	129XA8		SPARE 22" Lite Spoke Pneu Alum	NC
o	129XA9		SPARE 22" Lite Spoke Pneu Plastic	$144.00
o	129XA2		SPARE 24" Lite Spoke Pneu Alum	NC
o	129XA4		SPARE 24" Lite Spoke Pneu Plastic	$144.00

Rear Wheel Size

Ultra Lightweight Wheels only available with aluminum anondized handrim

Contact Built 4 Me for Spinergy color options

a	Option #		Definition	Price
o	129RW103		20" Mag Wheel 5 Spoke	NC
o	129RW104		22" Mag Wheel 5 Spoke	NC
o	129RW102		24" Mag Wheel 5 Spoke	NC
o	129RW90		20" Lite Spoke Wheel	$234.00
o	129RW91		22" Lite Spoke Wheel	$234.00
o	129RW92		24" Lite Spoke Wheel	NC
o	129RW93		25" Lite Spoke Wheel	NC
o	129RW94		26" Lite Spoke Wheel	$234.00
o	129RW06		22" Ultra Lightweight Spoke Wheel	$292.00
o	129RW10		24" Ultra Lightweight Spoke Wheel	$292.00
o	129RW28		22" Spinergy SPOX Black	$985.00
o	129RW95		24" Spinergy SPOX Black	$985.00
o	129RW96		25" Spinergy SPOX Black	$985.00
o	129RW97		26" Spinergy SPOX Black	$985.00
o	129RW32		22" Spinergy LX Black	$1,237.00
o	129RW98		24" Spinergy LX Black	$1,237.00
o	129RW99		25" Spinergy LX Black	$1,237.00
o	129RW40		26" Spinergy LX Black	$1,237.00
o	129RW50		24" Mountain	$292.00

Rear Wheel Tire

See Availability Matrix on PAGE 4.

a	Option #		Definition	Price
o	129RT1		Pneumatic	NC
o	129RT2		Pneumatic w/ Airless Insert ***N/A with Xtender***	$100.00
o	129RT3		Full Poly	NC
o	129RT38		SHOX Performance Poly	$110.00
o	129RT39		High-Pressure Clincher Black	$110.00
o	129RT21		Schwalbe Marathon Plus Evolution	$164.00
o	129RT27		Knobby V- Trak by Kenda	$110.00
o	129RT12		Mountain Tire ***Only available in 24"***	$110.00

Handrim

Plastic is not available on all wheels in all sizes. Natural Fits available in 22" 24" and 25". Natural Fit 5 spoke MAG only in 24". Handrim Selection N/A w/ Xtender.

See matrix figure 7 in appendix page 6 for reference.

a	Option #		Definition	Price
o	129HR1		Aluminum Anodized	NC
o	129HR4		Plastic Coated ***N/A w/ 25" & 26" Lite Spoke Wheels***	$129.00
o	129HR22		Natural Fit - Standard Grip	$549.00
o	129HR23		Natural Fit - Super Grip	$611.00
o	129HR39		Ellipse 3R	$550.00
o	129HR40		Carbolife - Gekko	$575.00
o	129HR41		Carbolife - Curve L Tetra Grip	$625.00
o	129HR28		Q-Grip Neoprene Coated	$475.00
o	129HR8		8 Oblique Projections ***N/A w/ SHOX tires, only available with 24" Lite Spoke***	$146.00
o	129HR7		Omit Handrim ***N/A w/ Ultra Lightweight Wheels***	NC

4. REAR WHEELS / SEAT HEIGHT (CONT)

Axle

a	Option #		Definition	Price
o	129AX1		Quick Release Steel	NC
o	129AX13		Quad Quick Release ***Only available on 24" Mountain rear wheels***	$182.00

Axle Plates

Selection not required for XTR Suspension

a	Option #		Definition	Price
o	129AL40		***Some limitations with STFH***	NC
o	129AL43		Standard Long (COG 0.5" to 4")	NC
o	129AL100		Axle plate amputee short	$272.00
o	129AL101		Axle plate amputee long	$272.00
o	129AL38		Suspension (COG -0.5 to +3) ***N/A with XTR Suspension Frame***	$654.00
o	129AL39		Xtender (COG 1"-4")	NC

Rear Suspension Rating

Selection not required for XTR Suspension. N/A with Xtender

a	Option #		Definition	Price
o	129SV04		From 0 weight to 100 lbs (Yellow)	NC
o	129SV05		From 101 lbs weight to 190 lbs (Red)	NC


Center of Gravity

Measure from front of back post to center of rear axle. XTR Suspension Frame will be factory set at 1.5", Xtender 1"- 4"

a	Option #		Definition	Price
o	129CC29		-1.5"	NC
o	129CC11		-1"	NC
o	129CC10		-.5"	NC
o	129CC12		0"	NC
o	129CC1		.5"	NC
o	129CC2		1"	NC
o	129CC3		1.5"	NC
o	129CC4		2"	NC
o	129CC5		2.5"	NC
o	129CC6		3"	NC
o	129CC7		3.5"	NC
o	129CC8		4"	NC

	Wheel				
Tire	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain
Pneumatic	20,22,24	20,22,24,26		22,24,26	
Pneumatic w/ Airless Insert	20,22,24	20,22,24		22,24	
Full Poly	20,22,24	20,22,24			
SHOX Performance Poly	20,22,24	20,22,24,25		24,25	
High-Pressure Clincher Blk		22,24,26	22,24	22,24,26	
Schwalbe Marathon		24,25,26		24,25,26	
Knobby V Trak by Kenda		24		24	
Mountain Tire					24
20"= 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590mm					

Figure 9



5. WHEEL LOCKS

Wheel Lock

Not all wheel locks available in all seat depths, wheel sizes and COG settings. XTR suspension only available w/ high mount push

a	Option #		Definition	Price
o	129WL46		High Mount Push	NC
o	129WL50		High Mount Pull	NC
o	129WL68		Ergo Scissor	NC
o	129WL64		Compact	$162.00
o	129WL10		Omit Wheel Lock	NC
o	129WL12		Do Not Mount	NC

Extension Handle

Only available w/ High Mount Push or Pull Wheel Locks

a	Option #		Definition	Price
o	129EH7		Extension handles (Pair)	$88.00

Grade Aid

Only avail w/ high mount push or pull. N/A w/ clincher, poly, SHOX, Schwalbe, or Xtender

a	Option #		Definition	Price
o	129GA1		Grade Aid	$175.00

6. FOOTREST OPTIONS

Footrest

Angle Adjustable Flip Back N/A w/ 19-20 wide & w/ 0" inset. Some Flip back configurations won't allow 2 " ground clearance and will need to raise up. Flip back not available 1" inset & 5"- 6" casters. See matrix in Appendix

a	Option #		Definition	Price
o	129F229		Aluminum Tubular	NC
o	129F274		Aluminum Tubular w/ Snap On Cover	$43.00
o	129F276		Angle Adjustable Flip Back	$247.00
o	129F51		Aluminum Platform Footrest	$110.00
o	129F278		High Mount Platform	$110.00

Adaptive Platform Spacer

Only available with Platform footrest & Angle Adj. Flip Back.

a	Option #		Definition	Price
o	129AD1		Adaptive Platform 2" Spacer	$110.00
o	129AD2		Adaptive Platform 4" Spacer	$110.00
o	129AD3		Adaptive Platform 6" Spacer	$110.00

Footrest Options

a	Option #		Definition	Price
o	129DI03		Leg Strap	$47.00
o	129DI04		Extra Leg Strap	$47.00
o	129DI05		Gel Hanger Cover	$123.00
o	129DI07		Neoprene Impact Guards	$105.00
o	129DI02		Wear Pads	$58.00
o	129DI06		Position Plates ***(N/A for Tubular Footrest)***	$58.00

a	Option #		Definition	Price

7. BACKREST

Back Type

o	129BT3		Non-Folding Angle Adjustable	NC
o	129BT82		Folding w/ Double Lock Down Angle Adj	NC

Height Adjustable Backrest

o	129B110		10 "- 14" Back Height	NC
			o 10 o 11 o 12 o 13 o 14	
o	129B127		12 "- 16" Back Height	NC
			o 12 o 13 o 14 o 15 o 16	
o	129B111		14 "- 18" Back Height	NC
			o 14 o 15 o 16 o 17 o 18	
o	129B137		16 "- 20" Back Height	NC
			o 16 o 17 o 18 o 19 o 20	

Push Handles

o	129PH6		Integral Push Handles	NC
o	129PH5		Fold-Down Push Handles ***N/A w/ 10-14" Back Height or Standard Upholstery.***	$168.00
o	129PH2		Omit Push Handles	NC

JAY J3 Back

o	129JB01		JAY J3 Back	See OF

Back Upholstery

o	129BU1		Standard Upholstery	NC
o	129BU88		Standard Upholstery EXO	$58.00
o	129BU70		Tension Adj Ballistic Nylon ****Reference 186BU70 for billing***	$367.00
o	129BU71		Tension Adj 3DX Vented ****Reference 186BU71 for billing***	$387.00
o	129BU89		Backrest Tension Adj EXO ****Reference XXXBU89 for billing***	$387.00
o	129BU14		Omit Upholstery ***(Use w/ JAY Backs)***	NC

8. ARMREST / SIDEGUARDS

Armrest - Padded Swing Away

Includes angle adj receiver hardware. Some limitations may exist w/ integral push handles, some RSHs and wheel sizes

o	129AA18		Padded Swing-Away STD (9" to 14")	$278.00
	Setting	o 9" o 10" o 11" o 12" o 13" o 14"		

Armrest - Height Adjustable

o	129AA11		Single Post Ht. Adjust. Std ***(8.5"-13.5")***	$339.00
o	129AA29		Single Post Ht. Adjust. Low ***(7"-11")***	$339.00

Armpad - For Height Adjustable

Recomend using classic style for tray tables

o	129AM23		Desk Length (10") ***Classic***	NC
o	129AM21		Desk Length (10") ***Contour***	NC
o	129AM24		Full Length (14") ***Classic***	NC
o	129AM22		Full Length (14") ***Contour***	NC

Side Guards

Some interferance my exist w/ padded SA arms at lowest settings

o	129SG33		Plastic - Standard	$169.00
o	129SG29		Plastic - Low	$169.00
o	129SG35		Aluminum ***N/A w/ push, pull or compact wheel lock. N/A w/ XTR suspension***	$286.00
o	129SG36		Aluminum w/ Fender ***N/A w/ push, pull or compact wheel locks. N/A w/ XTR suspension***	$321.00
o	129SG12		Carbon Fiber ***N/A w/ push, pull or compact wheel locks***	$537.00

Color For Aluminum Side Guard

o	129GE01		Black	NC
o	129GE02		Matching Frame - Standard Colors	NC
o	129GE03		Matching Frame - Personalized Color	$58.00

9. KEY OPTIONS AND ACCESSORIES

Anti Tips

Anti Tip Single Side N/A 12W-14W

Anti Tip Single Side N/A with XTR, Xtender and Suspension Axles.

o	129AT1		Anti Tips (Pair)	$121.00
o	129AT3		Omit Rear Anti Tip	NC
o	129AT14		Anti Tip Single Side - Flip Up	$259.00
			o Left o Right (Facing from rear)	

Pack Accessories

o	129A1		Black Utility Backpack	$112.00
o	129A2		Back Pack Black/Charcoal Kids	$105.00
o	1291A01		Seat Pouch Black	$83.00

Positioning Belt

o	129A826		Belt 1.5" Auto Buckle Universal	$89.00
o	129A802		Belt 2" Aircraft Buckle	$139.00
o	129A803		Belt 2" Aircraft Buckle Padded	$204.00
o	129A804		Positioning Belt Clamps Only	$77.00

Spoke Guards

Only available with Lite spoke. N/A w/ Xtender Power Assist

o	1294A01		Quickie Black Spoke Guard	$134.00
o	1294A02		Quickie Clear Spoke Guard	$134.00

Miscellaneous

o	129A2100		Crutch Holder ***Only available with Integral or Fold Down push handle. Not available with XTR Suspension***	$227.00
o	1298A22		Bolt-On Push Handles	$134.00
o	1298A21		Tool Kit	$63.00
o	1298A05		Gel-Pad swing away arm	$175.00
o	129A999		Quickie Water Bottle & Holder	$63.00

Upper Extremity Support Tray

o	129A29		X-Small (Suggested for 10"-12")	$298.00
o	129A30		Small (Suggested for 13"-14")	$298.00
o	129A31		Medium (Suggested for 15"-17")	$311.00
o	129A32		Large (Suggested for 18"-20")	$311.00

a	Option #	Definition	Price	a	Option #	Definition	Price

10. COLOR

Standard Frame Colors

o	129C211		Celestial Sapphire	NC
o	129C212		Desert Sage	NC
o	129C213		Desert Tan	NC
o	129C214		Gloss Black	NC
o	129C215		Matte Black	NC
o	129c216		Metallic Shadow	NC
o	129C217		Midnight Blue	NC
o	129C218		Mystic Purple	NC
o	129C219		Petal Pink	NC
o	129C220		Pink Power	NC
o	129C221		Red Pepper	NC
o	129C222		Royal Blue	NC
o	129C223		Salute Green	NC
o	129C224		Scarlet Noir	NC
o	129C225		Silver Ore	NC
o	129C226		Solar Orange	NC
o	129C227		Starlight White	NC
o	129C228		Storm Gray	NC
o	129C229		Turbo Green	NC
o	129C230		Velocity Yellow	NC

Personalized Frame Color

o	129C112		Desert Camouflage	$460.00
o	129C43		Mossy Oak Camouflage	$460.00
o	129C42		Stars & Stripes	$460.00
o	129C108		Zebra ***(Choose above base color) Starlight White std N/A w/ Matte Black, or Gloss Black)***	$460.00

Anodized Color Package

Package includes only forks, backcane bracket, axle sleeve clamp, footrest clamp, Axle plate plug.

Frog Legs, Xtender, suspension axle plates black only

o	129EB02		Black	NC
o	129EB06		Blue	$187.00
o	129EB05		Lime Green	$187.00
o	129EB04		Purple	$187.00
o	129EB03		Red	$187.00
o	129EB01		Ti Grey	$187.00
o	129EB13		Gold	$187.00

Built For Me

To inquire about any unique personalized Quickie products, contact:
Phone: (800) 263-3390, Email: cscanada@sunmed.com

built4me

ON-CHAIR BUILD-READY SPECIALIZED UNIQUE

11. APPENDIX

Figure 6

	Footrest-Caster Fork Matrix		Tubular/Flat			Platform			AA Flipback		
	Caster Size	Fork	Inset			Inset			Inset		
			0"	1"	2"	0"	1"	2"	0"	1"	2"
75° Frame Angle	3"	All Forks							N/A		
	4"	All Forks							N/A		
	5"	All Multi/Carb/Sngl	N/A						N/A	N/A	
		Froglegs	N/A			.5"			N/A	N/A	
	6"	Multi/Carb/Sngl (tall only)	N/A			2"			N/A	N/A	
		Froglegs	N/A	N/A		N/A	.5"		N/A	N/A	N/A
80° Frame Angle	3"	All Forks							N/A		
	4"	All Multi/Carb/Sngl							N/A		
		Froglegs	N/A						N/A	N/A	
	5"	All Multi/Carb/Sngl	N/A			1"			N/A	N/A	
		Froglegs	N/A			2"			N/A	N/A	
	6"	Multi/Carb/Sngl (tall only)	N/A	1"		N/A	1"		N/A	N/A	
		Froglegs	N/A	N/A		N/A	2"		N/A	N/A	N/A
85° Frame Angle	3"	All Forks							N/A		
	4"	All Multi/Carb/Sngl	N/A			.5"			N/A		
		Froglegs	N/A			.5"			N/A	N/A	
	5"	Multi/Carb (short only)	N/A	.5"		2"	.5"		N/A	N/A	
		All Wide, Sngl Sd, or Tall	N/A	1"		2"	.5"		N/A	N/A	
		Froglegs	N/A	1.5"		N/A	1.5"		N/A	N/A	
	6"	Multi/Carb/Sngl (tall only)	N/A	2"		N/A	2"		N/A	N/A	
		Froglegs	N/A	N/A	1.5"	N/A	N/A	1.5"	N/A	N/A	N/A

(blank) = Go N/A = No Go # = "Soft" Go with limited seat angle (FSH-RSH). Min seat angle listed in cell.

Figure 7

	Wheels					
Handrims	5 Spoke Mag	Lite Spoke	Ultra Lightwght	Spinergy LX	Spinergy SPOX	Mountain
Aluminum Anodized	20,22,24	20,22,24,25,26	22,24	22,24,25,26	22,24,25,26	24
Plastic Coat	20,22,24	22,24		24,26	24,26	
Natural Fit - Standard Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit - No Thumb Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit - Super Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - Std Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - Super Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Natural Fit LT - No Thumb Grip	24	22,24,25,26		22,24,25,26	22,24,25,26	
Surge - Oval w/Gription Strip		22,24,25,26		22,24,25,26	22,24,25,26	
Surge LT - Oval w/Gription Strip		22,24,25,26		22,24,25,26	22,24,25,26	
Q Grip Neoprene Coated		24,25		24,25	24,25	
Spinergy Flex Rim				24,25		
8 Oblique Projection		24				
Omit Handrim	20,22,24			22,24,25,26	22,24,25,26	24

SUNRISE MEDICAL

Sunrise Medical Canada Inc.
1000 Creditstone Road, Unit 2, Concord, ON L4K 4P8